



TECHNICAL DATASHEET

Item code	671	
Name	EXTRA ERGOLIFT – Folding motorized evacuation stair chair	
Dimensions [mm]	Length	840 ÷ 1500
	Width	550
	Height	1040 ÷ 1630
	Folded height	1040
	Folded thickness	290
	Seat height	500
Weight [kg]	32	
Max load [kg]	230	
Used materials	Aluminium alloy, ABS	
Wheels	n°2 Ø200 mm (fixed) n°2 Ø125 mm (rotating w/brake)	
Power source	Engine	Li-ion battery 28V 5Ah
	Console	nr. 2 batteries type AA 1,5V
Charging time	Max 1 h	
Declared use	Manual or electro-assisted transportation of seated patients with use-compatible pathologies who are temporarily or permanently disabled. The product is NOT fit for use inside rescue vehicles with patient on board	
Conformity	MDR 2017/745 - Class I Medical Device EN 1789 Certified	



2309

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it