



TECHNICAL DATASHEET

Item code	16150/G	16150/V
Name	ERGON - Variable Geometry Comfort Scoop Stretcher	
	Yellow/Blue	Green/Blue
Dimensions [mm]	Length	1650 ÷ 2010
	Folded length	1200
	Width	430
	Height	75
	Folded height	90
Used materials	Al alloy, Polypropylene, Thermoplastic elastomer	
Weight [kg]	10,5	
Max load [kg]	170	
Temperature [°C]	Working	-10 ÷ +70
	Storage	-50 ÷ +125
Relative humidity [%]	Working	5 ÷ 95
	Storage	
Characteristics	X-ray transparent in the middle zone	
Declared use	Transport, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible	
Conformity	MDD 93/42/CEE – Class I Medical Device EN 1865-1 – EN 1789 certified	



2008

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it