



## TECHNICAL DATASHEET

<b>Item code</b>	<b>670/BR</b>	
<b>Name</b>	<b>"EXTRA" STAIR CHAIR – Foldable track assisted stair chair with armrest EN 1865 certified</b>	
<b>Dimensions [mm]</b>	Length	810 ÷ 1390
	Width	550
	Height	1010
	Folded thickness	290
	Seat height	500
<b>Weight [kg]</b>	20,3	
<b>Max load [kg]</b>	230	
<b>Used materials</b>	Al alloy, ABS, S235JR	
<b>Wheels</b>	n°2 Ø200 mm (fixed) n°2 Ø125 mm (rotating w/brake)	
<b>Declared use</b>	Manual transportation of seated patients with use-compatible pathologies who are temporarily or permanently disabled. . The product is NOT fit for use inside rescue vehicles with patient on board	
<b>Conformity</b>	MDR 2017-745 – Class I Medical Device EN 1865-4 – EN 1789 certified	



2309

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)