



## TECHNICAL DATASHEET

<b>Item code</b>	<b>7077/4RG PROOF</b>
<b>Name</b>	<b>MERCURY 5 - Stainless steel certified high version 5 position self loading stretcher with 4 rotating wheels</b>
<b>Dimensions [mm]</b>	Length 1970
	Length with folded trolley 1670
	Width 570
	Height 1010 – 910 – 820 – 700 – 390
	Folded height 310
	Loading height 690 ÷ 710
<b>Used materials</b>	AISI 304, aluminium, HDPE
<b>Wheels</b>	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm
<b>Weight [kg]</b>	49,5
<b>Max load [kg]</b>	250
<b>Declared use</b>	Transporting, in the lying down position, of temporarily or permanently disabled patients with pathologies that allow this. The stretcher is suitable for transporting patients on vehicles appropriately prepared with fastening systems
<b>Conformity</b>	MDR2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified



2308

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)