



13010 Icarus

Multi level stretcher with side rails



FEATURES

Multi level stretcher ICARUS is built in stainless steel, has 4 diam. 125 mm wheels (2 with brakes), side rails and foot rest. Four gas pistons allow a perfect adjustment to all necessary heights for a practical and professional work. Available in following versions with:

- ABS lying part for a perfect cleaning and disinfection (Art. 13010)
- stainless steel tank (Art 13012).

MAIN FEATURES

Size	Adults
Usage	Transport
Multilevel	
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Min height (cm)	42
Max height (cm)	88
Width (cm)	62
Length (cm)	191
Weight (kg)	29,5
Loading capacity (kg)	170
Material	Stainless steel ABS

IMAGES



SEE ALSO



Cronos 13018
Multi level stretcher with back rest
and Trendelenburg



Icarus 13012
Multi level stretcher with steel tank

ACCESSORIES



1880-001
Mattress for Icarus



1880-002
Mattress pillow to be used with
Icarus stretcher



607/B
Blue belt in 2 pieces with plastic
buckle



607/Arneb
Orange belt in two pieces with
plastic buckle and Me.Ber. band
height 50mm



1880-001
Mattress for Icarus