



# 9130 - ME.BER. BABY FIX


Fixing system with belts for transport of children inside the ambulance



## FEATURES

Me.Ber. Baby Fix is a device for transport of children, adaptable to all model of stretcher by mean of 6 belts. Suitable for newborn and children up to 130 cm and 35 kg, allows a safe transport thanks to the restraint system at 5 points. Supplied with cover art. 12034. **Certified UNI EN 1789.** Length 74 x Head width 39 x Feet width 20 cm

## MAIN FEATURES

Pediatric	
Length (cm)	74
Width (cm)	39
Weight (kg)	1,3
Usage	Transport

---

## SEE ALSO

---



**9128**

**Me.Ber. Vacuum Baby Fix**

Vacuum mattress for fixing system  
with belts for transport of children  
inside the ambulance

---

---

## ACCESSORIES

---



**12034**

Red empty carrying bag for "Baby  
Fix" and "Vacuum Baby Fix"

---