



630/N Maxima

Black painted anti trauma scoop stretcher



FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.












Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

- Effective hooking and unhooking device to open/close the stretcher.
- Adjustable length with automatic locking device.
- Bag with n. 3 quick unhooking safety belts to immobilize the patient (603/MEB).
- Easy arrangement inside the ambulance when folded.
- Completely built with light-weight and sturdy aluminium.






MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	
Material	Painted aluminium

SEE ALSO

			
Maxima 630 Grey anodized anti trauma scoop stretcher	Maxima 630/A Orange painted anti trauma scoop stretcher	Maxima 630/B Blue painted anti trauma scoop stretcher	Maxima 630/G Yellow painted anti trauma scoop stretcher
			
Maxima 630/N Black painted anti trauma scoop stretcher	Maxima 630/N Black painted anti trauma scoop stretcher	Maxima 630/V Green painted anti trauma scoop stretcher	Maxima 631/G Yellow painted anti trauma scoop stretcher with head/trunk in one piece
			
Maxima 631 Grey anodized anti trauma scoop stretcher with head/trunk in 1 pc			

ACCESSORIES

			
616 Certified vertical fastening system for scoop stretchers	616/A Certified vertical fixing system for scoop stretchers	614 Certified horizontal fastening system for scoop stretchers	630/CS Carrying case for scoop stretchers
			
603/MEB Nylon bag containing 3 belts Art. 607/Ameb			