



## 19118

350 mm ramp for stretcher supports



### FEATURES

350 mm ramp with raised sides and gas spring.

### MAIN FEATURES

Length (cm)	35
Width (cm)	62
Height (cm)	5,4
Weight (kg)	4,3

### ACCESSORIES



#### 19130

"MTX" single space stretcher support with lateral movement with ramp mm. 180 - EN 1789 certified



#### 19130-955

"MTX-955" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 955 - EN1789 certified



#### 19130-956

"MTX-956" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 956 - EN1789 certified



#### 19136

"MX" fixed single space stretcher support with ramp of mm 180 - EN1789 certified



#### 19136-955

"MX-955" fixed single space stretcher support with ramp of mm. 180, equipped with certified fixing system art. 955 - EN1789 certified



#### 19136-956

"MX-956" fixed single space stretcher support with ramp of mm. 180, equipped with en1789 certified fixing system art. 956 - EN1789 certified

---

## ACCESSORY OF

---



**19130**

"MTX" single space stretcher  
support with lateral movement  
with ramp mm. 180 - EN 1789  
certified

---