



6200/V Maxima

Scoop stretcher for veterinary use



FEATURES

Completely constructed in aluminium, the structure is a break-apart design, dividing into two parts, allowing the stretcher to be positioned under the patient without having to remove it. It is equipped with two quick-release buckles.

MAIN FEATURES

Size	Children
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	130
Weight (kg)	7,8
Loading capacity (kg)	170

SEE ALSO



Maxima 6200
Pediatric anti trauma scoop
stretcher



6300/B
Vacuum mattress for veterinary use



6100
Long folding stretcher



6150
Short folding stretcher for
veterinary use
