



## 962

EN 1789 certified fixing system for stretcher/chair "EVOX"



### FEATURES

Fixing system for self loading stretcher/chair EVOX art. 7112 Proof is entirely made of stainless steel to ensure corrosion resistance and high durability. Due to its sturdy structure is able to withstand shocks of up to 10G in all directions, as regulated by the UNI EN 1789 norm.



Device certified according to the European harmonized safety standards UNI EN 1789.

### ACCESSORY OF



**Evox Proof 7112**  
"EVOX" certified chair cot with  
mattress



**Evox Proof 7112I**  
"EVOX" certified self-loading  
stretcher convertible into chair  
complete with mattress