



12092

Battery MILWAUKEE M28 RED 5,0 Ampere for Extra Ergolift chair



FEATURES

28V battery

Autonomy: 18 ramps by 20 steps with a load of 130 kg

Charging time: max 1 hour

ACCESSORY OF



Extra Ergolift 671
Powered stair chair



Extra Ergolift FX 675
Fixed frame powered Stair Chair
suitable for transportation



Extra Ergolift 671
Powered stair chair



Extra Ergolift FX 675
Fixed frame powered Stair Chair
suitable for transportation
