



8432-001

Pedestal for art. 8432



FEATURES

Cabinet pedestal for art. 8432. Back holes for wall hanging, lower holes for fixing it to the ground and interlocking system for cabinet to be assembled to the pedestal. Made of carbon steel enamelled by epoxy polyester white powder RAL9010 heated at 180°.

MAIN FEATURES

Width (cm)	25
Height (cm)	100
Depth (cm.)	15

SEE ALSO



8396-001

Pedestal for art. 8396 – 8398

ACCESSORY OF



8432

Teca di colore bianco mod.
standard per defibrillatore
