



# 9114 Chesterfield


Vacuum mattress with independent chambers with central welding



## FEATURES

CHESTERFIELD vacuum mattress with independent chambers with central welding.  
Usable with X-ray, CT and MRI.

## MAIN FEATURES

Trauma	
Length (cm)	210
Width (cm)	100
Weight (kg)	7,8
Loading capacity (kg)	150

## IMAGES



---

SEE ALSO

---



**Snake Snow 894/SNOW**  
Independent chambers vacuum  
mattress with anti-slip system for  
snow



**"Flake" 891**  
Independent chambers vacuum  
mattress



**891/SG**  
Flake basic independent chambers  
vacuum mattress



**Snake 894**  
Indipendent chambers vacuum  
matress



**Naxos 9104**  
"NAXOS" orange and blue vacuum  
mattress without chambers