



1322

Rescue knapsack ATENE



FEATURES

ATENE is a compact, functional and light knapsack. Internally, the central compartment has Velcro® stripes for holding six standard coloured bags. Equipped with a frontal pocket and two small double zip side pockets, this backpack is completed by three reflective high visibility strips and frontal and lateral elastic bands to accomodate different kind of medical material.

MAIN FEATURES

Width (cm)	40
Height (cm)	50
Depth (cm.)	24
Weight (kg)	1,6

IMAGES



SEE ALSO



1390
Red cordura backpack "Yankee"



Rescuemed 1343
Rescue knapsack



1325/RG
Compact professional royal and
grey rescue backpack Antartide



1325/RB
Compact professional rescue
backpack Antartide



1326/RG
Compact professional rescue
backpack Artide



1362
"Ambumed" rescue knapsack