



## 633 Axim

Grey anodized anti trauma scoop stretcher with head/trunk in one piece



### FEATURES

Thanks to its features the scoop stretcher has an important role in first aid as it can be put under the patient without removing him. The shape of the stretcher also enables the perfect immobilization of the patient during his transport.

- Once folded easy arrangement inside the ambulance. Supplied with a bag with 3 belts.

•



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

### MAIN FEATURES

Length (cm)	168
Max length (cm)	211
Folded length (cm)	120
Width (cm)	42
Height (cm)	6
Folded height (cm)	9
Weight (kg)	8,3
Loading capacity (kg)	170