




## 12200

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit



### FEATURES

The bariatric wings are a practical accessory that allows you to convert the self-loading stretcher into a practical bariatric recovery tool. Thanks to the small footprint (and a weight of only 14 kg) they can be easily transported by ambulance to be used when needed. Their use allows to face in safety, and in the most adequate way possible, every recovery scenario in the presence of bariatric patients. The kit includes two wings that extend the support surface of the stretcher by over 40 cm; assembly and disassembly on the self-loading stretcher are extremely simple and intuitive and do not require the use of tools. The bariatric wings are equipped with mattresses (art. 12170 for right bariatric wing and art. 12172 for left bariatric wing), belts with reel and assembly kit art. 12080-001. They can be purchased in kits including both the right and the left (art. 12200) or separately (art. 12202 Right bariatric wing and 12204 left bariatric wing).


For Frog and MBS-01 ambulance cots use art. [12084](#).

Dimensions:  
1365x205x85 mm (a wing)

### MAIN FEATURES

Length (cm)	136,5
Width (cm)	20,5

### IMAGES



---

SEE ALSO

---



**12084**

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

---

ACCESSORIES

---



**12170**

Black mattress for right bariatric wing for self-loading stretchers



**12172**

Black mattress for left bariatric wing for self-loading stretchers

---

STANDARD ACCESSORIES

---



**12080-001**

Monting kit for bariatric platform on self loading stretcher

ACCESSORY OF



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with variable height and four swivel castors



**Mercury 7070 PROOF**  
Certified stainless steel variable height stretcher



**Mercury 7070/4RG PROOF**  
Stainless steel variable height stretcher certified with 4 rotating wheels



**Mercury 7072 PROOF**  
Certified stainless steel variable height stretcher high version



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite Cinque 7095/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Lite Cinque 7097/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher - high version