



## 1740

---

Rotating wall aneroid sphygmomanometer



---

### FEATURES

Round sphygmomanometer. Easy to read, it is ideal for ambulances and ambulatories. 65° rotating in both directions, roomy cuff basket for storage of accessories, coiled tube extendible to a length of 3 m, latex bulb and chrome-plated air release valve, highly legible dial-plate diam. 147 mm. Specially tempered copper-beryllium diaphragm. Calibrated nylon Velcro® cuff with 2-tube latex bladder.

---

### IMAGES



---

## ACCESSORIES

---



**1740-001**

Pediatric cuff for rotating wall  
aneroid sphygmomanometer Art.  
1740.

---