



616

Certified vertical fastening system for scoop stretchers



FEATURES

Aluminium support with safe quick locking device for wall mount or vertical storage of a scoop stretcher. The stretcher can be folded and then locked.

With our item 616/A it is possible to fix the scoop stretcher folded.



Device certified according to the European harmonized standards UNI EN 1789

MAIN FEATURES

Length (cm)	44
Width (cm)	2,5
Height (cm)	7,5
Weight (kg)	0,5

SEE ALSO



Maxima 630
Grey anodized anti trauma scoop
stretcher



616/A
Certified vertical fixing system for
scoop stretchers



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630/V
Green painted anti trauma scoop
stretcher



Maxima 630/G
Yellow painted anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 630/A
Orange painted anti trauma scoop
stretcher