



# 630/V Maxima

---

Green painted anti trauma scoop stretcher



---

## FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
Type	
10G Certification	<input checked="" type="checkbox"/>
Material	Painted aluminium
Trauma	<input checked="" type="checkbox"/>
Alpine	<input checked="" type="checkbox"/>

SEE ALSO



**Maxima 630/N**  
Black painted anti trauma scoop  
stretcher



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher



**Maxima 630/G**  
Yellow painted anti trauma scoop  
stretcher



**Maxima 630/B**  
Blue painted anti trauma scoop  
stretcher



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher

---

ACCESSORIES

---



**630/CS**

Carrying case for scoop stretchers



**616/A**

Certified vertical fixing system for  
scoop stretchers



**616**

Certified vertical fastening system  
for scoop stretchers



**614**

Certified horizontal fastening  
system for scoop stretchers



**603/MEB**

Nylon bag containing 3 belts Art.  
607/Arneb