




644 Fixed stretcher

Stiff fixed stretcher for storage



FEATURES

Fixed rigid stretcher for storage, easily piled on top of same-type stretchers for storage optimization. Supplied with 2 quick-release straps for patient immobilization (art. 604/MEB).

- Light and handy
- Canvas in orange nylon, fire proof, impermeable, tear-resistant and easy to clean
- Sturdy aluminium frame
- Non-foldable

MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Height (cm)	17
Width (cm)	48
Max length (cm)	186
Weight (kg)	4,9
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Rigid stretcher 644/S
Stiff fixed stretcher for storage with
backrest

ACCESSORIES



604/MEB
Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon

STANDARD ACCESSORIES



604/MEB
Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon
