



## 7072 PROOF Mercury

Certified stainless steel variable height stretcher high version



### FEATURES

Designed to take different configurations, it is provided with a multi level system using 2 sliding pistons on each axle. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "Fixing Systems"). Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

#### Features:

- Ergonomic controls (for rescue vehicle loading operations and for configuration adjustments of the patient support surface)
- Adjustable backrest with lockable gas spring for inclinations from 0° to 75°
- Fully folding 600 mm sides
- HDPE buffers for loading onto the rescue vehicle (on both front and back legs)
- 2 fixed front wheels, Ø 200 mm with parking brakes
- 2 rear castors, 360° Ø 200 mm self centring
- 4 loading wheels of which 2 are Ø 100 mm and 2 are Ø 125 mm (enhanced loading system)
- Equipped with mattress Art. 7013 and belts (2 Art. 608-A/MEB o 608-N/MEB and 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB)

This version has the same features of the standard model with the only difference that it is possible to load the stretcher on ambulances with higher loading height.




Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
10G Certification	
Multilevel	
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Type	Self loading stretcher
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	43
Loading capacity (kg)	250
Material	Stainless steel ABS HDPE (High-density polyethylene)

---

## IMAGES

---



---

SEE ALSO

---



**Mercury 7070/4RG PROOF**  
Stainless steel variable height  
stretcher certified with 4 rotating  
wheels



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable  
height stretcher high version with  
4 rotating wheels



**Mercury 7070 PROOF**  
Certified stainless steel variable  
height stretcher

---

## ACCESSORIES



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**7014/A**  
4 pcs. Skay fire-proof electrowelded mattress



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12018N**  
Telescopic foldable phlebotoclysis stand for self loading stretchers



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**12070/BG**  
"COVER OX 2-3 BG" - professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG" - professional cylinder holder for 5 lt oxygen bottles



**12072/BR**  
"COVER OX 5 BR" - professional cylinder holder for 5 lt oxygen bottles



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**7013/A**  
3 pcs. Skay fire-proof electrowelded mattress



**7013/N**  
3 pcs. Skay fire-proof electrowelded mattress



**12054**  
Oximeb Maxi - Oxygen cylinder holder



**707**  
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



**12084-001**  
Cover for bariatric wings set



**12170**  
Black mattress for right bariatric wing for self-loading stretchers



**12172**  
Black mattress for left bariatric wing for self-loading stretchers



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**12025**

Black multi usage object holder for automatic loading stretcher



**12200**

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

---

## STANDARD ACCESSORIES

---



**12025**

Black multi usage object holder for automatic loading stretcher



**7013/A**

3 pcs. Skay fire-proof electrowelded mattress



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)