



7240 PROOF Frog

“FROG PLUS” self loading stretcher with trendelenburg certified



FEATURES

Made with electro polished STAINLESS STEEL. Supplied with 4 wheels Ø 200 mm (2 fixed with brake and 2 rotating self centring). The inclination of the legs is designed to give the stretcher better stability during use. The trendelenburg position is operated by a gas piston and can be locked in any position from 0° to 20°.

- Rigid platform made of reinforced thermoplastic material.
- Reclining backrest 0°- 80°, operated by a gas piston, can be locked in any position.
- Folding sides.
- Vertical legs protected by plastic buffers.
- Casters are self-centring so that the stretcher can be positioned in the ambulance without them having to be turned manually.
- Front trolley can be lowered to reduce overall dimensions of the stretcher.

In coupling with the proper fixation Art. 955 Proof, it is 10G resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB)



Device certified according to the European harmonized standards UNI EN1865-UNI EN1789

MAIN FEATURES

Usage	Transport
10G Certification	
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Type	
Height (cm)	92
Max loading height (cm)	66
Folded height (cm)	32,5
Min loading height (cm)	60
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	166
Weight (kg)	44
Loading capacity (kg)	250
Material	Stainless steel
Filter feature	

SEE ALSO



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium
stretcher with trendelenburg
certified



Frog 7210 PROOF
"FROG" yellow self loading stretcher
with trendelenburg certified



Frog 7240/INC PROOF
"FROG PLUS INC PROOF" self
loading trolley for incubators
predisposed for certified fixing
system



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg

ACCESSORIES



936

Phlebotomy stand for automatic loading stretcher Winner and Frog



939

Instrument support for stretchers



938

Oximeb oxygen cylinder holder



12080

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



955 PROOF

EN 1789 fixing system for Proof roll-in stretchers



12080-001

Monting kit for bariatric platform on self loading stretcher



12025

Black multi usage object holder for automatic loading stretcher



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



7013/A

3 pcs. Skyak fire-proof electrowelded mattress



7013/N

3 pcs. Skyak fire-proof electrowelded mattress



12084

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



12154

"BEAR PODS" Bariatric containment system for self-loading stretchers



12106

Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



12018N

Telescopic foldable phlebotomy



12070/BG

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles

stand for self loading stretchers



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR
"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles

STANDARD ACCESSORIES



7013/A
3 pcs. Skay fire-proof electrowelded
mattress



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)