



8020/N

Black flat phonendoscope



FEATURES

Flat phonendoscope with anodized stamped and lathed aluminium head diameter 46 mm. Threaded metal ring for membrane blocking. Metal chromium plated binaural, length 150 mm. Y rubber tube in one piece, diameter 5 mm. length 60 cm. Single printed carton box. Available in 3 different colors.

MAIN FEATURES

Size	Adults
Cardio	<input checked="" type="checkbox"/>

SEE ALSO



8025/N

Black stethophonendoscope



8025/G

Yellow stethophonendoscope



8025/R

Red stethophonendoscope



8020/R

Red flat phonendoscope



8020/G

Yellow flat phonendoscope



8030/B

Blue stethophonendoscope
"Rappaport"

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE