



8025/R

Red stethophonendoscope



FEATURES

Stethophonendoscope with head diameter 47/31 mm. Aluminium anodized head, chromium plated metal binaural, Y rubber tube in one piece. Rotating piper-fitting for both functions. Available in 3 different colors.

MAIN FEATURES

Size	Adults
Cardio	<input checked="" type="checkbox"/>

SEE ALSO



8020/R
Red flat phonendoscope



8030/B
Blue stethophonendoscope
"Rappaport"



8025/G
Yellow stethophonendoscope



8025/N
Black stethophonendoscope



8025/R
Red stethophonendoscope



8020/G
Yellow flat phonendoscope