



870 Sirius

Rigid cervical collar in 1 piece size TALL



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	

SEE ALSO



Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Sirius 876

Cervical collar set 6 pieces with carrying case



Antares 856

Cervical collar set 6 pieces with carrying case



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848