



894 Snake


Indipendent chambers vacuum mattress



FEATURES

Vacuum mattress with 14 indipendent chambers for an excellent and fast immobilisation. Extremely light and sturdy, 8 ergonomics handles, restraining belts and a perfect patient containment, give to the mattress unique solutions. Usable between - 30°C and + 70°C. An unmatched vacuum durability without any deflection, **the UNI EN 1865 certificate**, according to the european harmonized standards, **complete the features of this device.**

MAIN FEATURES

| | |
|-----------------------|---|
| Size | Adults |
| Usage | Immobilization |
| Length (cm) | 210 |
| Width (cm) | 100 |
| Weight (kg) | 8 |
| Trauma |  |
| Material | Polyurethane |
| Loading capacity (kg) | 150 |

SEE ALSO



Chesterfield 9114

Vacuum mattress with
independent chambers with
central welding

ACCESSORIES



888

Aluminium aspiration pump



890-001

Carrying case for vacuum mattress

STANDARD ACCESSORIES



894-001

Aspiration valve for vacuum splints
and mattress
