



914 PROOF Winner

Certified self loading stretcher with Trendelenburg and Fowler



FEATURES

Art. 914 PROOF is similar to model 914, the main difference is given by the mean of a rigid steel frame which allows, in coupling with the appropriate fixing system, to absorb those stresses caused by sudden accelerations or decelerations of the rescue vehicle. Complete with mattress and belts.

This article has been designed to offer a practical and safe instrument to the operator in case of emergency. The stretcher, built with anodized aluminium, is strong and handy thanks to its four wheels diam. 150, two fixed with brakes and two rotating. Important feature of this article is the loading height regulation system. This stretcher in fact is supplied with two sturdy perforated plates that allow to change the little front wheels height in order to load the stretcher on ambulances of various model and with different inside fittings.

It also has the possibility of having Trendelenburg and Fowler positions becoming this way a more complete and functional support for every emergency situation. Back rest and lying part made with thermoplastic material reinforced to allow cardiac massage.


Easy to clean and to disinfect. Adjustable back rest 0° - 80°, moved by a gas piston that allow to fix the back rest safely in any position. Trendelenburg and Fowler positions moved by gas pistons. Side protections completely overturnable. Blocking device for folded legs with a system that avoid the compression of the blocking bar. Safety lever for the legs positioned near the unlocking device allowing thus the release operation for the loading, keeping the hands on the stretcher. Vertical legs protected by nylon wedges. Self centring device mounted on rotating wheels.

This system automatically blocks the back wheels in the central position during the loading of the stretcher on the ambulance avoiding to turn the wheels manually.



Device certified according to the European harmonized safety standards UNI EN1865 - UNI EN 1789.

MAIN FEATURES

Usage	Transport
Lying surface	With Trendelenburg and Fowler
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	Self loading stretcher
10G Certification	
Height (cm)	87,5
Width (cm)	56
Weight (kg)	40
Loading capacity (kg)	170
Folded height (cm)	26
Max loading height (cm)	65
Min loading height (cm)	62
Length (cm)	194
Length with folded trolley (cm)	155
Material	Anodized aluminium

SEE ALSO



Winner 910 PROOF
Self loading stretcher - certified

ACCESSORIES



7014/A
4 pcs. Skay fire-proof electrowelded
mattress



955 PROOF
EN 1789 fixing system for Proof roll-
in stretchers



936
Phleboclysis stand for automatic
loading stretcher Winner and Frog



939
Instrument support for stretchers



938
Oximeb oxygen cylinder holder



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)



12025
Black multi usage object holder for
automatic loading stretcher



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



7013/A
3 pcs. Skay fire-proof electrowelded
mattress



7013/N
3 pcs. Skay fire-proof electrowelded
mattress



12018N
Telescopic foldable phleboclysis
stand for self loading stretchers



12070/BG
"COVER OX 2-3 BG"- professional
cylinder holder for 2 and 3 lt oxygen
bottles



12070/BR
COVER OX2-3 BR professional
cylinder holder with belt and
orange MeBer logo for 2 and 3 lt
oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



12054

Oximeb Maxi - Oxygen cylinder holder

STANDARD ACCESSORIES



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



7013/A

3 pcs. Skay fire-proof electrowelded mattress