



960

"Extra fixing system" - fixing system for chair Extra



FEATURES

Fixing device designed to ensure a permanent fixation of the Extra chair inside the ambulance patient compartment. It is made of corrosion resistant and easy to clean materials (stainless and thermoplastic). Extra fixing system enables to lock the whole Extra range in a simple, fast and effective way. The operator can easily lock and unlock the chair, thanks to the innovative and exclusive monobloc configuration.

Device certified according to the European harmonized standards UNI EN 1789.



MAIN FEATURES

Width (cm)	23
Height (cm)	42
Depth (cm.)	20
Weight (kg)	5,8

IMAGES



SEE ALSO



975
"Ergolift fixing system" - Fixing
system for chair Extra Ergolift

ACCESSORY OF



Extra 670/BR
Track assisted stair chair with arms



679
Narrow foldable sedan chair



Extra Lite Standard 676
Foldable carrying chair with 4
wheels, footrest and arms



678 Extra Lite Basic
Extra Lite Basic, foldable carrying
chair