



## 15244

Folding Scangliss



### FEATURES

Foldable transfer board to use under the patient allowing a safe transfer from stretcher to bed. It is translucent, a-magnetic and it is divided in two sections. This device consists of a rigid polyethylene core, padded with an amorphous material of 1,5 cm thickness, all wrapped in a sliding nylon sheet. Scan Gliss is CE certified according to Directive 93/42 EEC Medical Devices.

**Maximum load allowed: 180 kg.**

### MAIN FEATURES

Length (cm)	180
Width (cm)	50
Height (cm)	2
Weight (kg)	4
Loading capacity (kg)	180
Material	Polyethylene

### SEE ALSO



**1398**

White/blue absorbing and transfer sheet cm. 220x100



**15180**

Big Gurney

---

ACCESSORIES

---



**15244-001**  
Nylon cover for Scangliss

---