

Product catalogue 2023

**MEBER<sup>+</sup>**



Products, technologies and services  
for first aid and medical emergency


[www.meber.it](http://www.meber.it)





# INDEX

Category	Page
Stretchers	14
Ambulance Cots and Trolleys	34
Carry and Stair Chairs	71
Accessories	90
Immobilization	110
Diagnostics	136
Resuscitation	150
First Aid	162
Fastening Systems	186
Restraint Systems	192
Supports	196
Veterinary	203



**Discover MeBer products on-line**

Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website. You might be required to install a free application to navigate.

Alternatively, go to:

<https://www.meber.it/en/prodotti.aspx>



# MEBER

PRODUCTS,  
TECHNOLOGY AND  
SERVICE FOR **FIRST  
RESPONSE** AND  
**EMERGENCY AID**

## AT YOUR SIDE IN EVERY STEP OF THE WAY

---

Founded in 1987 in Parma, MeBer has grown to become a leading player in the field of design and production of medical equipment for emergency aid, rescue and first response. Ever since its origin MeBer stood out for its company culture, based on passion, professionalism, quality and innovation: essential characteristics to become a loyal partner for all operators.

**MEBER<sup>+</sup>**  
IN CASE OF EMERGENCY





IN CASE OF EMERGENCY



## MISSION AND CORPORATE CULTURE

6

In MeBer, planning and design are an out-and-out mission for all staff members: helping and supporting the rescue operators in order to achieve excellent performances, in complete safety, is our true mission.

## MEBER A RELIABLE PARTNER, WORLDWIDE

MeBer counts on a vast sales network on the domestic sector as well as a widespread distribution network made up of distributors and re-sellers in more than 60 countries.

MeBer is an ever-present member of the most important Trade Fairs such as REAS (Montichiari-Italy), Medica (Düsseldorf – Germany), Arab Health (Dubai – Arab Emirates), Rettmobil (Fulda – Germany), FIME (Miami – USA) and actively participates in a variety of local events and congresses.





## INNOVATION, TECHNOLOGY AND CUSTOMER CARE

MeBer is a young and dynamic company, strongly oriented towards the needs of the market and a constant improvement of the quality of its products. Thanks to the efficiency of a qualified team and by means of a dedicated communication with emergency operators, MeBer invests a major portion of its resources in the research and development of new advanced technological devices, certified according to the most recent European Standards (UNI EN 1865 – UNI EN 1789). The research of innovative materials as well as the application of state-of-the-art mechanisms have characterized MeBer's range over the years.

## TOTAL QUALITY MANAGEMENT

MeBer has a whole department dedicated to the implementation of the quality management system of its products, following every phase of their life-cycle, starting from their design up to the after-sales service.

This led the company to achieve the European Standards UNI EN ISO 9001:2015 and UNI CEI EN ISO 13485:2016.

Every step of the production, from the design phase up to the development and manufacturing, are carried out internally with the aid of state-of-the-art machines and equipment as every phase undergoes strict quality monitoring.

The rigid quality process continues up to the after sales service, namely customer service and technical assistance (ordinary and extra-ordinary) that the company oversees to optimize performance and reliability of the products.

# THE BEST PARTNER...

**PRODUCTS,  
TECHNOLOGIES AND  
SERVICES FOR FIRST  
AID AND MEDICAL  
EMERGENCY**



8

## MEBER ON-SITE (Italy) Directly at your premises

One of our technicians will visit your premises and carry out all necessary activities which might include:

- ✓ Ordinary maintenance
- ✓ Extra-ordinary maintenance
- ✓ Restoration
- ✓ Updates (if applicable)

## MEBER SERVICE NETWORK (Italy and World) Through authorized Service Centres

The products are shipped directly to Me.Ber or to an authorized Service Centre of your choice. We carry out all necessary tests and maintenance. The devices are then shipped back to you in the shortest possible time.

# DEDICATED TO MEDICAL EMERGENCY OPERATORS



**MEBER**, THE ONLY PARTNER THAT  
CARES ABOUT YOUR DEVICES AND  
KEEPS THEM EFFICIENT

A professional after sales service carried out by the manufacturer is a crucial support for any structure: it is a requirement in most Italian regions if the brand wishes to obtain accreditation.

Furthermore, the law states that it is compulsory for every structure to ensure the correct efficiency of their own devices. Proper maintenance is a fundamental aspect of the medical equipment sector, as it ensures the safety of rescue operators and patients. It also constitutes legal protection to the Coordinators and Administrators in the event of an accident.



# MEBER CUSTOMER CARE

**MeBer's customer care comprises all of the exclusive after-sales services offered by the company.**

The service is available either on-site (i.e. at the customer's premises), at MeBer's headquarters or at one of the authorized service centres, depending on the type of maintenance required and on the geographical area.







# QUALITY

## 10G PROOF PRODUCTS



### MEBER: SHOCK-PROOF EMERGENCY!

MeBer has put an emphasis on the safety and reliability of the products ever since the beginning of its journey, pioneering the birth of European standards and regulations as well as anticipating the needs of the emergency sector as we know it. The philosophy behind the company's path so far is driven by the hard-rooted conviction that a safe and regulated medical device shall safeguard both patients and operators, reducing the risk of accidents and consequently protecting the operators from moral and legal responsibilities. MeBer has performed 10G tests since 2001: today, the vast majority of the products in the range are certified!

### EUROPEAN REGULATIONS FOR A SAFE TRANSPORT

The European Standards that regulate the transport of stretchers and rescue devices on the ambulances are EN 1865 and EN 1789, which integrates and completes the former. European Standard EN 1865, known as UNI EN 1865 in Italy, has been in force since 1999 and dictates rigid technical standards, configurations, maximum height and length, maximum weight, minimum load, etc. to ensure a safe transport of the medical devices during transport.

UNI EN 1789 completed the regulations, imposing the requirement of 10G PROOF testing of the medical devices and relative fastening and locking systems: the devices transported by an ambulance must endure, in the event of an accident, an acceleration equal to 10 times the weight of the device. Patient's transport devices must also ensure that they can hold a 75 KG individual (average body weight).

Specifically, the regulation requires that the afore-mentioned force is applied in all directions and that no damage is sustained by either patient or operators on board of the vehicle.

# WHAT'S NEW ON THIS CATALOGUE



## NEW PRODUCTS AND UPDATES

<p>p. 44</p> <p>Art. 7075/4RG PROOF 7077/4RG PROOF</p>	<p>p. 46</p> <p>Art. 7095/4RG PROOF 7097/4RG PROOF</p>	<p>p. 84</p> <p>Art. 668</p>	<p>p. 96</p> <p>Art. 12084 - 12086 - 12088 12170 - 12172</p>	<p>p. 97</p> <p>Art. 12174 - 12176</p>	<p>p. 103</p> <p>Art. 707</p>
<p>p. 105</p> <p>Art. 12084-001</p>	<p>p. 142</p> <p>Art. 8377</p>	<p>p. 143</p> <p>Art. 8540</p>	<p>p. 156</p> <p>Art. 2411 - 2413 - 2414 - 2415</p>	<p>p. 158</p> <p>Art. 1233 - 1234</p>	<p>p. 164</p> <p>Art. 1319</p>
<p>p. 164</p> <p>Art. 1318</p>	<p>p. 166</p> <p>Art. 1322</p>	<p>p. 177</p> <p>Art. 11162</p>	<p>p. 181</p> <p>Art. 1982</p>	<p>p. 187</p> <p>Art. 977</p>	<p>p. 195</p> <p>Art. 600/N</p>
<p>p. 195</p> <p>Art. 601/N</p>	<p>p. 200</p> <p>Art. 19170 - 19172</p>	<p>p. 202</p> <p>Art. 19115</p>	<p>p. 202</p> <p>Art. 19123</p>	<p>p. 202</p> <p>Art. 787L</p>	

## AN OVERVIEW OF THE CATALOGUE...

### APPLICATION SYMBOLS

To understand if a product suits  
your needs, at a glance!



Pediatric

**NEONATAL AND/OR PEDIATRIC  
DEVICE**



Bariatric

**PRODUCT FOR BARIATRIC  
PATIENTS**



Cardio

**DEVICE FOR THE TREATMENT OF  
CARDIOLOGICAL EMERGENCIES**



Trauma


**DEVICE FOR THE TREATMENT OF  
TRAUMATIC EMERGENCIES**



Alpine

**PRODUCT FIT FOR ALPINE  
RESCUE**

The new products are  
**easily spotted!** Search for  
the "New Product" Logo  
to easily identify what's  
new!



# STRETCHERS

Stretchers

14







A stretcher is often the first rescue and transfer device of a patient. Ranging from the traditional and ultra-tested scoop stretchers in the 630 SERIES to the revolutionary ERGON, Me.Ber's stretchers are the archetype of reliability at an international level. The constant exchange of information with the rescue operators is the key of the

success of these products. Thanks to the choice of materials, their simplicity of use and ergonomics, they represent a benchmark in the medical equipment field. Me.Ber offers a very rich line of stretchers for any operating environment, ranging from adults to pediatric and bariatric patients.

Discover MeBer products on-line. Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website. You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/stretchers.aspx>



# STRETCHERS

Article Code	Type
16150/V	ERGON variable geometry comfort stretcher
16150/G	ERGON variable geometry comfort stretcher
16160/V	ERGON Junior variable geometry comfort stretcher
16160/G	ERGON Junior variable geometry comfort stretcher
631	MAXIMA grey anodized anti trauma scoop stretcher with head/trunk in 1 pc
631/G	MAXIMA yellow painted anti trauma scoop stretcher with head/trunk in one piece
630	MAXIMA grey painted anti trauma scoop stretcher
630/G	MAXIMA yellow painted anti trauma scoop stretcher
630/V	MAXIMA green painted anti trauma scoop stretcher
630/B	MAXIMA blue painted anti trauma scoop stretcher
630/A	MAXIMA orange painted anti trauma scoop stretcher
630/N	MAXIMA black painted anti trauma scoop stretcher
6200	Pediatric anti trauma scoop stretcher
638	Folding rescue stretcher with adjustable back rest
640	Folding rescue stretcher
639	Fixed rescue stretcher with comfort position and backrest
643	Duo fold rescue stretcher
646	Folding rescue stretcher
644	Stiff fixed stretcher for stowage
644/S	Stiff fixed stretcher for stowage with back rest
660	4 pieces folding stretcher with telescopic handles
662	2 pieces folding stretcher with telescopic handles
1550	Military folding stretcher with green canvas
1510 PROOF	Stiff stretcher with nylon cover
1515 PROOF	Stiff stretcher with rigid cover
13010	ICARUS multi level stretcher with side protection
13012	ICARUS multi level stretcher with stainless steel tank
13018	CRONOS multi level stretcher with back rest and Trendelenburg
16100	TOBOGA stretcher with belts and foot support
16102	TOBOGA detachable stretcher with belts and foot support

Materials:  Stainless steel  Anodized aluminium  Painted Aluminium  ABS

## COMPARATION CHART

Colour	10G Certification	Material	Height	Width	Min length	Max length	Folded length	Thickness folded	Weight	Load
	✓	HC <sup>3</sup>	7,5 cm	43 cm	165 cm	201 cm	120 cm	9 cm	10,5 kg	170 kg
	✓	HC <sup>3</sup>	7,5 cm	43 cm	165 cm	201 cm	120 cm	9 cm	10,5 kg	170 kg
		HC <sup>3</sup>	7,5 cm	43 cm	129 cm	153 cm	85 cm	9 cm	8,4 kg	170 kg
		HC <sup>3</sup>	7,5 cm	43 cm	129 cm	153 cm	85 cm	9 cm	8,4 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	9,5 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	9,5 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
	✓		7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
			7 cm	42 cm	130 cm	-	-	-	7,8 kg	170 kg
			14 cm	49 cm	187 cm	-	98 cm	17 cm	9,2 kg	170kg
			14 cm	48 cm	186 cm	-	92 cm	17 cm	6,5 kg	170 kg
			14 cm	49 cm	189 cm	-	-	-	12 kg	170 kg
			12 cm	49/12 cm	195 cm	-	97 cm	17 cm	5 kg	170 kg
			14 cm	50 cm	201 cm	-	101 cm	17 cm	6,2 kg	170 kg
			17 cm	48 cm	186 cm	-	-	-	4,9 kg	170 kg
			17 cm	48 cm	186 cm	-	-	-	5,9 kg	170 kg
			18 cm	57/11 cm	192 cm	216 cm	94 cm	21 cm	9,3 kg	170 kg
			18 cm	57/11 cm	192 cm	216 cm	192 cm	-	9,3 kg	170 kg
			17 cm	58/11 cm	198 cm	231 cm	198 cm	-	9 kg	170 kg
	✓		16 cm	57 cm	196 cm	233 cm	-	-	15 kg	170 kg
	✓		16 cm	57 cm	196 cm	233 cm	-	-	20,8 kg	170 kg
			42/88 cm	62 cm	191 cm	-	-	-	29,5 kg	170 kg
			51/96 cm	62 cm	198 cm	-	-	-	39 kg	170 kg
			48/86 cm	59 cm	198 cm	-	-	-	37 kg	170 kg
			17 cm	61 cm	216 cm	-	-	-	15 kg	300 kg
			17 cm	61 cm	216 cm	-	108/114 cm	-	15 kg	300 kg

# ERGON

VARIABLE GEOMETRY COMFORT STRETCHER



Art.  
16150/G



Alpine



Pediatric



Trauma

Art.  
16150/V



Art. 16150/V / Art. 16150/G



## ERGON

### Variable geometry comfort stretcher

ERGON was designed by rescue operators and experts of the medical emergency sector to offer outstanding performances, never achieved before. This stretcher ensures maximum comfort, especially in bed sore areas, thanks to a special material (the exclusive MeBer HC3 – High Comfort Composite Cell) consisting in a soft material moulded on a high resistance plastic polymer.



## Technical data

Art.	Colour	Material	Height	Width	Min length	Max length	Folded length	Thickness folded	Weight	Load
16150/V		HC <sup>3</sup>	7,5 cm	43 cm	165 cm	201 cm	120 cm	9 cm	10,5 kg	170 kg
16150/G		HC <sup>3</sup>	7,5 cm	43 cm	165 cm	201 cm	120 cm	9 cm	10,5 kg	170 kg


# # flexilock

## OPERATING PROCEDURES



A technological revolution: the new jointed locking system Flexilock allows to separate the stretcher in two halves and facilitates rescue operations even in extreme conditions.

Beside fitting all standard rescue techniques, ERGON allows to easily perform rescue operations in extreme environments and when the patient is in critical conditions.



## # HC<sup>3</sup> High Comfort Composite Cell

Hard plastic reinforcements were added to the corners to avoid damaging of the soft material during rescue operations or storage on the vehicles. These reinforcements also grant better sliding characteristics inside of the storage compartments. All indications on how to operate the Flexilock locking device are highlighted with a high contrast colour, in order to make them visible in dark and low-visibility situations.



Beside patient's comfort, ERGON is designed to assist the rescue operators thanks to unique solutions.

The handles are upholstered in HC<sup>3</sup> and present a 20° inclination, thus allowing a better grip. Dedicated eyelets and supports were integrated for the fastening of straps and carabiners. When split, the distance between the two halves measures 25 mm: this ensures excellent support to the spinal column, in contrast with regular scoop stretchers which present an 80 mm distance. The

shape of the blade is also larger than usual and is easily placed between the patient and the ground. The employment of HC<sup>3</sup> material, besides conferring excellent structural stability, ensures great thermal insulation. This feature prevents hypothermia and allows for a longer stay of the patient on the ERGON stretcher, thanks to its design

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

# ERGON


## ACCESSORIES



**Art.  
9150**

### **X-LOCK Ergonomic head immobilizer**

One-piece head immobilizer with integrated head strap, for a fast head immobilization directly on the stretcher. Made with laminated polyesters and bi-density polyurethane foams.



**Art.  
9156**

### **X-FIX Restraint system with straps and carabiners**

Standard restraint straps with carabiners: an innovative system designed to ease and accelerate any rescue operation, ensuring high comfort and providing stability to the spinal column.

**Art.  
16150-001**



Winch bag for transfers  
cm 124x32x16

**Art.  
Art. 16150-002**



Storage bag cm  
122x47x8

**Art.  
9150-001**



Blue carrying case for  
X-LOCK ergonomic  
head immobilizer

**Art.  
621**



X-DOUBLE padded  
adjustable head  
immobilizer. ORANGE  
and BLUE, specific for  
ERGON stretchers

**Art.  
572**



Blue strap in  
2 pieces with  
metallic buckle  
(155 cm length)

**Art.  
574**

Blue bag  
with 3 belts  
art. 572

**Art.  
959**



Wall mount  
system

Accessories

Factory-equipped  
with:

**Art.  
584**



Black bag with 3 orange straps in 2  
pieces and quick-release plastic buckle

# ERGON JUNIOR

PEDIATRIC VARIABLE  
GEOMETRY  
COMFORT  
STRETCHER



Pediatric

Art.  
16160/G

Art.  
16160/V

Art. 16160/V / Art. 16160/G



## ERGON JUNIOR

Pediatric variable geometry comfort stretcher

**ERGON share its advanced technology with the junior model. All benefits and main features of the adult stretcher are preserved in its Junior version.** The Flexilock locking system was moved from the centre to the side, in order to allow operations without invading the head-space. The large blade allows for maximum comfort even to pediatric patients.

21

### Technical data

Art.	Colour	Material	Height	Width	Min length	Max length	Folded length	Thickness folded	Weight	Load
16160/V		HC <sup>3</sup>	7,5 cm	43 cm	129 cm	153 cm	85 cm	9 cm	8,4 kg	170 kg
16160/G		HC <sup>3</sup>	7,5 cm	43 cm	129 cm	153 cm	85 cm	9 cm	8,4 kg	170 kg



# MAXIMA

ANTI TRAUMA SCOOP STRETCHER

United  
torso / head  
sections

Separate  
torso/head  
sections

Art. 630

Art. 631



Alpine



Trauma

Art. 6200



Art. 630

## MAXIMA

Grey anodized anti trauma scoop stretcher

Thanks to its features, this stretcher can be placed under the patient with ease, without lifting the body. Built with carefully designed materials, MAXIMA ensures solidity, reliability and durability.

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Art. 6200

## MAXIMA

Pediatric anti trauma scoop stretcher



Pediatric

Fit for pediatric rescue, reduced size.

Art. 631

## MAXIMA

Grey anodized anti trauma scoop stretcher with head/torso in 1 piece

Scoop stretcher: the difference between this model and Art. 630 is that the head and torso areas are united in 1 piece. They share the same size.

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

- ✓ Efficient splitting system with locking device
- ✓ Adjustable length with automatic locking device
- ✓ Bag with 3 straps (quick release) for the immobilization of the patient (Art. 603 / MEB)
- ✓ Easily arranged inside the ambulance when folded
- ✓ Built with light and sturdy aluminium
- ✓ Carrying case available on request (Art. 630/CS).



Art. 630/N

Art. 630/B



Art. 630/A



Art. 630/V



Art. 630/G



## MAXIMA

### Painted anti trauma scoop stretcher

Keeping the very same features of art. 630, thanks to a special paint (designed for optimal durability) this stretcher can be used to differentiate operators or a specific area of operation. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



Art. 631/G

## MAXIMA

### Yellow painted anti trauma scoop stretcher with head/torso in one piece

Keeping the same features of Art. 631, this stretcher helps differentiating operators or a specific area of operation thanks to its special yellow paint (designed for optimal durability). **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789**

#### Accessories

Art. 630/CS



Carrying case for scoop stretchers

Art. 624



**DOUBLE**  
Head  
immobilizer for  
scoop stretchers

Art. 623



**DOUBLE PLUS**  
Adjustable head  
immobilizer for scoop  
stretchers


#### Fastening Systems

Art. 614



Certified horizontal fastening system for scoop stretchers

Art. 616  
Art. 616/A



Certified vertical fastening system for scoop stretchers

Art.	Colour	Material	Height	Width	Min length	Max length	Folded length	Thickness folded	Weight	Load
630		Anodized aluminium	7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
631		Anodized aluminium	7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	9,5 kg	170 kg
631/G		Painted Aluminium	7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	9,5 kg	170 kg
630/G 630/V 630/B 630/A 630/N		Painted Aluminium	7 cm	42 cm	167 cm	202 cm	120 cm	9 cm	10 kg	170 kg
6200		Anodized aluminium	7 cm	42 cm	130 cm	-	-	-	7,8 kg	170 kg

# RESCUE STRETCHERS



Art. 639

## Fixed rescue stretcher

### Fixed rescue stretcher with comfort position and backrest

Rescue stretcher ideal for ambulances and rescue vehicles in general. It features a reclining backrest, operated by a gas spring, safe and lockable in any position.

This item also presents a leg raiser for comfort, making this stretcher a very versatile device.

The wheels and feet of the stretcher are turned inward, to decrease its volume.

Built in aluminium and fire-resistant, impermeable, tear-resistant PVC-coated nylon canvas.

It is supplied with a bag containing 2 straps for patient immobilization (Art. 604/MEB)..




Art. 640

## Collapsible rescue stretcher

Suitable as a secondary stretcher inside of an ambulance. No backrest.

It is supplied with a bag containing 2 straps for patient immobilization (Art. 604/MEB).



Art. 638

## Collapsible rescue stretcher

### With adjustable backrest

Excellent as a secondary stretcher inside of an ambulance. It is supplied with a bag containing 2 straps (Art. 604 / MEB).

Supplied with:

Art. 604/MEB

Orange bag with 2 straps 605 / MEB



Fastening System

Art. 964

Front fastening system for rescue stretcher



Art. 950/4

Rear fastening system



Art.	Colour and Material	Height	Width	Min length	Folded length	Thickness folded	Weight	Load
638	Anodized aluminium	14 cm	49 cm	187 cm	98 cm	17 cm	9,2 kg	170 kg
639		14 cm	49 cm	189 cm	-	-	12 kg	170 kg
640		14 cm	48 cm	186 cm	92 cm	17 cm	6,5 kg	170 kg



Art. 646

## Collapsible rescue stretcher

Foldable in 2 pieces



Designed for the transport also of heavy patients, this stretcher can be folded in two pieces. The canvas is built in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Straps are not included (see "Restraint systems" Art. 604/MEB).

Art. 643

## Collapsible rescue stretcher

Foldable in 4 pieces

Built in tubular anodized aluminium, canvas in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Ideal for storage inside of any emergency vehicle. Straps are not included (see "Restraint systems" Art. 604/MEB).

Art.	Colour	Material	Height	Width	Max length	Folded length	Thickness folded	Weight	Load
643		Anodized aluminium	12 cm	49/12 cm	195 cm	97 cm	17 cm	5 kg	170 kg
646		Anodized aluminium	14 cm	50 cm	201 cm	101 cm	17 cm	6,2 kg	170 kg



# RESCUE STRETCHERS



Art. 644

Piled  
Stretchers

Art. 644/S

Art. 644

## Stiff fixed stretcher

Suited for storage

Fixed rigid stretcher for storage, easily piled on top of same-type stretchers for storage optimization.  
Supplied with 2 quick-release straps for patient immobilization (Art. 604/MEB).


Art. 644/S

## Stiff fixed stretcher

Suited for storage, with backrest

Fixed rigid stretcher for storage, featuring a 3 position adjustable backrest. Supplied with 2 quick-release straps for patient immobilization (Art. 604/MEB)

- ✓ Light and handy
- ✓ Canvas in orange nylon, fire proof, impermeable, tear-resistant and easy to clean
- ✓ Sturdy aluminium frame
- ✓ Non-foldable

Art.	Colour	Material	Height	Width	Min length	Weight	Load
644		Anodized aluminium	17 cm	48 cm	186 cm	4,9 kg	170 kg
644/S			17 cm	48 cm	186 cm	5,9 kg	170 kg

Stretcher  
foldable in 4  
pieces



Art. 660



Art. 662



Folded  
stretcher

Art. 660

## Collapsible stretcher

Foldable in 4 pieces, telescopic handles

This 4 pieces foldable stretcher is the most complete and versatile in the range. The sturdy tubular aluminium frame makes it light and robust. The stretcher features 4 telescopic handles that allow for a larger support surface.

The canvas is in fire-resistant, impermeable, tear resistant, easily cleanable grey nylon.

Supplied with 2 quick-release straps for patient immobilization (Art. 605/Nmeb).

Art. 662

## Collapsible stretcher

Foldable in 2 pieces, telescopic handles

Foldable only width-wise, this stretcher shares the same features of Art. 660. Supplied with 2 quick-release straps for patient immobilization (Art. 605/Nmeb).

Accessories



Art. 662/CS



Carrying case for collapsible stretcher (foldable in 2 pieces) with telescopic handles.

Art. 660/CS



Carrying case for collapsible stretcher (foldable in 4 pieces) with telescopic handles.

Supplied  
with



Art. 605/Nmeb



2 single-piece strap belts with plastic buckle and MeBer band

Art.	Colour	Material	Height	Width	Min length	Max length	Folded length	Thickness folded	Weight	Load
660		Anodized aluminium	18 cm	57/11 cm	192 cm	216 cm	94 cm	21 cm	9,3 kg	170 kg
662		Anodized aluminium	18 cm	57/11 cm	192 cm	216 cm	192 cm	-	9,3 kg	170 kg

# RESCUE STRETCHERS

Art. 1550




Folded  
stretcher

Art. 1550

## Military Stretcher

Collapsible stretcher with green canvas

Width-collapsible stretcher built in tubular aluminium and green fireproof military canvas. It features telescopic handles and strap belts for patient immobilization. Ideal for military applications.

Art.	Colour	Material	Height	Width	Min length	Max length	Folded length	Weight	Load
1550		Anodized aluminium	17 cm	58/11 cm	198 cm	231 cm	198 cm	9 kg	170 kg
1510 PROOF		Anodized aluminium	16 cm	57 cm	196 cm	233 cm	-	15 kg	170 kg
1515 PROOF		Anodized aluminium	16 cm	57 cm	196 cm	233 cm	-	20,8 kg	170 kg



# FIXED STRETCHERS

FOR AMBULANCE  
TROLLEYS

**Art.  
7060/4RG +  
1510 PROOF**

Supplied with straps only in combination with a trolley: 2 pieces of art. 608-A/MEB or 608-N/MEB, 1 piece of art. 612/AN-2/MEB or 612/N-2/MEB



**Art.  
7060/4RG +  
1515 PROOF**




**Art. 1510 PROOF**

## Fixed stretcher with Nylon cover

Fixed stretcher built in anodized aluminium, canvas in fire-resistant and easily cleanable nylon. This model is suitable for our ambulance trolleys (TITANUS) or as a hand-carried stretcher.

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

- ✓ 4 extractable and adjustable handles
- ✓ 4 wheels (100 mm Ø)
- ✓ Adjustable backrest
- ✓ Supplied with straps (Art. 604/MEB)
- ✓ Foldable side rails



**Art. 1515 PROOF**

## Fixed stretcher with rigid bed

Fixed stretcher built in anodized aluminium with painted rigid backrest and bed. Resistant and easy to clean. Trendelenburg and Fowler positions. Suitable for our ambulance trolleys (TITANUS) or as a hand-carried stretcher. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

- ✓ 4 extractable and adjustable handles
- ✓ 4 wheels (100 mm Ø)
- ✓ Adjustable backrest
- ✓ Supplied with: 1 piece of Art. 7014, 2 pieces of Art. 608-A/MEB or 608-N/MEB, 3 pieces of Art. 605C, 1 piece of Art. 612/AN-2/MEB or 612/N-2/MEB
- ✓ Foldable side rails

# ICARUS

MULTI-LEVEL STRETCHER  
WITH SIDE RAILS



Art. 13010



Art. 13012

## ICARUS

Multi-level stretcher with  
steel tank



Icarus +  
mattress  
(1880-001)  
and cushion  
(1880-002)



Reclined  
Icarus

### Accessories

Art. 607/A-MEB



Orange straps in 2  
pieces with plastic  
buckle

Art. 607/B



Blue straps in 2  
pieces with plastic  
buckle

Art. 1880-001



Skai mattress for  
ICARUS

Art. 1880-002



Skai Cushion for  
ICARUS

Art. F1813



Cover for  
ICARUS

Art. 13010



## ICARUS

### Multi-level stretcher with side rails

Built in stainless steel, our ICARUS multi-level stretcher features 4 wheels (Ø 125 mm, two with brakes), side rails and foot rest. 4 gas springs ensure precise adjustment to any required height, for a professional and practical operation.

Available with ABS loading bed for optimal cleanliness and disinfection (Art. 13010) and steel tank (Art. 13012). Both mattress and cushion are made of a fire-resistant, waterproof and easily cleanable material.

### Technical data

Art.	Colour	Material	Height	Width	Length	Weight	Load
13010		Stainless steel	42/88 cm	62 cm	191 cm	29,5 kg	170 kg
13012		Stainless steel	51/96 cm	62 cm	198 cm	39 kg	170 kg

# CRONOS

MULTI-LEVEL STRETCHER WITH BACK REST AND TRENDELENBURG



Art. 13018



Cronos set in low position



Cronos from the side

## Accessories

Art. 602/MEB

Orange bag with 3 straps (plastic buckle)



Art. 7013/A



Art. 7013/N




Art. 13018

## Technical data

## CRONOS

### Multi-level stretcher with back rest and trendelenburg

CRONOS is built in electropolished stainless steel, with foldable bed and side rails. The backrest and Trendelenburg are easily adjustable by means of gas springs, enabling virtually infinite positions. The 4 gas springs positioned below the loading bed allow perfect alignment at any height. Lastly, 4 wheels (Ø 200mm), 2 rotating wheels with brakes and 2 fixed swivel wheels, making CRONOS a handy and versatile stretcher.

Art.	Colour	Material	Height	Width	Length	Weight	Load
13018		Stainless steel	48/86 cm	59 cm	198 cm	37 kg	170 kg



# TOBOGA

TOBOGA STRETCHERS WITH STRAPS  
AND FOOT REST



Alpine



Trauma

Art. 16100



Stretchers

## Accessories



Art. 16100-005

Certified lifting sling for  
Toboga stretchers



Art. 16100-002

Flotation system for  
Toboga stretcher



32

Art. 16100

## TOBOGA

**Stretcher with straps and foot rest.**

Consisting of an ABS modified polymer in the outer shell, supported by a strong aluminium frame. 12 hand grips ensure a very easy handling of the stretcher. Chemical-, UV-, rust- and corrosion resistant. Robust brass grommets, positioned on the sides, enable lifting the stretcher by means of art. 16100-005. One layer of EVA foam ensures high comfort to the patient. Supplied with 4 straps and adjustable foot rest.

## Technical data

Art.	Colour	Material	Height	Width	Length	Folded length	Weight	Load
16100		Anodized aluminium and ABS	17 cm	61 cm	216 cm	-	15 kg	300 kg
16102		Anodized aluminium and ABS	17 cm	61 cm	216 cm	108/114 cm	15 kg	300 kg



Alpine



Trauma

Art. 16102



Accessori



Art. 16100-005

Certified lifting sling for  
Toboga stretchers



Art. 16102-001

Carry bag for detachable  
Toboga stretcher  
Size: 121 x 64 x 24 cm



Art. 16102

## TOBOGA

**Detachable stretcher with straps and foot rest**

This stretcher shares the same features as art. 16100.  
This model is made up of two halves that enable easy storage  
and carrying. The two halves are kept together by strong alumi-  
nium rods locked together with steel pivots.



# AMBULANCE COTS & TROL



Ambulance cots and trolleys



# TROLLEYS


Ambulance cots and trolleys are the emblem of emergency rescue and first aid. They play the most delicate role on an emergency vehicle and must ensure a safe transport in order to avoid the worsening of a patient's condition. Furthermore, they must also safeguard the rescuers. Many precautions must be taken into consideration to ensure that the rescuers stay injury-free:

ergonomics, weight reduction and many other measures that prevent back injuries. Besides, all "10G PROOF" models are certified to ensure resistance to 10G acceleration forces, which are applied to the device during its testing phase to prove that it doesn't deform nor detach in case of a road crash.



Discover MeBer products on-line  
Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website.  
You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/self-loading-stretchers.aspx>



# AMBULANCE COTS & TROL

	Item code	Name	Type	Colour	10G certificate	Material	Adjustable height	Trendelenburg	Fowler	Intermediate height
Self Loading Trolleys	7060 PROOF	TITANUS	Low version	■	✓	◆	✓			61 cm
	7062 PROOF		High version	■	✓	◆	✓			61 cm
	7060/4RG PROOF	TITANUS 4	4 rotating wheels, low version	■	✓	◆	✓			61 cm
	7062/4RG PROOF		4 rotating wheels, high version	■	✓	◆	✓			61 cm
	7240/INC	FROG INC	Plus, for incubators	■		◆				
	7240/INC PROOF		CERTIFIED plus for incubators	■	✓	◆				
	7260/INC PROOF		CERTIFIED Lite for incubators	□	✓	◆				
	910/GIN	WINNER GIN	For incubators	■		◆				
	910/GIN PROOF		CERTIFIED for incubators	■	✓	◆				
Self Loading Stretchers	7034		Fixed trolley for incubator transport	■		◆				
	7112 PROOF	EVOX	CERTIFIED Chair/Cot	■	✓	◆				
	7080/4RG PROOF	MERCURY LITE 4	4 rotating wheels, low version	■	✓	◆	✓	✓		72 cm
	7082/4RG PROOF		4 rotating wheels, high version	■	✓	◆	✓	✓		75 cm
	7070 PROOF	MERCURY	Low version	■	✓	◆	✓	✓		72 cm
	7072 PROOF		High version	■	✓	◆	✓	✓		75 cm
	7070/4RG PROOF	MERCURY 4	4 rotating wheels, low version	■	✓	◆	✓	✓		72 cm
	7072/4RG PROOF		4 rotating wheels, high version	■	✓	◆	✓	✓		75 cm
	7075/4RG PROOF	MERCURY 5	4 rotating wheels, low version	■	✓	◆	✓	✓		95 - 87 - 79 - 62 - 39 cm
	7077/4RG PROOF		4 rotating wheels, high version	■	✓	◆	✓	✓		100 - 90 - 81 - 65 - 39 cm
	7095/4RG PROOF	MERCURY 5 LITE	4 rotating wheels, low version	■	✓	◆				94,5 - 87 - 79 - 68,5 - 39 cm
	7097/4RG PROOF		4 rotating wheels, high version	■	✓	◆				100 - 90 - 81 - 68,5 - 39 cm
	7300	MEBER-X		■		◆	✓	✓	✓	
	7210	FROG	With Trendelenburg	■		◆		✓		
	7210 PROOF		CERTIFIED with Trendelenburg	■	✓	◆		✓		
	7240	FROG PLUS	Plus with Trendelenburg	■		◆		✓		
	7240 PROOF		CERTIFIED Plus with Trendelenburg	■	✓	◆		✓		
	7260 PROOF	FROG LITE	CERTIFIED with Trendelenburg	■	✓	◆		✓		
	910	WINNER	Ambulance cot	■		◆				
	910 PROOF		CERTIFIED ambulance cot	■	✓	◆				
	914		Trendelenburg and Fowler	■		◆		✓	✓	
	914 PROOF		CERTIFIED with Trendelenburg and Fowler	■	✓	◆		✓	✓	
	7122/V	GECKO	Entry-level cot, green	■		◆				
	7122/AZ		Entry-level cot, light-blue	■		◆				
	7270	MBS	Ambulance cot for bariatric transport "MBS-01 MeBer Bariatric Stretcher"	■	✓	◆				

Materials: ◆ Stainless steel ◆ Painted steel ◆ Anodized Aluminium ◆ Painted Aluminium

## COMPARISON TABLE

# TROLLEYS

Max height	Max loading height	Min Loading Height	Folded height	Folded height (with stretcher)	Width	Rear axle width	Length	Min length	Folded length	Chair position min length	Weight load capacity	Loading capacity
84 cm	64 cm	60 cm	25 cm		44 cm	56 cm	197 cm				34,5 kg	220 kg
89 cm	72 cm	68 cm	25 cm		44 cm	56 cm	197 cm		200 cm		35 kg	220 kg
84 cm	64 cm	60 cm	25 cm		44 cm	56 cm	197 cm				36,5 kg	220 kg
89 cm	72 cm	68 cm	25 cm		44 cm	56 cm	197 cm		200 cm		37 kg	220 kg
82 cm	66 cm	60 cm	27 cm		56 cm		195 cm	169 cm			37 kg	250 kg
82 cm	66 cm	60 cm	27 cm		56 cm		195 cm	166 cm			38 kg	250 kg
83 cm	66 cm	60 cm	26 cm		56 cm		196 cm	167 cm			27,5 kg	250 kg
80 cm	65 cm	45 cm	26 cm		56 cm		194 cm	162 cm			28 kg	170 kg
80 cm	65 cm	62 cm	26 cm		56 cm		194 cm	155 cm			31 kg	170 kg
74 cm					57 cm		90 cm				30,5 kg	
80 cm	67 cm		27 cm		60 cm		193 cm		198 cm	123 cm	44,5 kg	220 kg
95 cm	66 cm	64 cm	31 cm		57 cm		197 cm	165 cm			40 kg	250 kg
99 cm	71 cm	69 cm	31 cm		57 cm		197 cm	165 cm			40,5 kg	250 kg
95 cm	66 cm	64 cm	31 cm		57 cm		197 cm	165 cm			42 kg	250 kg
99 cm	71 cm	69 cm	31 cm		57 cm		197 cm	165 cm			43 kg	250 kg
95 cm	66 cm	64 cm	31 cm		57 cm		197 cm	165 cm			44 kg	250 kg
99 cm	71 cm	69 cm	31 cm		57 cm		197 cm	165 cm			45 kg	250 kg
	66 cm	64 cm	31 cm		57 cm		197 cm	165 cm			48 kg	250 kg
	71 cm		31 cm		57 cm		197 cm	165 cm			49,5 kg	250 kg
	66 cm	64 cm	31 cm		57 cm		197 cm	165 cm			42 kg	250 kg
	71 cm	69 cm	31 cm		57 cm		197 cm	165 cm			41 kg	250 kg
123 cm	85,5 cm	0 cm			60 cm		197 cm	166 cm			45,9 kg	320 kg
92 cm	66 cm	60 cm	32,5 cm		56 cm		195 cm	169 cm			43 kg	250 kg
92 cm	66 cm	60 cm	32,5 cm		56 cm		195 cm	166 cm			44 kg	250 kg
92 cm	66 cm	60 cm	32,5 cm		56 cm		195 cm	169 cm			43 kg	250 kg
92 cm	66 cm	60 cm	32,5 cm		56 cm		195 cm	166 cm			44 kg	250 kg
91,5 cm	66 cm	60 cm	31 cm		56 cm		196 cm	167 cm			35,1 kg	250 kg
87,5 cm	65 cm	45 cm	26 cm		56 cm		194 cm	162 cm			34 kg	170 kg
87,5 cm	65 cm	62 cm	26 cm		56 cm		194 cm	155 cm			37 kg	170 kg
87,5 cm	65 cm	45 cm	26 cm		56 cm		194 cm	162 cm			37 kg	170 kg
87,5 cm	65 cm	62 cm	26 cm		56 cm		194 cm	155 cm			40 kg	170 kg
90,5 cm	62 cm		30 cm		56 cm		197 cm				37 kg	170 kg
90,5 cm	62 cm		30 cm		56 cm		197 cm				37 kg	170 kg
92 cm	60 cm		25 cm		56 cm		197 cm				47,5 kg	350 kg

# MERCURY

## LITE<sup>4</sup>

AMBULANCE COT IN ALUMINIUM  
WITH VARIABLE HEIGHT AND 4  
SWIVEL CASTORS



Bariatric



Pediatric



Trauma



Art. 7080/4RG Proof

### MERCURY LITE 4

Ambulance cot in aluminium with variable height and 4 swivel castors

Our "Mercury Lite" stretcher is available in the 4 swivel castors version, ensuring excellent manoeuvrability. This device is ideal for narrow spaces. The two front castors can be unblocked by means of an intuitive mechanism. They can revert to their original non-rotating position. The fastening of the stretcher inside of an emergency vehicle is ensured by our product Art. 955 Proof, certified UNI EN 1789 (see "Fastening Systems"). The stretcher is supplied with mattress, straps and object holder. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**


Art. 7082/4RG Proof

### MERCURY LITE 4

Ambulance cot in aluminium with variable height and 4 swivel castors – high version

Sharing the same features of the 7080/4RG Proof model, this version is easily loaded on higher emergency vehicles. Supplied with the same accessories as 7080/4RG. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**





## Ambulance cot Mercury Lite 4 accessories

See **ACCESSORIES**  
at page 90

### Technical Data

Art.	Colour	Material	Adjustable height	Trendelenburg	Intermediate height	Max height	Max loading height	Min loading height	Folded height	Width	Length	Min length	Weight	Loading capacity
7080/4RG PROOF	■	Painted aluminium	✓	✓	72 cm	95 cm	66 cm	64 cm	31 cm	57 cm	197 cm	165 cm	40 kg	250 kg
7082/4RG PROOF	■	Painted aluminium	✓	✓	75 cm	99 cm	71 cm	69 cm	31 cm	57 cm	197 cm	165 cm	40,5 kg	250 kg

# MERCURY

STAINLESS STEEL CERTIFIED AMBULANCE COT WITH VARIABLE HEIGHT



Bariatric



Pediatric



Trauma



## Art. 7070 Proof

### MERCURY

Stainless steel ambulance cot with variable height

Designed to fit a variety of configurations, this ambulance cot is provided with a multi-level system with 2 pistons on each axle. The fastening of the stretcher inside of an emergency vehicle is ensured by our product Art. 955 Proof, certified UNI EN 1789 (see "Fastening Systems").

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

## Art. 7072 Proof

### MERCURY

Stainless steel ambulance cot with variable height - high version

This model shares the same features of art. 7070, but also allows loading on higher emergency vehicles. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



### Unique Features

- ✓ Ergonomic control for loading operations and for set-up adjustments
- ✓ Adjustable backrest with lockable gas spring for 0° - 75° inclination
- ✓ Fully foldable 600 mm side rails
- ✓ HDPE bumpers (both on front and back legs)
- ✓ 2 fixed front castors Ø 200 mm with brakes
- ✓ 2 self-centring back swivel castors (360°) Ø 200 mm
- ✓ 2 load castors Ø 150 mm (easy loading system)
- ✓ Supplied with mattress (Art. 7013), 2 straps (Art. 608-A/MEB or 608-N/MEB) and 1 strap (Art. 612/AN-2/MEB or Art. 612/N-2/MEB)
- ✓ Supplied with multi-purpose object holder (Art. 12025).

Loaded  
position



Ambulance cots and trolleys



### Art. 955 PROOF

Fastening system for  
ambulance cots



## Mercury stretcher accessories

See **ACCESSORIES**  
at page 90

41

Art.	Colour	Material	Adjustable height	Trendelenburg	Intermediate height	Max height	Max loading height	Min loading height	Folded height	Width	Length	Min length	Weight	Loading capacity
7070 PROOF	■	Stainless steel	✓	✓	72 cm	95 cm	66 cm	64 cm	31 cm	57 cm	197 cm	165 cm	42 kg	250 kg
7072 PROOF	■	Stainless steel	✓	✓	75 cm	99 cm	71 cm	69 cm	31 cm	57 cm	197 cm	165 cm	43 kg	250 kg

# MERCURY4

VARIABLE HEIGHT STAINLESS STEEL CERTIFIED AMBULANCE COT WITH 4 SWIVEL CASTORS



Bariatric



Pediatric



Trauma



Ambulance cots and trolleys

42

Art. 7070/4RG Proof

## MERCURY4

Variable height stainless steel certified ambulance cot with 4 swivel castors

Our "Mercury" stretcher is available in the 4 swivel castors version, ensuring excellent manoeuvrability. This device is ideal for narrow spaces. The two front castors can be unblocked by means of an intuitive mechanism. They can be reverted to their original non-rotating position. The stretcher is supplied with mattress, straps and object holder.

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**


Art. 7072/4RG Proof

## MERCURY4

Variable height stainless steel certified ambulance cot with 4 swivel castors – high version

Sharing the same features of the 7070/4RG Proof model, this version is easily loaded on higher emergency vehicles. The stretcher is supplied with mattress, straps and object holder. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**





## Mercury 4 Ambulance Cot accessories

See **ACCESSORIES**  
at page 90

43

### Technical Data

Art.	Colour	Material	Adjustable height	Trendelenburg	Intermediate height	Max height	Max loading height	Min loading height	Folded height	Width	Length	Min length	Weight	Loading capacity
7070/4RG PROOF	■	Stainless steel	✓	✓	72 cm	95 cm	66 cm	64 cm	31 cm	57 cm	197 cm	165 cm	44 kg	250 kg
7072/4RG PROOF	■	Stainless steel	✓	✓	75 cm	99 cm	71 cm	69 cm	31 cm	57 cm	197 cm	165 cm	45 kg	250 kg

# MERCURY CINQUE

STAINLESS STEEL  
CERTIFIED AMBULANCE  
COT, 5 HEIGHTS, 4 SWIVEL  
CASTORS



Bariatric



Pediatric



Trauma



Art. 7075/4RG Proof

## MERCURY Cinque

Stainless steel certified ambulance cot, 5 heights, 4 swivel castors

Mercury Cinque stems from one of our most successful projects in the emergency sector, surely our most valued product, the Mercury 4RG. Mercury Cinque aims at further extending the operational capacity of the original model thanks to the introduction of the 5 level system, namely the 5 height feature. However, Mercury Cinque is much more than sheer evolution: it is the only stretcher with 5 heights and foldable front trolley and the only ambulance cot that can grant a load capacity of 250 kg on all 5 levels. Supplied with mattress, straps and a multi-purpose object holder. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**





Art. 7077/4RG Proof

## MERCURY Cinque

Stainless steel certified ambulance cot, 5 heights, 4 swivel castors - high version

Sharing the same features of the 7075/4RG Proof model, this version is easily loaded on higher emergency vehicles. The stretcher is supplied with mattress and straps. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



Ambulance cots and trolleys

## Mercury Cinque accessories

See **ACCESSORIES**  
at page 90

45

Art.	Colour	Material	Adjustable height	Trendelenburg	Height	Max loading height	Min loading height	Folded height	Width	Length	Folded trolley length	Weight	Loading capacity
7075/4RG PROOF	■	Stainless steel	✓	✓	95 - 87 - 79 - 62 - 39 cm	66 cm	64 cm	31 cm	57 cm	197 cm	167 cm	48 kg	250 kg
7077/4RG PROOF	■	Stainless steel	✓	✓	100 - 90 - 81 - 65 - 39 cm	71 cm	69 cm	31 cm	57 cm	197 cm	167 cm	49,5 kg	250 kg

# MERCURY

# CINQUELITE



Bariatric



Pediatric



Trauma


NEW  
PRODUCT

Art. 7095/4RG Proof

## MERCURY Cinque LITE

Certified Aluminium ambulance cot, 5 heights, 4 swivel castors

Entirely built in aluminium, this stretcher varies from the Mercury Cinque Art. 7075/4RG Proof and Art. 7077/4RG Proof (both in stainless steel) because of its light weight. This ambulance cot is provided with a multi-level system with 2 pistons on each axle. The aluminium tops are resistant, durable and easily cleanable. The fastening of the stretcher inside of an emergency vehicle is ensured by our product Art. 955 Proof, certified UNI EN 1789 (see "Fastening Systems"). Supplied with mattress, straps and a multi-purpose object holder. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

NEW  
PRODUCT

Art. 7097/4RG Proof



## MERCURY Cinque LITE

Certified Aluminium ambulance cot, 5 heights, 4 swivel castors – high version

Sharing the same features of the 7095/4RG Proof model, this version is easily loaded on higher emergency vehicles. The stretcher is supplied with mattress, straps and multi-purpose object holder. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**







Ambulance cots and trolleys

## Mercury Cinque Lite accessories

See **ACCESSORIES**  
at page 90

### Technical Data

Art.	Colour	Material	Adjustable height	Trendelenburg	Height	Max loading height	Min loading height	Folded height	Width	Length	Folded trolley length	Weight	Loading capacity
7095/4RG PROOF	■	Painted aluminium	✓	✓	94,5 - 87 - 79 - 68,5 - 39 cm	65 cm	63 cm	31 cm	57 cm	197 cm	167 cm	41 kg	250 kg
7097/4RG PROOF	■	Painted aluminium	✓	✓	100 - 90 - 81 - 68,5 - 39 cm	70 cm	68 cm	31 cm	57 cm	197 cm	167 cm	42 kg	250 kg

# MEBER-X

MULTI-LEVEL STRETCHER WITH MATTRESS AND STRAPS



Bariatric



Pediatric



Trauma

Ambulance cots and trolleys

48

Art. 7300

## MEBER-X

Multi-level stretcher with mattress and straps

The multi-level stretcher MEBER-X, with its X-shaped frame, belongs to our new White Line: it has been designed to facilitate loading, unloading, translation, lifting and dropping of a patient. Its main feature is the possibility to be adjusted at 9 different heights, supporting transfers and manoeuvrability under all conditions. Built in a special aluminium alloy, this stretcher has a 320 kg load capacity, ensuring maximum safety for the transport of bariatric patients. Further characteristics of this device include: 4 swivel castors with 200 mm diameter, 2 with brakes; loading bed in aluminium, easily cleanable; backrest inclination from 0° to 80°, lockable in any position; Fowler and Trendelenburg positions; retractable loading trolley. **Supplied with: Mattress Art. 7002, Straps Art. 608-A/MEB or Art. 608-N/MEB (2 pieces) and Strap Art. 612/AN-2/MEB or Art. 612/N-2/MEB (1 piece).**

Rotating front castors release lever



Rotating front castors

Ground position with extracted loading trolley




Fowler position



Adjustable headrest with pillow



Retracted trolley



Ambulance cots and trolleys

## MEBER-X Multi-level Stretcher accessories

See **ACCESSORIES**  
at page 90

49

### Technical Data

Art.	Colour	Material	Adjustable height	Trendelenburg	Fowler	Multi-level	Min height	Max height	Max loading height	Min loading height	Width	Length	Min length	Weight	Loading capacity
7300	■	Painted aluminium	✓	✓	✓	✓	40,5 cm	123 cm	85,5 cm	0 cm	60 cm	197 cm	166 cm	45,9 kg	320 kg

# MEBER-X

MULTI-LEVEL STRETCHER WITH  
MATTRESS AND STRAPS





Supplied  
with

**Art.  
7002**

Electro-welded  
4 pieces black  
mattress



**Art.  
612/N-2/MEB**

Black chest belt with  
4-buckles system



**Art. 608-N/MEB**

2 x 2-pieces black  
straps with metal  
buckle



**Art.  
980**

Cot fastening system composed of front and side fixing  
and rear safety hook



**Art.  
12018N**

P-120 - Telescopic Infusion Pole for ambulance cots



**Art.  
12120**

Lower - Black objects holder for MEBER-X  
(Size: 92x43 cm)



**Art.  
12122**

Upper - Black objects holder for MEBER-X  
(Size: 39x52 cm)




**Art.  
12100/N**

Black Pillow for MEBER-X headrest



**Art.  
12128**

Adjustable length system for MEBER-X headrest



Accessories

Ambulance cots and trolleys

51

## Technical Data

Art.		Intermediate Height 1	Intermediate Height 2	Intermediate Height 3	Intermediate Height 4	Intermediate Height 5	Intermediate Height 6	Intermediate Height 7	Intermediate Height 8	Intermediate Height 9
	Backrest	123 cm	114,5 cm	106 cm	96,5 cm	87 cm	77,5 cm	67 cm	55,5 cm	40,5 cm
7300	Foot-rest	91 cm	90,5 cm	88 cm	85 cm	81 cm	76 cm	69 cm	59 cm	43 cm
	Loading Area	84 cm	76 cm	67 cm	58 cm	48 cm	38,5 cm	29 cm	17 cm	0 cm

# FROG LITE PROOF

ALUMINIUM  
AMBULANCE  
COT WITH  
TRENDLENBURG



Bariatric



Pediatric



Trauma



Art. 7260 Proof

## FROG Lite Proof

Certified Aluminium Ambulance cot with Trendelenburg

Entirely built in robust painted aluminium, this extremely light ambulance cot is also provided with Trendelenburg. 4 Ø 200 mm castors in polyurethane (2 fixed with brakes and 2 self-aligning swivels) offer a high level of comfort on uneven floors and surfaces. This device is easy to clean and disinfect. Together with the fastening system 955 proof, the device can absorb the inertial stresses caused by sudden acceleration or deceleration of the rescue vehicle. Supplied with mattress and straps (Art. 7013, 2 pieces of Art. 608-A/MEB or 608-N/MEB and 1 strap Art. 612/AN-2/MEB or Art. 612/N-2/MEB).

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



Adjustable Loading bed

### Unique features

- ✓ Painted aluminium backrest
- ✓ Reclining backrest 0°-80°, lockable in any position
- ✓ Foldable side rails
- ✓ Vertical legs with bumpers
- ✓ Swivel castors with self-aligning system
- ✓ Adjustable Front trolley



Ground position



Trendelenburg Position

Supplied with



#### Art. 7013

3 pieces  
fire-proof  
electro-welded  
Skai mattress



#### Art. 608-A/MEB - 608-N/MEB



Orange or black strap in 2 pieces  
with metal buckle and Me.Ber. logo  
(155 cm length)

#### Art.

612/N-2/MEB  
612/AN-2/MEB



Orange/black adjustable chest belt  
with 4-buckles systems for ambulance  
cots or chairs

## Frog Lite Proof accesories

See **ACCESSORIES**  
at page 90

### Technical Data

Art.	Colour	Material	Height	Max loading height	Min loading height	Folded height	Width	Length	Min length	Weight	Loading capacity
7260 PROOF	■	Painted aluminium	91,5 cm	66 cm	60 cm	31 cm	56 cm	196 cm	167 cm	35,1 kg	250 kg

# FROG

YELLOW  
AMBULANCE  
COT WITH  
TRENDLENBURG

# TRENDLENBURG



Bariatric



Pediatric



Trauma

Art. 7210

## FROG

Yellow ambulance cot with Trendelenburg

Entirely built in powder-coated steel, this device is extremely resistant and reliable. The remarkable smoothness of the 4 Ø 200 mm castors (2 fixed with brakes and 2 self-aligning swivel castors) ensure high comfort to the patient. This model is equipped with Trendelenburg, from 0° to 20°.

Frog  
accessories

See **ACCESSORIES**  
at page 90

Art.	Colour	Material	Trendelenburg	Height	Max loading height	Min loading height	Folded height	Width	Length	Length with folded trolley	Weight	Loading capacity
7210	■	Painted steel	✓	92 cm	66 cm	60 cm	32,5 cm	56 cm	195 cm	169 cm	43 kg	250 kg
7210 PROOF	■	Painted steel	✓	92 cm	66 cm	60 cm	32,5 cm	56 cm	195 cm	166 cm	44 kg	250 kg



## Unique Features



- ✓ Loading bed in reinforced thermoplastic
- ✓ Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- ✓ Foldable side rails
- ✓ Vertical legs with cot bumpers
- ✓ Swivel castors with self-aligning system to ensure positioning of the cot on the emergency vehicle without manual turning of the castors
- ✓ Foldable Front trolley for volume reduction
- ✓ Trendelenburg.



Art. 7210PROOF

### FROG Proof

Certified yellow Ambulance cot with Trendelenburg

Sharing the very same characteristics of Art. 7210, this model is supplied with mattress and straps (Art. 7013, 2 straps Art. 608-A/MEB or 608-N/MEB and 1 strap Art. 612/AN-2/MEB or Art. 612/N-2/MEB). Combined with the fastening system 955 proof, this device can endure 10G stress. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

### Frog Proof accessories

See **ACCESSORIES**  
at page 90

Ambulance cots and trolleys

55

Supplied with



#### Art. 7013

3 pieces  
fire-proof  
electro-welded  
Skai mattress



#### Art. 608-A/MEB - 608-N/MEB

Orange or black strap in 2 pieces  
with metal buckle and Me.Ber. logo  
(155 cm length)



#### Art. 612/N-2/MEB 612/AN-2/MEB

Orange/black adjustable chest belt  
with 4-buckles systems for ambulance  
cots or chairs



Product Catalogue 2023

# FROG PLUS

## TRENDELENBURG

STAINLESS STEEL AMBULANCE COT WITH TRENDELENBURG



Bariatric



Pediatric



Trauma



Art. 7240

### FROG Plus

Ambulance cot with Trendelenburg

Entirely built in electro-polished stainless steel, this device features 4 Ø 200 mm castors (2 fixed with brakes and 2 self-aligning swivel castors). The leg's inclination was studied to ensure excellent stability to the stretcher. This model is equipped with Trendelenburg, from 0° to 20°, lockable in any position by means of a gas spring.

Frog Plus

See **ACCESSORIES**  
at page 90

Art.	Colour	Material	Trendelenburg	Height	Max loading height	Min loading height	Folded height	Width	Length	Length with folded trolley	Weight	Loading capacity
7240	■	Stainless steel	✓	92 cm	66 cm	60 cm	32,5 cm	56 cm	195 cm	169 cm	43 kg	250 kg
7240 PROOF	■	Stainless steel	✓	92 cm	66 cm	60 cm	32,5 cm	56 cm	195 cm	166 cm	44 kg	250 kg

## Unique Features

- ✓ Loading bed in reinforced thermoplastic
- ✓ Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- ✓ Foldable side rails
- ✓ Vertical legs with cot bumpers
- ✓ Swivel castors with self-aligning system to ensure positioning of the cot on the emergency vehicle without manual turning of the castors
- ✓ Foldable Front trolley for volume reduction



Art. 7240PROOF

## FROG Plus Proof

Certified Ambulance Cot with Trendelenburg

In addition to stability, comfort and all other features of Art. 7240, this model introduces excellent structure robustness, ensuring high resistance to stress on the emergency vehicle during transport. Supplied with mattress and straps (Art. 7013, n. 2 straps Art. 608-A/MEB or 608-N/MEB and 1 strap Art. 612/AN-2/MEB or Art. 612/N-2/MEB). **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



## Frog Plus Proof accessories

See **ACCESSORIES**  
at page 90

Ambulance cots and trolleys

57

Supplied with



### Art. 7013

Electro-welded 3 pieces fire-proof skai mattress



### Art. 608-A/MEB - 608-N/MEB



Black or Orange belt in 2 pieces with metal buckle and Meber Logo (155 cm)

### Art. 612/N-2/MEB 612/AN-2/MEB



Orange/black adjustable chest belt with 4-buckles systems for ambulance cots or chairs

Product Catalogue 2023

# FROG INC

INCUBATOR TRANSPORTER

Art. 7240/INC



Pediatric



Art. 7240/  
INC PROOF



Art. 7240/INC

## FROG Plus INC

Incubator transporter

This device belongs to the "FROG PLUS" line, sharing the very same features of the other models, such as reliability, safety and durability. Built in electropolished stainless steel, this device has no backrest, no side rails and no rigid loading bed. It is ideal for incubator transport on emergency vehicles.

Art. 7240/INC PROOF

## FROG Plus INC Proof

Incubator transporter "FROG PLUS" with fastening system configuration (10 G)

Integrating the configuration for the certified fastening system Art. 955 Proof, this model is ideal for incubator transport on emergency vehicles. Without backrest, side rails nor rigid loading bed, this device belongs to the "FROG PLUS" line.



Art.  
7240/INC PROOF  
+ 12048



Art.  
7240/INC PROOF  
+ 12050



Art.	Colour	Material	Height	Max loading height	Min loading height	Folded height	Width	Length	Length with folded trolley	Weight	Loading capacity
7240/INC	■	Stainless steel	82 cm	66 cm	60 cm	27 cm	56 cm	195 cm	169 cm	37 kg	250 kg
7240/INC PROOF	■	Stainless steel	82 cm	66 cm	60 cm	27 cm	56 cm	195 cm	166 cm	38 kg	250 kg
7260/INC PROOF	□	Painted aluminium	83 cm	66 cm	60 cm	26 cm	56 cm	196 cm	166 cm	27,5 kg	250 kg



#### Art. 7260/INC PROOF

### FROG Lite INC Proof

Incubator transporter with fastening system configuration (10 G)

This device belongs to the "FROG LITE" line, sharing the very same features of the other models, such as reliability, safety and durability. Built in painted aluminium, without backrest, side rails nor rigid loading bed. It is ideal for incubator transport on emergency vehicles.



#### Art. 12048



Certified 10G interface for incubator housing


#### Art. 12050



Extended interface for incubator housing



Art.  
7260/INC PROOF  
+ 12048



Art.  
7260/INC PROOF  
+ 12050

# WINNER

## GIN

AMBULANCE COT  
FOR INCUBATOR  
TRANSPORT

Art. 910/GIN



Pediatric



Art. 910/  
GIN PROOF



Art. 910/GIN

### WINNER GIN

Incubator transporter "WINNER"

Designed for the housing and transport of an incubator on emergency vehicles, this device belongs to the "Winner" line. Entirely built in anodized aluminium.

Art. 910/GIN PROOF

### WINNER GIN Proof

Incubator transporter "WINNER" with configuration for certified fastening system

Entirely built in anodized aluminium, this device shares the same characteristics of Art. 910/GIN and additionally features the configuration for the certified fastening system Art. 955 Proof.

Art.	Colour	Material	Height	Max loading height	Min loading height	Folded height	Width	Length	Length with folded trolley	Weight	Loading capacity
910/GIN	■	Anodized aluminium	80 cm	65 cm	45 cm	26 cm	56 cm	194 cm	162 cm	34 kg	170 kg
910/GIN PROOF	■	Anodized aluminium	80 cm	65 cm	62 cm	26 cm	56 cm	194 cm	155 cm	37 kg	170 kg

# FIXED TROLLEY FOR INCUBATOR HOUSING



Pediatric

Art. 7034

## Fixed trolley for incubator housing

Fixed height trolley with stainless steel frame, 4 castors  $\varnothing$  200 mm (2 with brakes) for incubator transport inside hospitals. Reduced size for optimal manoeuvrability. Supplied with universal incubator interface.



Pediatric

Accessories



Art. 12048



Certified 10G  
interface for  
incubator housing

Ambulance cots and trolleys

61

Art.	Colour	Material	Height	Width	Length	Weight
7034		Stainless steel	74 cm	57 cm	90 cm	30,5 kg

# TITANUS

STAINLESS STEEL VARIABLE HEIGHT CERTIFIED AMBULANCE TROLLEY



Ambulance cots and trolleys

62

Art. 7060 PROOF

## TITANUS

Stainless steel variable height certified ambulance trolley

Titanus represents the ideal solution among ambulance trolleys with variable height. The anchoring system for the anchoring of a stretcher to a TITANUS trolley holds the UNI EN 1789 certification. The variable height, as well as the leg folding mechanism for loading and unloading, is handled by a 2 piston system on each axle. 200 mm castors, together with a unique leg structure, confer stability, manoeuvrability and handiness to the TITANUS. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Art. 7062 PROOF

## TITANUS

Stainless steel variable height certified ambulance trolley – high version

Sharing the same features of the 7060 Proof model, this version is easily loaded on higher emergency vehicles. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**





Intermediate  
Position



Ground  
Position

Ambulance cots and trolleys

Accessories

Art. 1515 PROOF



Art. 1510 PROOF



Art. 956 PROOF



Fastening system for TITANUS  
trolleys (Certified UNI EN 1789)

Art.	Colour	Material	Adjustable height	Intermediate height	Max height	Max loading height	Min loading height	Folded height	Width	Rear axle width	Length	Folded length	Weight	Loading capacity
7060 PROOF	■	Stainless steel	✓	61 cm	84 cm	64 cm	60 cm	25 cm	44 cm	56 cm	197 cm	-	34,5 kg	220 kg
7062 PROOF	■	Stainless steel	✓	61 cm	89 cm	72 cm	68 cm	25 cm	44 cm	56 cm	197 cm	200 cm	35 kg	220 kg

# TITANUS<sup>4</sup>

STAINLESS STEEL VARIABLE HEIGHT CERTIFIED AMBULANCE TROLLEY  
WITH 4 SWIVEL CASTORS



Ambulance cots and trolleys

64

Art. 7060/4RG PROOF

## TITANUS<sup>4</sup>

Stainless steel variable height certified ambulance trolley with 4 swivel castors

The 4 swivel castors option is available for the TITANUS trolley: the TITANUS 4 shares the same features of Art. 7060 Proof. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Art. 7062/4RG PROOF

## TITANUS<sup>4</sup>

Stainless steel variable height certified ambulance trolley with 4 swivel castors – high version

Sharing the same features of the 7060/4RG Proof model, this version is easily loaded on higher emergency vehicles. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**



Intermediate  
Position

Art. 7060/4RG +  
1515 PROOF

Art. 7060/4RG +  
1510 PROOF



Front castors  
release lever

Front swivel  
castors

## Technical Data

Art.	Colour	Material	Adjustable height	Intermediate height	Max height	Max loading height	Min loading height	Folded height	Width	Rear axle width	Length	Folded length	Weight	Loading capacity
7060/4RG PROOF	■	Stainless steel	✓	61 cm	84 cm	64 cm	60 cm	25 cm	44 cm	56 cm	197 cm	-	36,5 kg	220 kg
7062/4RG PROOF	■	Stainless steel	✓	61 cm	89 cm	72 cm	68 cm	25 cm	44 cm	56 cm	197 cm	200 cm	37 kg	220 kg

# EVOX PROOF

CERTIFIED  
CHAIR COT

Unlimited  
Adjustments



Art. 7112 PROOF

## EVOX PROOF

EVOX Certified Chair Cot with mattress

The EVOX Chair Cot allows to place the patient in the most suitable position during rescue and transport. By converting into chair, this device ensures access to narrow entrances, tight aisles and elevators, where a regular Ambulance Cot won't fit. A gas piston system regulates the unlimited adjustments up to the sitting position, relieving the operator from any efforts to move the patient. The frame is entirely built in electro-polished stainless steel. The chair cot features Ø 200mm polyurethane castors, safety command to prevent accidental closure. The mattress is split in 4 separate sections (one for each portion of the cot's loading bed). All of these features, together with the cot's captivating design, are strong points of this device, which is experiencing increasing demand on the market. Its innovative load system, with its 4 swivel castors, facilitates the loading on a variety of supports and many different heights. Thanks to the mushroom head-lock together with the fastening system Art. 962, **the EVOX Chair Cot Art. 7112 PROOF is capable of enduring 10G stress in all directions, as per UNI EN 1789.** The ability to accommodate patients in the most appropriate position (depending on the seriousness of the injury) makes the EVOX chair cot a perfect fit for emergency waiting rooms, where the patients may be requested to queue for a long time. **Supplied with 1 mattress Art. 7009/A or 7009/N and 2 straps Art. 608-A/MEB or 608-N/MEB.**



Supplied with\*

**Art. 7009/N - Art. 7009/A**



Electro-welded mattress, orange or black, for EVOX Chair Cot

**Art. 608-A/MEB - Art. 608-N/MEB**



Orange or black strap in 2 pieces with metal buckle and Me.Ber. logo (155 cm length)

(\*) Supplied with 1 mattress and 2 straps (to be choosen among available colours)



Chair Cot without Mattress

Ambulance cots and trolleys

Accessories

**Art. 12016E**



Foldable telescopic Infusion Pole for EVOX Proof Chair Cot.

**Art. 12100/N**



Black Electro-welded headrest with MeBer logo for electro-welded mattress

**Art. 962**



EVOX PROOF Chair Cot fastening system (EN 1789)

Technical Data

Art.	Colour	Material	Height	Max loading height	Folded height	Width	Length	Folded length	Chair position min length	Weight	Loading capacity
7112 PROOF	■	Stainless steel	80 cm	67 cm	27 cm	60 cm	193 cm	198 cm	123 cm	44,5 kg	220 kg

# GECKO

ENTRY-LEVEL  
AMBULANCE  
COT

Art. 7122/V



Art. 7122/AZ

Art. 7122/V - Art. 7122/AZ

## GECKO

Green/light blue entry-level ambulance cot

Designed to optimize quality/ratio, this stretcher is entirely built in high-resistance painted steel. This model features all characteristics of regular ambulance cots: adjustable backrest, 4 Ø 200 mm castors (2 with brakes and 2 self-aligning swivel castors), foldable side rails. Supplied with straps.

## Gecko accessories

See **ACCESSORIES**  
at page 90

Ambulance cots and trolleys

### Accessories

Art.  
7016/A



One-piece electro-welded  
orange mattress in PVC

Art.  
965



Front fastening system for  
ambulance cot

Art.  
950/3

(H. 148 mm)



Rear fastening system for  
ambulance cot



Art.  
603/MEB



Supplied with straps

Art.	Colour	Material	Height	Max loading height	Folded height	Width	Length	Weight	Loading capacity
7122/V		Acciaio verniciato	90,5 cm	62 cm	30 cm	56 cm	197 cm	37 kg	170 kg

# MBS-01


## AMBULANCE COT FOR BARIATRIC TRANSPORT



Art. 7270 +  
Extension Wings art.12084 +  
Mattress art.12170 e art. 12172 +  
Mattress art.7013/N



Art. 7270 +  
Top art. 12080 +  
Mattress art. 7015 +  
Strap roller art. 560



Bariatric

### MBS-01 ambulance cot accessories

See **ACCESSORIES**  
at page 90

Art. 7270

## MBS-01

Ambulance cot for bariatric transport

Specifically designed for **bariatric transport**, the MBS-01 is built in **high resistance stainless steel**. The **reinforced** frame ensures a load **capacity of 350 kg**, which is **higher** than all cots in MeBer's range and features **2 front Ø 125x25 mm castors and 2 rear Ø 100x25 mm castors** for maximum safety during loading and unloading. The model features 4 Ø 200x50 mm castors, 2 with brakes and 2 self-aligning swivel castors, ensuring maximum comfort to the patient. The rigid loading bed is built in reinforced thermoplastic, making it easily cleanable. The device features a reclining backrest (0°-80°) operated by a gas spring, lockable in any position, and foldable side rails. The legs are protected by cot bumpers. Supplied with mattress and straps (Art. 7013, 2 straps Art. 608-A/MEB or 608-N/MEB and 1 strap Art. 612/AN-2/MEB or

612/N-2MEB). The higher load capacity makes this stretcher a **perfect fit for the MBP bariatric loading bed** (art. 12080).

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Accessories

### Art. 12080

MBP wide bariatric loading  
bed for ambulance cots

**Specify the cot model when  
ordering this accessory**



Art.	Colour	Material	Height	Max loading height	Folded height	Width	Length	Weight	Loading capacity
7270	■	Stainless steel	92 cm	60 cm	25 cm	56 cm	197 cm	47,5 kg	350 kg

# WINNER

AMBULANCE  
COT

**Art. 910PROOF**

**Art. 910**

**Art. 914PROOF**

**Art. 914**

**Art. 910 - Art. 910PROOF - Art. 914 - Art. 914PROOF**

## Winner

910 - Ambulance cot

914 - Ambulance cot with Trendelenburg and Fowler position

Entirely built in anodized aluminium, these models feature 4 150 mm castors (rear, self-aligning), adjustable backrest 0°-80°, foldable side rails, legs protected by cot bumpers. Sturdy and manoeuvrable.

**Also available in the "PROOF" certified version (UNI EN 1865 – UNI EN 1789).**

Art.	Colour	Material	Trendelenburg	Height	Max loading height	Min loading height	Folded height	Length with folded trolley	Width	Length	Weight	Loading capacity
910	■	Anodized aluminium		87,5 cm	65 cm	45 cm	26 cm	162 cm	56 cm	194 cm	34 kg	170 kg
910PROOF	■	Anodized aluminium		87,5 cm	65 cm	62 cm	26 cm	155 cm	56 cm	194 cm	37 kg	170 kg
914	■	Anodized aluminium	✓	87,5 cm	65 cm	45 cm	26 cm	162 cm	56 cm	194 cm	37 kg	170 kg
914PROOF	■	Anodized aluminium	✓	87,5 cm	65 cm	62 cm	26 cm	155 cm	56 cm	194 cm	40 kg	170 kg



# CARRY AND STAIR CHAIRS



MeBer has a whole line of carry- and stair-chairs. These devices are vital on emergency vehicles of all kinds, because they allow access to areas where cots can't fit (narrow spaces, elevators, tight stairs, high buildings). These models were developed thanks to the experience and support of rescue operators from all over the world: our range is highly

appreciated for the quality of the materials, for their movement accuracy and above all for their 10G certification (UNI EN 1865 – UNI EN 1789) which ensures safety when installed on emergency vehicles. The Stair-Chair EXTRA re-defined the market standards thanks to unique innovations implemented and to the use of extremely light materials.

Discover MeBer Chairs on-line  
Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website.  
You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/sedan-chairs.aspx>



# CARRY AND STAIR CHAIRS

# R



Article code	Name	Model
671	EXTRA ERGOLIFT	Powered Stair Chair
675	EXTRA ERGOLIFT FX	Fixed Frame Powered Stair Chair
670/BR	EXTRA	Track Assisted Stair Chair
676	EXTRA LITE	STANDARD Folding Stair Chair
678		BASIC Folding Stair Chair
679		ZENA Small Folding Stair Chair
668	KUMA	Black Powered Stair Chair
658/Y	QUICK	Folding black and yellow stair chair with 2 wheels
659/Y		Folding stair chair with footrest and armrest, 4 wheels
650	MAYA	Stair Chair
653		Stair Chair with 4 wheels, footrest and armrest
683	PROXIMA PRO	Padded stair chair with footrest and armrest

**Material:** Stainless steel Painted steel Anodized aluminium Painted aluminium

## COMPARATION CHART



Colour	10G Certification	Material	Seat height	Height	Folded height	Width	Min length	Max length	Folded thickness	Weight	Loading capacity
■	✓	Painted aluminium	50 cm	104 cm	-	55 cm	84 cm	150 cm	29 cm	29 kg	230 kg
■		Painted aluminium	50 cm	140 cm	-	55 cm	84 cm	150 cm	-	34 kg	230 kg
■	✓	Painted aluminium	50 cm	101 cm	-	55 cm	81 cm	139 cm	29 cm	17,2 kg	230 kg
■	✓	Painted aluminium	50 cm	101 cm	-	55 cm	78 cm	139 cm	29 cm	10,5 kg	200 kg
■	✓	Painted aluminium	50 cm	101 cm	-	55 cm	76 cm	139 cm	29 cm	9,9 kg	200 kg
■		Painted aluminium	56 cm	102 cm	-	47/50 cm	65 cm	122 cm	29 cm	13 kg	230 kg
■		Painted aluminium	50 cm	108 cm	-	52 cm	79 cm		24 cm	34 kg	181 kg
■	✓	◆	49 cm	90 cm	-	51 cm	71 cm	126 cm	23 cm	10 kg	170 kg
■	✓	◆	50 cm	90 cm	99 cm	51 cm	71 cm	126 cm	23 cm	14 kg	170 kg
■	✓	◆	51 cm	91 cm	-	52 cm	74 cm	128 cm	24 cm	10 kg	170 kg
■	✓	◆	52 cm	90 cm	98 cm	52 cm	75 cm	131 cm	24 cm	14 kg	170 kg
■	✓	◆	50 cm	118 cm	-	57 cm	90 cm	163 cm	-	21 kg	170 kg

Carry and stair chairs

73

# EXTRA ERGO **LIFT**

POWERED  
STAIR CHAIR



Carry and stair chairs

74

Art. 671

## EXTRA Ergolift

Powered Stair Chair

Extra Ergolift is the powered version of our Extra stair chair: it shares the very same structural characteristics of the Extra, such as a light aluminium frame, front and rear telescopic handgrips, coated armrest, front footrest, long and high-grip sliding tracks (allowing to balance the patient between 3 steps on 4 wheels), 2 rear Ø 200 mm castors and 2 swivel Ø 125 mm front ones. The real innovation of the Extra Ergolift lies in the electric motor that powers the tracks, thus facilitating the climb and descent of the patient on

the stairs and reducing the operator's physical effort. The motor is activated by means of a wireless console. The long-running battery ensures completion of the descent in case of sudden discharge. Extra Ergolift has a load capacity of 230 kg and may be used under adverse weather conditions. The device weighs 29 kg.

**Device certified according to the European harmonized safety standard UNI EN 1789.**



**Art. 12016-S**

Foldable Telescopic Infusion Pole for Extra and Extra Ergo-lift Chairs



**Art. 975**



Extra Ergolift fastening system



**Art. 671-001**

Black breathable dust sheet for Extra Ergolift chairs



**Art. 670-004**

Ø 200 mm castor's clamp with brake for Extra chair



**Art. 670-001**

Headrest for Extra Stair chair



**Art. 566**

Black bag with 3 two-pieces black straps and MeBer logo in orange for Extra stair chair



**Art. 12090**

MILWAUKEE Battery Charger M28C for M28RED re-chargeable battery for Extra Ergolift Stair Chair



**Art. 12092**

MILWAUKEE M28RED Lithium Battery 5.0 Ampere for Extra Ergolift Stair Chair




Supplied with



Battery housing support

Wireless console



Art.	Material	Length	Width	Height	Seat height	Weight	Loading capacity	Battery	Folded thickness	Autonomy	Charging time
671	Painted aluminium	84-150 cm	55 cm	104 cm	50 cm	29 kg	230 kg	28 V	29 cm	18 ramps by 20 steps with a load of 130 kg	max 1 hour

# EXTRA ERGO LIFT

## FX

FIXED FRAME POWERED STAIR CHAIR



Chair lock  
release lever



Art. 675

## EXTRA Ergolift FX


Fixed Frame Powered Stair Chair

Extra Ergolift FX is the fixed frame version of the powered stair chair Extra Ergolift. The two models share the same structural characteristics, such as light frame, telescopic front handgrips, wireless console, coated armrest, long and high-grip sliding tracks, footrest, small and simple electric motor. What characterizes this model is that the rescued patient can be carried to the final destination without transshipments. As a matter of fact, this device is homologated for transport on emergency vehicles as well as assistive vehicles.

Some of the main features include: high and stuffed backrest, for a safe and comfortable seating position.

Chest belt with 4-buckles system, 200 mm rear castors, telescopic grab bar, long-running battery. The EXTRA Ergolift FX has a 230 kg load capacity and may be used in adverse weather conditions. The device weighs 34 kg. 28V battery, duration of approximately 18 flights with 20 steps, carrying 130 kg. Fully charged in max 1 hour.

**Device certified according to the European harmonized safety standard UNI EN 1789.** (with fastening system Art. 988)



### Accessories



**Art. 988**



Fastening system for Extra Ergolift FX chair

### Supplied with



**Art. 12090**



MILWAUKEE Battery Charger M28C for M28RED re-chargeable battery for Extra Ergolift Stair Chair

**Art. 12092**



MILWAUKEE M28RED Lithium Battery 5.0 Ampere for Extra Ergolift Stair Chair

### Supplied with



**Art. 610**



Chest-belt with 4-buckles system for Extra Ergolift FX

**Art. 12118**



Mattress

**Art. 566**



Black bag with 3 two-pieces straps, black background with orange MeBer Logo for Extra Ergolift



Art.	Material	Length	Width	Height	Seat height	Weight	Loading capacity	Battery	Autonomy	Charging time
675	Painted aluminium	84-150 cm	55 cm	140 cm	50 cm	34 kg	230 kg	28 V	18 ramps by 20 steps with a load of 130 kg	max 1 hour

# EXTRA

STAIR  
CHAIR



Art. 670/BR

## EXTRA

Stair Chair


EXTRA is a track-assisted stair chair by MeBer. Entirely built in aluminium, this device is one of the lightest models available on the market. The inclination of both the backrest and the rear handgrips, as well as the special configuration of the telescopic handles, ensure a perfect weight distribution, facilitating the operator's work. The long high-grip tracks allow to balance the chair on 3 steps, making the EXTRA an extraordinarily safe and manoeuvrable chair during the descent of a seated patient. During transport on flat surfaces, the rear 200 mm castors and the front

125 mm ones ensure high comfort and handiness. Coming with 2 coated armrests and a telescopic grab bar, the comfort of this chair is unparalleled. To complete the features of this outstanding device, the EXTRA introduces a foldable footrest with magnetic blocking system and a few innovative bio-mechanical solutions, way ahead of the regular stair chairs. As an optional, the two rear castors can be supplied with brakes. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789**





Closed chair



Supplied  
with



**Art. 566**



Black bag with 3 two-pieces black straps and MeBer logo in orange for Extra stair chair

Fastening  
Systems



**Art. 975**



ERGOLIFT FASTENING SYSTEM. Fastening system for Extra Ergolift, Extra (foldable) and Zena



**Art. 960**



EXTRAFASTENING SYSTEM. Fastening system for EXTRA stair chair, certified UNI EN 1789



Accessories



**Art. 670-001**

Headrest and immobilizer for Extra stair chair



**Art. 670-002**

Black breathable dust sheet for Extra Ergolift chairs



**Art. 670-004**

Ø 200 mm castor's clamp with brake for Extra chair



**Art. 12016-S**

Foldable Telescopic Infusion Pole for Extra and Extra Ergolift Chairs



Art.	Material	Seat height	Height	Width	Min length	Max length	Folded thickness	Weight	Loading capacity
670/BR	Painted aluminium	50 cm	101 cm	55 cm	81 cm	139 cm	29 cm	17,2 kg	230 kg

# EXTRA LITE

FOLDABLE  
STAIR CHAIR



Art. 678



Art. 676

Art. 678

## EXTRA LITE Basic

Stair Chair

EXTRA LITE basic is MeBer's stair chair entirely built in light aluminium, currently among the lightest models on the market. The rear handgrips and the unique front handles' design ensure a perfect weight distribution, facilitating the operator's job. The backrest and seat positioning offer extraordinary comfort to the patient. During transfer on flat surfaces, the rear 200 mm castors and the front 125 mm ones offer great manoeuvrability. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.** (with fastening system Art. 988)

Art. 676

## EXTRA LITE Standard

Foldable 4 wheels Stair Chair with footrest and armrest

The EXTRA LITE standard shares the same features of Art. 678. Foldable footrest with magnetic lock and the armrests complete this device. **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.** (with fastening system Art. 988)



Supplied  
with

**Art. 566**



Black bag with 3 two-pieces black straps and MeBer logo in orange for Extra stair chair

Fastening  
Systems

**Art. 975**



ERGOLIFT FASTENING SYSTEM. Fastening system for Extra Ergolift, Extra (foldable) and Zena



**Art. 960**



EXTRA FASTENING SYSTEM. Fastening system for EXTRA stair chair, certified UNI EN 1789



Accessories

**Art. 670-002**



Black breathable dust sheet for Extra Ergolift chairs

**Art. 670-004**



Ø 200 mm castor's clamp with brake for Extra chair

## Technical data

Art.	Material	Seat height	Height	Width	Min length	Max length	Folded thickness	Weight	Loading capacity
676	Painted aluminium	50 cm	101 cm	55 cm	78 cm	139 cm	29 cm	10,5 kg	200 kg
678	Painted aluminium	50 cm	101 cm	55 cm	78 cm	139 cm	29 cm	9,9 kg	200 kg

# EXTRA LITE

## ZENA

SMALL FOLDING  
STAIR CHAIR



Art. 679

## EXTRA LITE ZENA

Small Folding Stair Chair, 4 wheels, with footrest, armrest and aluminium tops

82

Extra Lite Zena was designed to be employed in tight and narrow spaces, such as public transport vehicles (Buses, trains, airplanes, etc.), elevators and tight staircases. Entirely built in aluminium, this device is easily manoeuvrable and extremely light, among the lightest in Me.Ber's range. Seat and backrest are coated in expanded EVA, to ensure maximum comfort and easy cleaning. Foldable, featuring telescopic front handles and foldable rear handgrips, this Stair Chair is perfect when you need to reduce volume; the ZENA Stair Chair can be installed on an emergency vehicle by means of the **fastening system Art. 960, certified EN 1789**. Further characteristics of this device include 2 front 125 mm castors with brakes, 2 200 mm rear castors, coated armrests and foldable footrest. Weight: 13 kg

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.** (with fastening system Art. 988)

Art.	Material	Seat height	Height	Min length	Max length	Min width	Max width	Folded thickness	Weight	Loading capacity
679	Aluminium alloy	56 cm	102 cm	65 cm	122 cm	47 cm	50 cm	29 cm	13 kg	230 kg





Old buildings



Trains and buses




Airplanes



Service Elevators



Folded Chair



Volume comparison with a regular stair chair

Accessories

Art. 171J  
(Left)



Telescopic handgrips for Extra Stair Chairs

Art. 172J  
(Right)

Art. 612/AN-2/MEB  
Art. 612/N-2/MEB



Orange/black adjustable chest belt with 4-buckles systems for ambulance cots or chairs

Art. 842L



Rear handles lock and release system for Zena

Fastening Systems

Art. 960



EXTRA FASTENING SYSTEM. Fastening system for EXTRA stair chair, certified UNI EN 1789



Art. 975



ERGOLIFT FASTENING SYSTEM. Fastening system for Extra Ergolift, Extra (foldable) and Zena



Supplied with

Art. 566

Black bag with 3 two-pieces black straps and MeBer logo in orange for Extra stair chair



# KUMA

POWERED STAIR CHAIR



NEW  
PRODUCT

Art. 668

## KUMA

Black powered stair chair

Kuma is the entry level model in MeBer's line of powered stair chairs. Frame in aluminium, lithium-ion battery powered, front handgrips, ergonomic armrests, telescopic grab bar (with a small console for track speed control) and 4 castors (2 fixed rear Ø200 mm with brakes and 2 front swivels castors Ø125 mm). The electric motor operates the tracks, allowing an effortless ascent and descent of the stairs. Kuma has a 181 kg load capacity. Supplied with battery, battery charger, seat cushion, headrest, chest belt, ventral belt. Practical and easy to manoeuvre, the KUMA is extremely easy to use. Colour: black.

Art.	Material	Seat height	Height	Length	Width	Folded thickness	Weight	Loading capacity
668	Painted aluminium	50 cm	108 - 164 cm	79 cm	52 cm	24 cm	34 kg	181 kg



### Accessories

Art. 760M

Battery

NEW  
PRODUCT

Art. 761M

Battery charger

NEW  
PRODUCT

### Fastening Systems

Art. 977

NEW  
PRODUCT

Fastening system for KUMA Powered Stair Chair



# PROXIMA

## PRO

PADDED  
STAIR CHAIR




Art. 683

### PROXIMA Pro

Padded Stair Chair, built in aluminium and certified UNI EN 1865

Entirely built in aluminium, this stair chair is notable for its light weight. Padded fire-proof seat and backrest, coated in Skai to ensure easy cleaning and comfort. Equipped with 4 Ø 150 mm castors (2 rear swivel castors with brakes and 2 fixed front castors), 2 rear handgrips, 4 telescopic handles (2 front and 2 rear, Ø 35 mm for ease of movement) and foldable coated armrests for maximum comfort. Supplied with chest belt. Dedicated fastening system Art. 988 UNI EN 1789 available as accessory (see "Fastening Systems").

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Art.	Colour	Material	Seat height	Height	Width	Min length	Max length	Weight	Loading capacity
683		Painted aluminium	50 cm	118 cm	57 cm	90 cm	163 cm	21 kg	170 kg


Art. 988



Fastening System for Extra  
Ergolift FX and Proxima Pro

# QUICK659

FOLDABLE STAIR CHAIR



Art. 659/Y

## QUICK

Black/Yellow Foldable Stair Chair with 4 castors, footrest and armrest

The QUICK is a crucial device for rescue operators. Equipped with 4 castors (2 fixed Ø 125 mm rear castors and 2 swivel Ø 100 mm front castors with brakes), foldable armrests, foldable footrest, 4 rear handgrips and 2 front telescopic handles. The backrest and seat are built in easily-cleanable thermoplastics. Supplied with black bag containing 3 straps (Art. 590/N). **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

86

Art.	Colour	Material	Alt. seduta	Height	Alt. chiusa	Width	Min length	Max length	Folded thickness	Weight	Loading capacity
659/Y		Acciaio verniciato	50 cm	90 cm	99 cm	51 cm	71 cm	126 cm	23 cm	14 kg	170 kg
658/Y		Acciaio verniciato	49 cm	90 cm	-	51 cm	71 cm	126 cm	23 cm	10 kg	170 kg



# QUICK658

FOLDABLE STAIR CHAIR



Art. 658/Y

## QUICK

Black/Yellow Foldable Stair Chair with 2 castors

Designed for a safe transport. Easily closed/opened. Thanks to a unique mechanism, the chair cannot be closed by mistake. Equipped with 2 Ø 125 mm castors, two rear handgrips and 2 front telescopic handles, this device is extremely manoeuvrable. The backrest and the seat are built in easily-cleanable thermoplastics. Supplied with black bag containing 3 straps (Art. 590/N). **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Carry and stair chairs

87

Supplied  
with

Art. 590/N



Black bag containing 3  
black straps

Fastening  
Systems

Art. 618



Fastening System for Stair chair  
(See FASTENING SYSTEMS)

Accessories

Art. 612/AN-2/MEB Art. 612/N-2/MEB



Orange/black adjustable chest  
belt with 4-buckles systems for  
ambulance cots or chairs

Product Catalogue 2023

# MAYA653

FOLDABLE STAIR CHAIR



Art. 653

## MAYA

Foldable Stair Chair with 4 castors, footrest and armrest

Designed to help rescue operators during rescue procedures. Entirely built in anodized aluminium. Equipped with 4 castors (2 swivel castors with brakes Ø 100 mm and 2 fixed castors Ø 125 mm), a unique open/close system, 4 rear handgrips and 2 front telescopic handles, backrest and seat in thermoplastics, footrest. Supplied with black bag containing 3 straps (Art. 590/N). **Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Art.	Colour	Material	Seat height	Height	Folded height	Width	Min length	Max length	Folded thickness	Weight	Loading capacity
653	■	Anodized aluminium	52 cm	90 cm	98 cm	52 cm	75 cm	131 cm	24 cm	14 kg	170 kg
650	■	Anodized aluminium	51 cm	91 cm	-	52 cm	74 cm	128 cm	24 cm	10 kg	170 kg

# MAYA650

FOLDABLE STAIR CHAIR



Art. 650

## MAYA

Foldable Stair Chair


Designed to serve as a crucial device for patients' transport, this model is entirely built in anodized aluminium. Equipped with 4 handgrips (2 rear handgrips and 2 front telescoping handles), 2 fixed castors Ø 125 mm. Intuitive open/close mechanism to avoid accidental closure of the chair. Supplied with black bag containing 3 straps (Art. 590/N).

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Carry and stair chairs

89


**Art. 590/N**



Supplied with

Black bag containing 3 black straps


**Art. 618**



Fastening Systems

Fastening System for Stair chair (See FASTENING SYSTEMS)

**Art. 612/AN-2/MEB** **Art. 612/N-2/MEB**



Accessories

Orange/black adjustable chest belt with 4-buckles systems for ambulance cots or chairs

# 4

## Accessories

[illegible]





Alternatively, go to:  
<https://www.meber.it/en/products/accessories.aspx>

\* Compatible with Art. 12174 and Art. 12176

# MATTRESSES

ELECTRO-WELDED MATTRESSES FOR AMBULANCE COTS



Art. 7014/A

## SKAI MATTRESS

Electro-welded 4 pieces fire-proof skai mattress

Ambulance cot mattress in fire-proof skai and padding. Entirely electro-welded and water-resistant, this item is suitable for all MeBer ambulance cots (EVOX excluded).

**Size: cm 186 x 47 x 9**



Art. 7013/A

Art. 7013/N



Art. 7013

## SKAI MATTRESS

Electro-welded 3 pieces fire-proof skai mattress

Same features of Art. 7014 but suitable for all MeBer ambulance cots with Trendelenburg, including CRONOS (EVOX and Art. 1515 excluded). **Size: cm 186 x 47 x 9**

Art. 1880-002

## SKAI CUSHION FOR ICARUS

Skai cushion for Icarus

Cushion for ICARUS. Inner fireproof padding, coated in anti-tear, water-resistant and fireproof skai. Double-stitched seams. **Size: cm 55x35x5**

Art. 1880-001

## MATTRESS FOR ICARUS

Mattress for Icarus

Mattress for Icarus: inner fireproof padding, coated in anti-tear, water-resistant and fireproof skai. Double-stitched seams. **Size: cm 180x55x5.**

Art. 7016/A

## PVC MATTRESS

Orange, PVC one-piece electro-welded mattress

**Size: cm 186 x 47 x 5**





**Art. 7009/A - Art. 7009/N**

## MATTRESS

Electro-welded mattress, orange or black, for EVOX Chair Cot

Electro-welded mattress for EVOX Chair cot in Skai and inner fire-proof padding. Equipped with straps. Split in two pieces to better adapt to the cot's various configurations. The sides are raised for a better containment of the patient.

**Dimensioni: cm 131x45x8 - 46x42x5**

**Art. 12100/N**

## CUSHION

Electro-welded, modeled black cushion with MeBer logo – for electro-welded mattress

Electro-welded Skai cushion with inner fire-proof padding for EVOX Chair Cot: it features eyelets for anchoring to the Mattresses Art. 7009 - 7013 - 7014. The unique shape makes it ideal for head containment and support.

**Max size: cm 25x37x12 (LxPxH)**



**Art. 7009-001**

## SET

Set with 4 black straps with VEL-CRO® for the anchoring of Art 7009 - 7013 - 7014 to the EVOX Chair Cot




**Art. 7002**

## MEBER-X MATTRESS

4-pieces electro-welded black mattress

Sharing the same features of Art. 7014 but suitable for all MeBer ambulance cots with Trendelenburg, including CRONOS (EVOX and Art. 1515 excluded).

**Size: cm 180 x 47 x 8**



**Art. 12128**

Adjustable length headrest system for MEBER-X Ambulance cot.

# BARIATRIC LI

ACCESSORIES FOR BARIATRIC TRANSPORT



Art. 12080

## MBP

Bariatric loading bed for FROG, MERCURY, GECKO and MBS-01 ambulance cots, supplied with straps and fastening system

The MBP loading bed is compatible with all Frog, Mercury and Gecko line ambulance cots and allows to transport bariatric patients in complete safety. The stainless steel structure and the aluminium loading bed confer extreme stability and lightness at the same time. Equipped with straps (with roller) on the sides and 4 hooks for the anchoring to the cot. **Size: cm 137,5 x 93,5 x 13. Weight: kg 18**

Art. 7070/4RG PROOF + 12080



Art. 12080-001

## ASSEMBLY KIT

Bariatric loading bed assembly kit

Bariatric loading bed MBP (art.12080) assembly kit for Mercury, Mercury Lite, Frog, Frog Plus, Frog Lite and Gecko ambulance cots. The kit includes all fastening accessories and is easily installed.



# LINE



Art. 7015

## MBP MATTRESS

Foldable Skai mattress for MBP loading bed

Mattress for bariatric loading bed (MBP) in expanded polyurethane and black Skai coating (fireproof, water-resistant, anti-tear). Foldable in two pieces for easy storage. Equipped with 6 VELCRO® bands for the anchoring to the loading bed.

**Size: cm 173x96x5,5**

Art. 7070/4RG PROOF  
+ 12080 + 7015



Accessories

Art. 7080/4RG PROOF  
+ 12154

Art. 12154

## BEAR PODS

"Bear Pods" Bariatric containment system for ambulance cots

The Bear Pods containment system is an accessory for the transportation of bariatric patients on ambulance cots. The Bariatric Equalizing Abdominal Restraint (BEAR) is a strong light weight nylon mesh fabric which is UV resistant, chemical resistant, easy to clean, and easy to store. The purpose of this device is to facilitate the transport of bariatric patients in total safety. Bear Pods fits any type of ambulance cot thanks to its fittable straps that must be anchored to the sides of the cot's frame.

**This product can withstand up to 300 kg.**



95



Bariatric



Product Catalogue 2023

# BARIATRIC

ACCESSORIES FOR BARIATRIC TRANSPORT



Bariatric

**NEW**  
PRODUCT

Art. 12084

Art. 7270 + 12084 +  
12170 + 12172

## BARIATRIC WINGS SET

Bariatric wings set for ambulance cots, complete with belts, belt roller and fastening kit

The set allows you to transform the stretcher in a bariatric rescue device quickly and independently. Thanks to the reduced volume (and a weight of merely 15 KG) the extensions can be kept on the ambulance to always be within reach. Their employ allows you to face any scenario of bariatric patient rescue, in total safety and with appropriate equipment. The kit includes two extension wings that increase the width of the stretcher by over 30 cm. The assembly and disassembly are very easy and intuitive and do not require any tools. When mounted on our MBS01 the total load capacity amounts to 335 kg.

**Size: 1365 x 205 x 85 mm (wing alone)**

**Size: 1365 x 205 x 134 mm (wing + closed fastening system, non-mounted)**

The bariatric wings are equipped with straps and belt roller, as well as a mounting kit for maximum safety.

The kits can be bought separately (Art. 12086 right wing, Art. 12088 left wing) or together (Art. 12084).

Cushions are not included in the kit (Art. 12170 black mattress for right bariatric wing and Art. 12172 black mattress for left bariatric wing).

# NE

## AVAILABLE SETS



Art. 7270 + 12084 +  
12170 + 12172 + 7013/N



Bariatric



### Art. 12084

Bariatric wings set for ambulance cot, complete with belts, belt roller and fastening kit.



Bariatric



### Art. 12086

Right bariatric wing for ambulance cot, complete with belts, belt roller and fastening kit.



Bariatric



### Art. 12088

Left bariatric wing for ambulance cot, complete with belts, belt roller and fastening kit.



Bariatric



### Art. 12170

Right mattress

Comfortable black mattress for right bariatric wing. It perfectly fits the support and ensures maximum comfort to the patient.



Bariatric



### Art. 12172

Left mattress

Comfortable black mattress for left bariatric wing. It perfectly fits the support and ensures maximum comfort to the patient.



Bariatric



### Art. 12174

Right interface kit for bariatric wings. Mercury Lite and Mercury stainless stell.



Bariatric




### Art. 12176

Left interface kit for bariatric wings. Mercury Lite and Mercury stainless stell.



# HI-COM

THERAPEUTIC GEL HIGH COMFORT MATTRESS




Art. 7018/A

## HI-COM

Orange mattress in slow-recovery expanded polyurethane and therapeutic gel

Mattress in fluid therapeutic gel composed of 3 different layers of expanded polyurethane (1 slow-recovery layer). Equipped with containment wings and two stiff areas: head and feet. The sacral/dorsal area features two therapeutic gel sheets (cm 40x40xH 150) to confer comfort. The three expanded polyurethane areas absorb vibrations caused by transport on the emergency vehicle, while the two therapeutic gel areas support the patient's weight, ensuring maximum comfort. Hand-made in waterproof, fireproof, anti-tear Skai, with double-stitched seams with raised sides. The gel sheets are glued to a unique "pouch" that ensures durability. Fit for any type of ambulance cot.

**Size: cm 179x47x13 – Weight kg 9,7**



Art. 7018/A-001

## CUSHION

Orange Skai cushion for mattress



Made in slow-recovery expanded polyurethane and coated in fireproof, water-resistant and anti-tear Skai, this cushion fits all of our ambulance cots. Equipped with two rubber bands, allowing positioning in the best position for the patient, this cushion prevents cervical pain and aids blood flow in the area. Comfort is ensured on any kind of transport. **Size: cm 48 x 35 x 7,5**



# AIR-COM

ANTI-BEDSORE MATTRESS FOR LONG-HAUL



Art. 12106 +  
Art. 12108 +  
Art. 7013/A +  
Art. 7070/4RG

**Art. 12106**

## AIR-COM

Anti-bedsore segmented mattress for ambulance cots

Overlay mattress in PVC, divided in 18 individual inner tubes. Specifically designed for patients at risk of developing bedsores. AIR-COM fits any mattress for ambulance cots.

The underlying part of the AIR-COM, made in skid-proof fabric, is equipped with 3 eyelets to facilitate anchoring to the cot.

Air compressor with pressure regulator is available (Art. 12108).

**Size: 180 x 46 x 10 cm. Load Capacity: 150 kg.**

**Art. 12108**


AIR-COM Mattress compressor



Accessories

99

**Alternate segments  
diagram of operation**



Product Catalogue 2023

# INFUSION POLES



Art. 12018N

## P-120

Telescopic Infusion Pole for ambulance cots

Foldable Telescopic Infusion pole in aluminium for ambulance cots

**Size:** H min cm 45 - H max cm 128  
1<sup>st</sup> step height: cm 70  
2<sup>nd</sup> step height: cm 98



Art. 12016-E

## EVOX INFUSION POLE

Foldable telescopic infusion pole for EVOX Proof Chair cot

**Size:**  
H min cm 43 - H max cm 72,5



Art. 12016-S

## EXTRA INFUSION POLE

Foldable telescopic infusion pole for EXTRA and EXTRA ERGOLIFT Stair Chairs

**Size:**  
H min cm 43 - H max cm 72,5



30 mm

Art. 936



46 mm

Art. 936/A



34 mm

Art. 12010

Art. 936 - 936/A - 12010

## INFUSION POLE

Art. 936 for Winner and Frog ambulance cots  
Art. 936/A for Art. 1510 Proof and Art. 1515 Proof  
Art. 12010 for Mercury and Gecko ambulance cots

Telescopic infusion pole for ambulance cots. The Pole can be installed anywhere on the cot thanks to a quick anchoring system.

**Size:** H min cm 47,5 - H max 130,5 cm




Art. 939

## Object Holder

Ambulance cot Object Holder

Support for medical equipment such as defibrillator or monitors to be anchored to the cot's frame. Supplied with straps for the fixing of the medical equipment. Foldable for a better storage.

**Size: cm 34x56x42h fully open, cm 11 depth when closed**




Pediatric

Art. 12050

## INCUBATOR INTERFACE

Extended interface for incubator housing




Pediatric

Art. 12048

## INCUBATOR INTERFACE

10G Interface for incubator housing



Accessories

101

Product Catalogue 2023



Art. 12140

## EVOLUTION

Automatic emergency hammer



Art. 12142

## CLASSIC GLOW

Compact emergency hammer

ACCESSORIES

# PATIENT TRANSFER SYSTEM



Art. 11506

## SMALL TRANSFER BOARD WITH COVER

Small rigid slide board for CPR with cover

Size: cm 59 x 32



Art. 11508

## TRANSFER BOARD WITH COVER

Rigid slide board for CPR with cover

Size: cm 76 x 32

Art. 15244

## FOLDABLE SCANGLISS

Foldable transfer board



Foldable transfer board, to be placed underneath a patient and facilitate sliding from cot to bed in complete safety. Split in two (foldable) sections, this device is radiolucent and a-magnetic.


Consisting in an inner polyethylene layer, padded in amorphous material (1,5 cm thick) and wrapped in a sliding nylon sheet. SCANGLISS is CE certified according to CEE 93/42 EEC Medical Devices. **Maximum Load: 180 kg.**

Size 180 x 50 cm

Art. 15244-001

## COVER

Spare nylon cover for Scangliss







Art. 707

## SLASH

Transfer sheet with 6 handles

Slash is the brand new transfer sheet by MeBer. This device facilitates patients' transfer from cot to bed and vice-versa. Equipped with 6 handles, this sheet features a fireproof, water-resistant and anti-tear Nylon. Fit for ambulance cots, the SLASH is easy-to-clean and has a small volume, to ensure easy handling during transfers.

Fit for all ambulance cots. **SLASH does not work as a transport means for injured patients.**

**Size: 188 x 50 cm**



Art. 12098

## BLACK CUSHION

Black cushion with backrest and seat for EXTRA Stair Chairs

Cushion for Stair Chair, designed to ensure maximum comfort during transport and protect the rigid thermoplastics from constant stress. Made in expanded polyurethane (slow-recovery), this cushion is coated in fire-proof, anti-tear Skai fabric. Removable and easy-to-clean. Specifically designed for Extra stair chairs, it can fit any chair in MeBer's range thanks to the Velcro® straps.



Art. 671-001 | Art. 670-002

## COVER

Black breathable dust sheet for Extra Ergolift chairs

**Size: cm 55 x 100 x 27 (LxHxD)**

## ACCESSORIES

# CARRYING BAGS AND CASES



Art. 16102-001

### TOBOGA CASE

Divisible case for Toboga

Size: cm 121x64x24



Art. 12144/N

### CARRYING BAG

Spinal board carrying bag

Spinal board carrying bag with pouch for Spider Fastrap and removable head immobilizer, colour black.

Size: cm 184x43



Art. 16150-001

### LIFTABLE CASE

Liftable case for transport

Designed to facilitate rescue operations, this case for the ERGON scoop stretcher ensures maximum operativity in hostile environment. Equipped with 4 shoulder straps. The ERGON stretcher (disassembled) fits with all of its accessories inside of the case, strongly reducing its volume. The EN566 hook in the centre ensures safe winch operation. **Size: cm 124x32x16**



Art. 16150-002

### STORAGE BAG

Storage bag

Measuring as little as 11 cm in depth, this bag is ideal for ERGON storage, together with its accessories. It was designed for storage purposes or equipped on emergency vehicles. Once opened on the ground, the bag works as an insulating layer, allowing to safely prepare medicines or operate electro-medical devices. **Size: cm 122x47x11**



Art. 12084-001

## CARRYING BAG

Bariatric wings carrying bag

Bariatric wings carrying bag (Art. 12086, Art. 12088), black.

Size: cm 136x18x21



Art. 660/CS - Art. 662/CS

## CARRYING BAGS

Carrying bags for stretchers

Grey nylon bags for Art. 660 / Art. 662.  
Equipped with handles and shoulder strap for transport.

Size art. 660/CS: cm 100 x 16 x 20

Size art. 662/CS: cm 213 x 13 x 20



Art. 660/CS



Art. 662/CS



Art. 630/CS

## CARRYING BAG

Scoop stretcher carrying bag

Nylon orange carrying bag for scoop stretchers (split in half) with handle.

Size: cm 130x47



Art. 9104-001

## CARRYING BAG

Blue carrying bag for NAXOS vacuum mattress

Size: cm 85xø36



# OBJECT HOLDER



Velcro straps for  
anchoring

Pouch for documents

Zip Pouch

Accessories

106

Art. 12025

## OBJECT HOLDER

Multi-purpose object holder for ambulance cots

Anti-tear nylon object holder, anchorable to any MeBer ambulance cot. Consisting in a closed zip pouch for small items and a bigger pouch for documents, the holder is anchored to the frame of the cot's trolley with Velcro® straps.

**Size: cm 50x37**





ACCESSORIES

# OBJECT HOLDER



Art. 12024

## CASE

Empty case for spider strap and head-immobilizer

Cordura case with handle for spider strap and head-immobilizer. It ensures volume optimization inside of the emergency vehicle and prevents the Velcro® in the head-immobilizer from decaying. Also fit for the storage of Art. 9128 and 9130.

**Size: cm 41x26x13**



Pediatric

Art. 12034

## CASE

Red case for BABY FIX and VACUUM BABY FIX

Red Nylon case with handle and double-zip for the transport of "Baby Fix" and "Vacuum Baby Fix" pediatric devices. Ideal for volume optimization.

**Size: cm 41x26x13**

# OXIMEB

NYLON TANK HOLDER

Inner padding to ensure protection against impact.  
Equipped with handle for transport and lateral  
pouch for tube and mask.

Art. 938

## OXIMEB

Tank holder for  
2 litres tank

Suitable for the housing of 2 litres tank. Thanks to its  
straps system, it can be mounted on a cot's side rail or  
on the backrest of a stair chair.

**Weight: 0,3 kg**



Accessories

108

Art. 12054

## OXIMEB MAXI

Tank holder for  
3 litres tank

Suitable for the housing of 3 litres tank. Thanks to its  
straps system, it can be mounted on a cot's side rail  
or on the backrest of a stair chair.

**Weight: 0,3 kg**



# COVER OX

## PROFESSIONAL TANK HOLDERS

**Made in high-resistance, anti-tear and easy-to-clean material, these covers are equipped with 2 rubber handles and two knurled plastic bases, allowing both horizontal and vertical operation.** They also feature a padded shoulder strap and 2 auxiliary carabiners to facilitate transport. The 6 eyelets placed on the sides constitute a perfect passage for the oxygen

tube. A big central pouch houses a variety of oxygen therapy masks as well as the oxygen extension tube.

Two belt loops with plastic hooks allow anchoring of these covers to any ambulance cot. Built with robust zips, this is a safe, easy-to-clean medical device, designed to last for a long time.

Art. 12072/BR - Art. 12072/BG

### COVER OX 5

Professional tank holder  
for 5 litres tank

Professional tank holder for 5 litres tank. Blue with yellow profiles and yellow Logo or blue with red profiles and orange Logo.

**Size: (LxHxD) 67 x 22 x 24 cm. Weight: 1,2 kg**



Art. 12070/BG - Art. 12070/BR

### COVER OX 2-3

Professional tank holder  
for 2 or 3 litres tank

Professional tank holder for 2 or 3 litres tank. Blue with yellow profiles and yellow Logo or blue with red profiles and orange Logo.

**Size: (LxHxD) 60 x 22 x 21 cm. Weight: 1,1 kg**




# IMMOBILIZATION

5

The management of an injured patient is always a critical phase of all emergency operations, especially the immobilization process. Thanks to high quality products, the operator may ensure a quick and efficient stabilization, avoiding further worsening of a patient's condition. Completing the classic Galaxy spinal board range, ever appreciated for its simplicity and build quality, many innovative products complete MeBer's range of immobilization devices, especially to cover all pediatric necessities.

Immobilization

110



Spinal boards

Head  
immobilization




# TION




Discover MeBer products on-line.  
Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website.  
You might be required to install a free application to navigate.

Oppure digita: <https://www.meber.it/en/products/immobilization.aspx>



Immobilization



**Vacuum  
mattresses**



**Splints**



**Cervical Collars**



**Spinal  
immobilizers**

# GALAXY

SPINAL BOARDS FOR ADULTS



Art. 9010  
Yellow

Art. 9012  
Orange

Art. 9016  
Red



Trauma

Art. 9010 - Art. 9012 - Art. 9016

## GALAXY

112

Yellow, orange and red Spinal Boards

Spinal boards in high-resistance plastic material. The GALAXY was designed paying special attention to the daily complications occurring in rescue operations. For a fast and practical fastening of the head immobilizer, the board features 4 eyelets on the top, as well as a hollow area on the support area of the board. This ensures perfect immobilization and avoids damage to the straps during ambulance transfers. The board features 20 handles for transport and is supplied with 3 quick-release straps (Art. 605/MEB). X-Ray translucent.

**Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.**

Accessories



Art. 12144/N  
Spinal board  
carrying case



Art. 690

## SPIDER

### Immobilization strap system for spinal board

This strap system ensures perfect immobilization of the patient. Thanks to its transversal sliding straps and VELCRO fastening, it suits any spinal board as well as any patient.



Art. 629

## JELLY

### Standard orange head immobilizer for spinal board

Adjustable to any spinal board model, this head immobilizer is made up of expanded close-cells material and coated in vinyl.



Art. 625

## XEROS

### Standard adjustable head immobilizer

Coated in antibacterial, liquids-proof vinyl. Expanded close-cells material inside. This head-immobilizer is unique thanks to the internal aluminium layer that allows to shape the device to the patient's head. The fastening takes place by means of a VELCRO system and a strap system to stabilize chin and forehead.



Art. 628

## GREEN

### Universal head immobilizer

Essential for head immobilization, this device is easily adjustable to spinal boards. Made in close-cells expanded material and coated in exclusive vinyl

Head Immobilizer mounting detail. Rear view.



Head Immobilizer mounting detail. Front view.



Under X-rays, the Galaxy boards can be spotted as a longitudinal shadow, allowing the performance of radiographic examinations without moving the patient. The positioning of the internal reinforcements of the boards was designed not to hinder the identification of organ injuries.

Art.	Colour	Length	Width	Thickness	Weight	Load
9010		183 cm	41 cm	5 cm	5,5 kg	150 kg
9012						
9016						



# DUCK

PEDIATRIC SPINAL BOARDS



Art. 9022  
Yellow

Art. 9024  
Orange



Pediatric



Trauma



Art. 9022 - Art. 9024

## DUCK

Pediatric spinal board, orange and yellow

Designed for the immobilization and transport of pediatric patients **up to 120 cm**, made in high resistance plastic material. The DUCK was designed paying special attention to the daily complications occurring in rescue operations. For a fast and practical fastening of the head immobilizer, the board features 4 eyelets on the top, as well as a hollow area on the support area of the board. This ensures perfect immobilization and avoids damage to the straps during ambulance transfers. The board features 10 handles for transport.

The DUCK allows the performance of radiographic examinations with X-Rays and can be supplied in 2 colours: yellow and orange. The straps are not included.

Art.	Colour	Length	Width	Thickness	Weight	Load
9022		125 cm	41 cm	5 cm	3,7 kg	100 kg
9024						





Hollow area on the head space



Art. 692

## SPIDER BABY Pediatric Spinal Board immobilization straps

Strap system for the immobilization of a pediatric patient on a DUCK spinal board. Thanks to the transversal sliding straps, this device can be easily adjusted to any patient.



Pediatric



Trauma



See page 107  
ACCESSORIES



Art. 627

## KIDS Pediatric padded head immobilizer

Padded head immobilizer for pediatric spinal boards, made in close-cells expanded material, coated in vinyl.



Pediatric



Trauma

# HEAD IMMOBILIZER



Art. 629

## JELLY

**Standard orange head immobilizer for spinal board**

Adjustable to any spinal board model, this head immobilizer is made up of expanded close-cells material and coated in vinyl.

Art. 625

## XEROS

**Standard adjustable head immobilizer**

Coated in antibacterial, liquids-proof vinyl. Expanded close-cells material inside. This head-immobilizer is unique thanks to the internal aluminium layer that allows to shape the device to the patient's head. The fastening takes place by means of a VELCRO® system and a strap system to stabilize chin and forehead.



Art. 628

## GREEN

**Universal head immobilizer**

Essential for head immobilization, this device is easily adjustable to spinal boards. Made in close-cells expanded material and coated in exclusive vinyl.



Art. 627

## KIDS

**Pediatric padded head immobilizer**

Padded head immobilizer for pediatric spinal boards, made in close-cells expanded material, coated in vinyl.



Pediatric





Art. 621

## X-DOUBLE

**Adjustable head immobilizer for ERGON Scoop Stretcher**

The X-DOUBLE is a two-piece detachable head immobilizer. The strap system of this immobilizer allows to make full use of the ERGON's Flexilock while immobilizing the patient's head for transport. Made in expanded close-cells material with an internal aluminium layer. Painted in non-absorbent and antibacterial finishing. Designed to be used with safety helmets.



Trauma



Art. 9150

## X-LOCK

**Ergonomic head immobilizer for ERGON stretcher**

Innovative head immobilizer, made up of a single element, integrated straps, ensuring a fast fastening of the head to the stretcher. Made in laminated polyester and bi-density expanded polyurethane.



Trauma



Art. 670-001

Orange/black head immobilizer for Extra Stair Chair



Trauma



Art. 623

## DOUBLE PLUS

**Adjustable head immobilizer for scoop stretchers**

DOUBLE PLUS is the first detachable and adjustable two-piece head immobilizer. The fastening system with straps allows the positioning of the immobilizer and the separation of the two halves of the stretcher, in order to safely load and transport the patient. Made in close-cells expanded material with an internal aluminium layer, painted in non-absorbent and anti-bacterial finish.



Trauma



Art. 624

## DOUBLE

**Head Immobilizer for scoop stretcher**

This head immobilizer shares the same characteristics of the DOUBLE PLUS as far as detachment and adjustability are concerned. Its uniqueness stands in the head immobilization system: instead of the mouldable aluminium layer, this model is supplied with the traditional lateral blocks. Supplied with padded straps, this device is made in close-cells expanded materials and painted in non-absorbent and anti-bacterial finish.



Trauma



IMMOBILIZATION

# SPINAL BOARDS

with pins

Art. 9040  
Yellow

Art. 9042  
Orange

Accessories



Art. 574

2-pieces blue belt with carabiners

Art. 9040 - Art. 9042

## GALAXY-P

Yellow and orange spinal board with pins

Spinal board with the same features as "Galaxy" Art. 9010, 9012 and equipped with 10 plastic pins for a quick-fastening of restrain systems (with carabiners). Supplied with 3 belts (Art. 605/MEB).

Art.	Colour	Lenght	Width	Thickness	Weight	Load
9040	Yellow	183 cm	41 cm	5 cm	5,8 kg	150 kg
9042	Orange					

Immobilization

118

Art. 9050

## NILIX

Yellow standard spinal board with pins

Basic spinal board in high density plastic material, equipped with 8 pins.

**Belts are not included.**

Art.	Lenght	Width	Thickness	Weight	Load
9050	183 cm	40 cm	5 cm	6,5 kg	150 kg





IMMOBILIZATION

# VACUUM MATTRESSES



Trauma

Art. 9104

## NAXOS

Single chamber vacuum mattress

Single chamber vacuum mattress, blue and orange.

Weight 6 kg - Size cm 200x80

Accessories



Art. 9104-001

Carrying bag for "NAXOS",  
single-chamber vacuum  
mattress



Outer  
view

Immobilization

119

Product Catalogue 2023

IMMOBILIZATION

# VACUUM MATTRESSES



Art. 894-001

Intake valve  
for vacuum  
mattresses  
and splints.



Art. 890-001

Carrying case for  
vacuum mattress.  
Size: 93x63x18 cm



Trauma

Art. 894

## SNAKE

Independent chambers vacuum mattress

Vacuum mattress with 14 independent chambers, for a quick and precise immobilization. Extremely light and resistant, this model features 8 ergonomic handles and restraining belts. Ensuring perfect containment, the "Snake" is a truly unique device. Working conditions between -30°C and +70°C. A lasting vacuum and no deformation, together with the UNI EN 1865 certification, top off an excellent device. Compatible with X-Rays, CT and MRI.

**Weight 8 kg - Size cm 210x100**





Trauma

Art. 895

## SNAKE PLUS

Independent chambers vacuum mattress with head immobilizer

Mattress with integrated head immobilizer for a quick head immobilization. Made of polyurethane, treated with flame retardant, certified according to the European harmonized safety standards UNI EN 1865. 14 independent chambers for a safe and precise immobilization. Working conditions between -30°C and +70°C. Extremely light and sturdy, equipped with 8 ergonomic padded handles and safety straps for an excellent patient containment. Compatible with X-Rays, CT and MRI.

**Weight 8,3 kg - Size cm 210x100**



Accessories



Art. 888

### Aluminium suction pump

Aluminium suction pump for vacuum mattresses with foot support.



Art. 896/DF

### ABS suction pump

ABS suction pump for vacuum mattresses, for inflation or deflation. Equipped with junctions.



Immobilization

121

IMMOBILIZATION

# VACUUM MATTRESSES



Art. 894-001

Intake valve  
for vacuum  
mattresses  
and splints.



Art. 9100-001

Carrying case for  
pediatric vacuum  
mattress.  
Size: 72x52x16 cm



Pediatric



Trauma

122

Art. 9100


## BUNNY

**Pediatric vacuum mattress with independent chambers**

Polyurethane mattress with flame retardant treatment, 14 independent chambers for a safe and precise immobilization. Working conditions between -30°C and +70°C. Extremely light and sturdy, equipped with 4 ergonomic padded handles and safety straps for an excellent patient containment. Compatible with X-Rays, CT and MRI.

**Weight 3 kg - Size cm 120x70**





Pediatric Trauma

Art. 9102

## BUNNY PLUS

**Pediatric vacuum mattress with independent chambers and head immobilizer**

Polyurethane mattress with flame retardant treatment, 14 independent chambers for a safe and precise immobilization. Working conditions between -30°C and +70°C. Extremely light and sturdy, equipped with 4 ergonomic padded handles and safety straps for an excellent patient containment. Compatible with X-Rays, CT and MRI.

**Weight 3,5 kg - Size cm 120x70**



### Accessories

Art. 888

#### Aluminium suction pump

Aluminium suction pump for vacuum mattresses with foot support.



Art. 896/DF

#### ABS suction pump

ABS suction pump for vacuum mattresses, for inflation or deflation. Equipped with junctions.



IMMOBILIZATION

# VACUUM MATTRESSES

Immobilization

124



Trauma

Art. 894/SNOW

## SNAKE SNOW

Independent chambers vacuum mattress with anti-sliding features

Vacuum mattress with 14 independent chambers, for a quick and precise immobilization. Equipped with anti-sliding features for rescue in snow-capped conditions. Extremely light and resistant, this model features 8 ergonomic handles and restraining belts. Ensuring perfect containment, the "Snake" is a truly unique device. Working conditions between -30°C and +70°C. A lasting vacuum and no deformation, together with the UNI EN 1865 certification, top off an excellent device. Compatible with X-Rays, CT and MRI. **Size: cm 210x100. Weight: 8,9 kg**

Accessories

Art. 888

### Aluminium suction pump

Aluminium suction pump for vacuum mattresses with foot support.



Art. 896/DF

### ABS suction pump

ABS suction pump for vacuum mattresses, for inflation or deflation. Equipped with junctions.





Trauma

**Art. 891**

## FLAKE

**Independent chambers vacuum mattress**

3-chambers vacuum mattress, for a quick and precise immobilization. Extremely light and resistant, this model features 8 ergonomic handles and restraining belts. Ensuring perfect containment, the "FLAKE" is a truly unique device. Long lasting vacuum and no deformation complete this model. **Weight 6,7 kg - Size cm 200x90**



### Accessories

**Art. 888**

#### Aluminium suction pump

Aluminium suction pump for vacuum mattresses with foot support.



**Art. 896/DF**

#### ABS suction pump

ABS suction pump for vacuum mattresses, for inflation or deflation. Equipped with junctions.



**Art. 894-001**

**Intake valve for vacuum mattresses and splints.**

Immobilization

125



# BABY FIX

PEDIATRIC RESTRAINT SYSTEM WITH VACUUM MATTRESS

Art. 9130

## MEBER BABY FIX

**Restraint system with straps for ambulance transport**

MeBer Baby Fix is a universal pediatric transport device, adjustable to any stretcher, with 6 straps. Suitable for infants and children up to 130 cm and 35 KG, this device ensures a safe transport thanks to its 5 anchor points. Supplied with case Art. 12034. Certified UNI EN 1789.

**Length: 74 cm; Head Width: 39 cm; Feet width: 20 cm**



Pediatric



Art. 9128

## MEBER VACUUM BABY FIX

**Small vacuum mattress for pediatric restraint system with straps (ambulance transport)**

MeBer Vacuum Baby Fix is a vacuum mattress, to be combined with Art. 9130. Extremely light and resistant. Thanks to its safety straps and to a precise and quick containment, this 9128 is a truly unique device. Certified UNI EN 1789, radiolucent, a-magnetic.

**Size: 79 x 45 cm**



Pediatric



Trauma

MeBer  
Exclusive



Accessories



Pediatric

Art. 12034

**Red Case for BABY FIX and VACUUM BABY FIX**



# RIGID SPLINTS



**Art. 430**  
Rigid splint set with  
carrying case



**Art. 425**  
Adult rigid splint for  
the lower limbs



Pediatric

**Art. 427**  
Rigid splint for upper limbs (adults) /  
lower limbs (pediatric patients)



Trauma

**Art. 425 - Art. 427 - Art. 430**

## ARCOS

### Rigid splints

Rigid splints in fire-proof and water-resistant PVC. This device ensures a correct limb immobilization.

Equipped with 4 VELCRO straps (different lengths), ensuring a safe and quick application.

This model is shaped for a firm lacing to ankles and upper leg and is equipped with integrated foot support. Radiolucent.

#### Available sizes:


**Adult arm / pediatric leg: length 52 cm – base 9 cm**

**Adult leg: length 74 cm – base 13 cm**

# SPLINTS



Art. 880



Art. 884



Art. 882



Art. 886

Accessories

Art. 889

Small aluminium suction pump



Art. 894-001

Intake valve for vacuum mattresses and splints.



Art. 880 - Art. 882 - Art. 884 - Art. 886



Trauma

## HALLEY

### Vacuum Splints

Made in high quality materials, the HALLEY ensures a perfect immobilization of the lower and upper limbs by means of vacuum suction. This device is equipped with an external intake valve to avoid annoyance to the patient. This device can be ordered as an individual unit or as a set. The set includes one arm vacuum splint, one leg vacuum splint, one forearm vacuum splint, a small suction pump and a nylon carrying case. All items are radiolucent.

Art. 880 Leg vacuum splint

Art. 882 Arm vacuum splint

Art. 884 Forearm vacuum splint

Art. 886 3 splints set with carrying case and suction pump

# SPLINTS



Art. 9160

## VACUUM SPLINT

**Leg vacuum splint with 16 independent chambers**

For a quick and precise immobilization. Equipped with internal valve and junctions.

**Compatible with X-Rays, CT and MRI.**



Art. 9162

## VACUUM SPLINT

**Arm vacuum splint with 14 independent chambers**

For a quick and precise immobilization. Equipped with internal valve and junctions.

**Compatible with X-Rays, CT and MRI.**

Art. 9140 - Art. 9142 - Art. 9144

## TRACTION SPLINT



Trauma

**Traction splint**

**Art. 9140 - Traction splint set for adult and pediatric application – with case**

Splint set for femoral traction, supplied with case. Equipped with 2 immobilization devices, adult and pediatric.

**Art. 9142 - Adult traction splint with case**


The device ensures a correct and aligned traction, thus reducing pain and preventing further injury during transport.

**Art. 9144 - Pediatric traction splint with case**

The device ensures a correct and aligned traction, thus reducing pain and preventing further injury during transport.



Pediatric



# SPLINTS



Art. 417



320 mm

Art. 412



400 mm

Art. 410



560 mm

Art. 415



610 mm

Art. 407



600 mm

Art. 405

Art. 405 - Art. 407 - Art. 410 - Art. 412 - Art. 415 - Art. 417

## NEOFIX

### Neoprene splint

Rigid flexible neoprene splints, equipped with an inner aluminium sheet to ensure perfect adaptation on a case by case basis. The aluminium sheet can be removed for cleaning. The neoprene material ensures constant body temperature and perfect shaping of the splint, as well as a being easily cleanable.

There are 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.



Trauma

Art. 405 Leg splint

Art. 407 Arm splint

Art. 410 Forearm splint

Art. 412 Wrist splint

Art. 415 Ankle/Elbow splint

Art. 417 Splint set (5 units) with carrying case



Trauma

Art. 722

## AIR SPLINT

### Set of 6 pneumatic splints with pump and carrying case

Plastic pneumatic splints including: carrying case, pump, hand/wrist splint, arm splint, forearm splint, foot splint and short leg splint, long leg splint





IMMOBILIZATION

# CERVICAL COLLARS



Trauma

Art. 836

## PATRIOT

Adjustable cervical collar



Trauma



Pediatric

Art. 837

## PATRIOT BABY

Adjustable pediatric cervical collar



Sharing the same features of its adult version, the PATRIOT BABY is an adjustable pediatric cervical collar with 3 different heights.



## IMMOBILIZATION

# CERVICAL COLLARS



## GEMINI

### Rigid 2 pcs cervical collar

The GEMINI cervical collar is a highly professional and indispensable device for rescue operators. Composed by 2 pieces, it fits in the front and rear parts of the neck, ensuring a perfect rachis immobilization. Equipped with a central hole for tracheotomy, it is easily adjustable thanks to the straps. Made in printed polyethylene, the GEMINI comes in various colourings to properly choose size. Internal covering in washable E.V.A., fixed with plastic rivets. Available in 6 different sizes.

Art.	Size	Set	Pieces per pack	
9213	Set 3 sizes	9204 + 9206 + 9208	3	With case
9214	Set 4 sizes	9204 + 9206 + 9208 + 9210	4	With case
9216	Set 6 sizes	9200 + 9202 + 9204 + 9206 + 9208 + 9210	6	With case



## ARTIS

### Rigid cervical collar in 1 piece

1 piece cervical collar equipped with a central hole for tracheotomy, it is easily adjustable thanks to the straps. Made in printed polyethylene, the ARTIS comes in various colourings to properly choose size. Internal covering in washable E.V.A., fixed with plastic rivets. Available in 6 different sizes.

Art.	Size	Set	Pieces per pack	
9233	Set 3 sizes	9224 + 9226 + 9228	3	With case
9234	Set 4 sizes	9220 + 9224 + 9226 + 9228	4	With case
9236	Set 6 sizes	9220 + 9222 + 9224 + 9226 + 9228 + 9230	6	With case



**Art. 842**  
**PEDIATRIC**



Pediatric

**Art. 840**  
**INFANT**



Pediatric

**Art. 846**  
**STANDARD**

**Art. 848**  
**MEDIUM**

**Art. 844**  
**STOUT**

**Art. 850**  
**TALL**

## ANTARES

### Rigid cervical collar in 2 pieces

The ANTARES cervical collar is a highly professional and indispensable device for rescue operators. Composed by 2 pieces, it fits in the front and rear parts of the neck, ensuring a perfect rachis immobilization. Equipped with a central hole for tracheotomy, it is easily adjustable thanks to the straps. Made in transparent polyethylene, the ANTARES comes in various colourings to properly choose size. Internal covering in washable E.V.A., fixed with plastic rivets. Available in 6 different sizes.

Art.	Size	Set	Pieces per pack	
853	Set 3 sizes	844 + 846 + 848	3	With case
854	Set 4 sizes	842 + 844 + 846 + 848	4	With case
856	Set 6 sizes	840 + 842 + 844 + 846 + 848 + 850	6	With case

## SIRIUS

### Rigid cervical collar in 1 piece

Collare cervicale rigido ad 1 pz realizzato in polietilene trasparente e rivestimento interno in E.V.A. lavabile fissato con rivetti in plastica. Dotato di chiusure regolabili a strappo colorate per contraddistinguere le varie misure, è provvisto di foro centrale per la tracheotomia. Disponibile in 6 diverse misure.



**Art. 862**  
**PEDIATRIC**



Pediatric

**Art. 860**  
**INFANT**



Pediatric

**Art. 866**  
**STANDARD**

**Art. 868**  
**MEDIUM**

**Art. 864**  
**STOUT**

**Art. 870**  
**TALL**

Art.	Size	Set	Pieces per pack	
873	Set 3 sizes	864 + 866 + 868	3	With case
874	Set 4 sizes	862 + 864 + 866 + 868	4	With case
876	Set 6 sizes	860 + 862 + 864 + 866 + 868 + 870	6	With case

IMMOBILIZATION

## EXTRICATION



Trauma



Immobilization

134

Art. 899

**MED****Spinal immobilizer**

The MED spinal immobilizer is ideal for the immobilization and extrication of a patient under critical conditions, such as road accidents. Thanks to its enveloping shape and structure, it ensures blocking of vertical and lateral movements of the head and spinal column, to minimize potential trauma during transfer. Equipped with coloured straps and belts, a padded cushion and two straps for forehead and chin. Supplied with carrying case. **Weight: 2,9 kg**





Trauma




Art. 9124/B

## SPINAL IMMOBILIZER

Standard blue spinal immobilizer

Essential device for a safe patient immobilization and extrication. Made of vinyl-coated nylon. Equipped with padded cushion, 2 straps and a carrying bag. **Weight: 2,9 kg**



Trauma



Pediatric

Art. 9120

## MED PED

Pediatric spinal immobilizer

The "MED PED" spinal immobilizer is ideal for the immobilization and extrication of a patient under critical conditions, such as road accidents. Thanks to its enveloping shape and structure, it ensures blocking of vertical and lateral movements of the head and spinal column, to minimize potential trauma during transfer. Equipped with coloured straps and belts, a padded cushion and two straps for forehead and chin. Supplied with carrying case. **Weight: 1,4 kg**



# DIAGNOSTIC

6

The diagnostics line from MeBer includes traditional devices such as phonendoscopes and sphygmomanometers as well as multiparameter monitors, pulse oximeters and TecnoHeart Automated Defibrillators. Beyond the obvious relevance in terms of vital signs' monitoring or, when it comes to defibrillators, the life-saving role they play, the diagnostics line is a crucial aspect for the patient's perception of an operator's degree of professionalism

Diagnostics

136

Product Catalogue 2023



Sphygmomanometers



Stethophonendoscopes



Defibrillators



Pulse Oximeters

S



Discover MeBer products on-line. Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website. You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/diagnostic.aspx>

# SPHYGMOMANOMETERS

Art. 8006

## SPHYGMOMANOMETER

Aneroid sphygmomanometer

Manometer, range 0-300 mm Hg. Accuracy  $\pm 3$  mm Hg.  
Adult Velcro® cuff and metal ring. Supplied with carrying case.



Art. 8008

## SPHYGMOMANOMETER

One hand aneroid sphygmomanometer

Antishock ABS case manometer Ø 60 mm. Range 0-300 mm Hg. Accuracy  $\pm 3$  mm Hg. Manometer with NON STOP PIN. Easy left-hand or right-hand operation. Adult nylon cuff with Velcro® strap and metal ring. Single corrugated latex ball for improved intake and chrome-plated vent valve. Supplied with carrying case with zip.



Art. 8012

## SPHYGMOMANOMETER

Red aneroid one hand sphygmomanometer (Plus)

Ø 48 mm. Precision sphygmomanometer with protective ring and extremely sensitive deflation valve. Anodised aluminium.







Art. 8014

## SPHYGMOMANOMETER

**Compact aneroid sphygmomanometer**

This handmade sphygmomanometer features a high quality finish. It features state-of-the-art components and the highest technological solutions. The German movement follows the highest quality standards.

The Polycarbonate glass is alcohol- and stress-proof. The drain valve is easy to activate and extremely sensitive. Equipped with natural rubber bladder. The nylon cuff is resistant and anti-allergic. The manometer, fixed on the cuff, is highly resistant to falls or casual blows. Made in Italy according to the EN 1060-1 and EN 1060-2.



Art. 8004

## SPHYGMOMANOMETER

**Coloured aneroid sphygmomanometer with case**

Range 0-300 mm Hg. Accuracy  $\pm 3$  mm Hg. Manometer with NON-STOP PIN.

Adult nylon cuff with Velcro® strap and metal ring. Single corrugated latex ball for improved intake and chrome-plated vent valve. Supplied with carrying case with zip.



Pediatric

Art. 8009

## SPHYGMOMANOMETER

**Shock proof aneroid sphygmomanometer for children**

Manometer range 0-300 mm/Hg. Accuracy  $\pm 3$  mm/Hg. Light and sturdy. Pediatric nylon cuff with Velcro® strap and metal ring. Supplied with carrying case with zip.



Pediatric

Art. 8011

## SPHYGMOMANOMETER

**Neonatal and pediatric sphygmomanometer SET**

Sphygmomanometer set equipped with 1 neonatal cuff and 1 pediatric cuff, supplied with coloured carrying case.

# WALL-MOUNT SPHYGMOMANOMETERS



Art. 1740

## SPHYGMOMANOMETER

### Rotating wall-mount aneroid sphygmomanometer

Round sphygmomanometer. Easy to read, it is ideal for ambulances and clinics. 65° rotation in both directions, roomy cuff basket for storage of accessories, coiled tube extendible to a length of 3 m, single latex ball for improved intake and chrome-plated vent valve, dial-plate diam. 147 mm. Specially tempered copper-beryllium diaphragm. Calibrated nylon cuff with Velcro® strap and 2-tube latex bladder.

#### Accessories



Pediatric

Art. 1740-001

Pediatric cuff for wall-mount aneroid sphygmomanometer  
Art. 1740

### Adult cuffs



**Art. 8140**  
1 tube

**Art. 8141**  
2 tubes

### Pediatric cuffs



**Art. 8142**  
1 tube

**Art. 8143**  
2 tubes



Pediatric



Pediatric

DIAGNOSTICS

# AUTOMATIC BLOOD PRESSURE MONITOR



Art. 8534

## AUTOMATIC BASIC

### Automatic blood pressure monitor with adult cuff

Easy to use, accurate and reliable medical device for blood pressure monitoring. This device complies with the requirements of the Medical Device Directive 93/42/EEC and offers different cuff sizes, depending on the circumference of the patient's arm. Optional «Easy» preformed cuffs are also available for a better fitting and comfort.

Dimensions	135,5 x 82 x 57 mm
Weight	340 g
Storage Temperature	-20/+55 °C - -4/+131 °F 15-95% relative max. humidity
Operating temperature	10 - 40 °C / 50 - 104 °F 15-95% relative max. humidity
Measurement range	20 - 280 mmHg – blood pressure / 40 - 200 beats per minute – pulse
Memories	last 30 measurements
Voltage source	4 x batteries 1,5 Volt; size AA Mains adapter DC 6V, 600 mA (optional)

#### Art. 8534-001

S-SIZE Pediatric cuff  
17-22 cm.



Pediatric



#### Art. 8534-002

M-L SIZE Standard cuff  
22-42 cm




#### Art. 8534-003

L-XL SIZE Obese cuff  
32-52 cm



Bariatric



#### Art. 8536

220V power cable for "AUTOMATIC BASIC"  
blood pressure monitor

Accessories



Diagnostics

141

Product Catalogue 2023

# PHONENDOSCOPES, STETHOPHONENDOSCOPES and **INFRARED THERMOMETER**




Art. 8030/B

## RAPPAPORT

### Blue stethophonendoscope

Professional stethophonendoscope. Thanks to its accessories, it is possible to perform different kinds of auscultation. Equipped with twin tube, metal head, rotating pavilions for the its double function, 3 cuffs in different sizes, 2 smooth olives, 2 rigid olives and 2 membranes.




Art. 8020

## PHONENDOSCOPE

### Flat phonendoscope

Flat phonendoscope with anodized stamped and lathed aluminium head diameter 46 mm. Threaded metal ring for membrane blocking. Metal chromium plated binaural, length 150 mm. Y rubber tube in one piece, diameter 5 mm. length 60 cm. Single carton box. Available in 3 different colours.



Art. 8025

## STETHOPHONENDOSCOPE

### Flat stethophonendoscope

Stethophonendoscope with 47/31 mm head diameter. Aluminium anodized head, Chromium plated metal binaural, Y rubber tube in one piece. Rotating piper-fitting for both functions. Available in 3 different colours.

Art. 8377

## INFRARED

### Infrared thermometer

Infrared thermometer for a contact-less, quick, precise and reliable body temperature measurement. Easy-to-use, intuitive and hygienical, this device ensures stable measurements, free from any interference. Certified 93/42/CEE-DM Class IIa



Temperature range	Body temperature:
	32÷43,9 °C / 89,6÷109,4 °F
	Surface temperature:
	0÷60 °C / 32÷140 °F
Terms of use	Room temperature:
	0÷40 °C / 32÷100 °F

Terms of use	Environment
	10÷40 °C / 50÷104 °F
Distance	Relative humidity ≤ 85%

Distance	3 ÷ 5 cm (1,2 ÷ 2 pollici)
----------	----------------------------

Response time	1 second
---------------	----------

Memory	last 32 measurements
--------	----------------------

Battery	2 x 1,5 V AAA type
---------	--------------------

Size	100 x 40 x 155 mm
------	-------------------

Weight	80 g (without battery)
--------	------------------------



# MONITOR

# TG



**NEW**  
PRODUCT

**Art. 8540**

## MONITOR TG

### Multiparameter monitor

The TG Monitor is a professional device, equipped with accessories and circuits with a double function: monitoring and diagnosis. This model ensures the collection of vital signs, which can be monitored individually or parallelly.

Features:

- Multiparameter monitoring, with value and tracing
- Battery or current powered – 100 : 240V ca, 50 : 60Hz power source. The battery ensures continuity of use in case of missing power source (up to 4 hours with its standard battery)
- Integrated printer
- Remote controller
- Parameter adjustment and setup
- 10 inch monitor

Supplied  
with



#### Art. 8540-001

SpO2 sensor for adults, re-usable



#### Art. 8540-002

Adult cuff 27-35 cm



#### Art. 8540-003

**NEW**  
PRODUCT

NIBP pediatric cuff 16-21,5 cm



#### Art. 8540-004

**NEW**  
PRODUCT

NIBP neonatal cuff 10-15 cm



#### Art. 8540-005

**NEW**  
PRODUCT

Rechargeable lithium battery for TG Monitor



#### Art. 8540-006

**NEW**  
PRODUCT

Thermal graph paper



Accessories



Diagnostics

5 parameters:

- ✓ Ecg - 5 channels / 2 tracks
- ✓ Systolic / diastolic / mean pressure
- ✓ Heartbeat
- ✓ Pulse oximeter
- ✓ Breath monitoring
- ✓ Double probe temperature

# TECNOHEART PLUS AED

SEMI-AUTOMATIC DEFIBRILLATOR

Dimensions	240 x 294 x 95 mm
Weight	2,65 kg
Length of electrodes	180 cm
Energy level	Adult 185-200 J ( $\pm 5\%$ ) Pediatric 45-50 J ( $\pm 5\%$ )
Battery life	200 shocks / 10 h monitoring
Battery	LiMnO <sub>2</sub> (15V-4200mAh)



Art. 8440

## TECNOHEART PLUS AED

Semi-automatic defibrillator


I defibrillatori AED oggi in commercio sono prevalentemente semiautomatici, pertanto vengono pilotati da comandi vocali e l'attivazione di scariche avviene con comando manuale azionato dall'operatore. Un'azione semplice e guidata.

### Istruzioni vocali guidate

Seguendo i messaggi di testo e le istruzioni vocali, chiunque è in grado di utilizzare questo dispositivo salvavita.

Tecnoheart Plus main features

- ✓ LiMnO<sub>2</sub> long life battery: 5 years, 200 discharges
- ✓ Critical components self-test: to keep its efficiency
- ✓ Impedance calculation: to verify the integrity of contact between the AED and the patient
- ✓ Heart on software: useful to archive data concerning AED usage and for firmware updates
- ✓ Pre-connected electrodes: to reduce time interventions
- ✓ IRDA port: to allow pc-defibrillator communication
- ✓ Variable energy: adult (<25kg) from 185 to 200 J/child (<25kg) from 45 to 50 J.




Tecnoheart Plus: technology

- ✓ Biphasic technology, the most advanced system for shocks delivery
- ✓ First shock efficacy: 100%, certified
- ✓ Mean shock per patient: 2.3
- ✓ Time required for an effective defibrillation: 55 secs

# TECNOHEART TRAINER

SEMI-AUTOMATIC DEFIBRILLATOR

Dimensions	314 x 259 x 109 mm
Weight	1,9 kg
Power	DC 10,8V / 2200mAh DV 15V / adaptor
Static power	< 5µA
Maximum power	< 500mA
Battery	Rechargeable 3-cell lithium battery



Art. 8442

## TECNOHEART PLUS TRAINER

### Semi-automatic defibrillator

The DAE training system is a basic first-response training simulation, for a correct application of the defibrillation procedures by the operators.

#### The DAE Trainer:

- ✓ simulates the provision of electric shocks and ensures a safe use during the training, since it doesn't deliver high voltages,
- ✓ includes 10 preconfigured scenarios of training that simulate realistic episodes of sudden cardiac arrest
- ✓ can be controlled with an infrared remote control that allows scenarios modification during the training so you can test the students' reaction.
- ✓ thanks to its remote control, you can also simulate the shockable rhythm and the non-shockable rhythm.

The PC software allows the configuration of custom scenarios. To connect the device to the PC serial port, you need a standard serial cable D-sub 1: 1 9-pin.

#### Operational conditions

- ✓ Operation: from 0 to 43 °C (from 32 to 109.4 °F)
- ✓ Relative humidity: from 5 to 95% (non-condensing)
- ✓ Altitude: from 0 to 4.475 m

#### Storage conditions:

- ✓ Storage: from -20 to 60 °C (from -4 to 140 °F)
- ✓ Relative humidity: from 5 to 95% (non-condensing)
- ✓ Altitude: from 0 to 12.192 m

#### Remote control:

- ✓ Remote control: DC 3V / Button batteries
- ✓ Static power: < 5uA
- ✓ Maximum power: < 10mA

#### List of components:

- ✓ N° 1 Remote control
- ✓ N° 1 Training electrode for adults (reusable up to 100 times)
- ✓ N° 1 Operator's manual
- ✓ N° 1 Carrying case (Nylon)
- ✓ N° 1 DC adapter
- ✓ N° 1 CD Software DAE Trainer (Optional)

Diagnostics

145

## DIAGNOSTICS

# DEFIBRILLATOR DISPLAY CABINET

Art. 8398

## THERMOREGULATED DISPLAY CABINET


### Thermoregulated defibrillator display cabinet

Made of P02 carbon steel with a thickness of 15/10, watertight welds, rabbit gaskets, water- and dust-proof door with IP44 degree of protection and snap closure with pull adjustment and potential padlock or seal addition (not included). Hinges in durable, tear-resistant, corrosion-resistant metal alloy. Shelf height-adjustable interior with grating for hot air recirculation. Comfortable techno-polymer handle for easy opening. Internal liquid crystals temperature detection, positioned at the front.

Electric resistance thermoregulation system with thermal protection 80 °C, 200W-230W power, room thermostat + 5 °C / + 20 °C, insulated with cable.

**NOTE: THIS CABINET IS NOT SUITABLE FOR MARITIME LOCATIONS**

Outer size	425 x 480 x 215 mm (L x H x P)
Internal size	420 mm (L) - 390 mm (H)
Door light	365 mm (L) - 365 mm (H)
Weight	14,50 kg



Thermostat display



Art. 8396

## ALARMED THERMOREGULATED CABINET

Alarmed thermoregulated cabinet for defibrillator, equipped with power unit

Main features as per Art. 8398

Equipped with anti-theft magnetic alarm, this cabinet emits a loud noise when drawn without permission. Siren powered by 3 "LR44" batteries. Alarm signal volume 100 dB.

**NOTE: THIS CABINET IS NOT SUITABLE FOR MARITIME LOCATIONS**





Art. 8432

## STANDARD CABINET

White standard cabinet for  
AED defibrillators



Built in galvanized P02 carbon steel sheet; white RAL9010 epoxy-polyester powder coating; processed in 180° C thermal tunnel. 1 window door, hinges with galvanized pin, opening 180°; lockable seal closure in nickel-plated zamak; pre-drilled for attachment to the wall and on a special floor stand (Art. 8432-001).

Outer size	452 x 425 x 160 mm (L x H x P)
------------	--------------------------------

Weight	6 kg
--------	------

Art. 8432-001

## FLOOR STAND FOR ART. 8432

Floor stand for white  
standard cabinet Art. 8432

Floor stand for white standard cabinet Art. 8432, equipped with holes for floor fastening and interlocking + pre-drilled system to facilitate any wiring and cable passages in safety. Built in galvanized carbon steel sheet P02. Powder coated epoxy-polyester RAL9010 white, processed in thermal tunnel at 180 °C.



Size	250 x 1000 x 150 mm
------	---------------------

Base size	300 x 150 mm
-----------	--------------

Art. 8396-001

## FLOOR STAND FOR ART. 8396 AND 8398

Floor stand for yellow cabinet  
Art. 8396 and 8398

Floor stand for yellow defibrillator cabinet (Art. 8396 and 8398), equipped with holes for floor fastening. Built in galvanized carbon steel sheet P02 7/10. Powder coated epoxy-polyester.



Size	270 x 800 x 120 mm
------	--------------------

Base size	300 x 150 mm
-----------	--------------

# VITAL TEST

HANDHELD PULSE OXIMETER



Art. 11200

## VITAL TEST

Handheld pulse oximeter with adult sensor

148

Vital Test is a simple and professional device: it can be used in intensive care units, surgeries, dentist's surgeries and emergency rooms, ensuring immediate, accurate and non-invasive monitoring. Vital Test has a pleasant and

ergonomic shape, it is easy to grip and to be transported. Vital Test allows SpO2 measurement, heart-rate detection, waveform visualization and also monitor the last 10 minutes trend value of these two parameters. Supplied with rechargeable battery.

Supplied  
with



Art. 11206

Adult sensor for  
pulse oximeter  
(Vital Test)



Art. 11210

Vital Test carrying  
case



### Accessories



Art. 11212

Neonatal sensor  
for Vital Test Pulse  
oximeter



Pediatric



Art. 11208

Pediatric sensor  
for Vital Test Pulse  
oximeter



Pediatric



Art. 11204

Wall-mount Pulse  
oximeter support

# OXYTEST and SATFINGER II

PORTABLE PULSE OXIMETERS

Supplied  
with



**Art. 11226**

Adult sensor for  
OXYTEST pulse  
oximeter



**Art. 11220**

## OXYTEST

Portable oximeter equipped with adult sensor, USB cable and rechargeable lithium battery.

The Oxytest is a light, easy-to-use pulse oximeter, perfect for transport. Thanks to its sensors Art. 11224 and Art. 11222 (not included), the Oxytest can fit pediatric patients and new-borns. It is powered by means of a rechargeable lithium battery.

- ✓ OLED display, 4 directions to monitor arterial oxygen saturation, pulse rate and pulse strength
- ✓ Battery voltage indicator
- ✓ Standby or automatic sleep mode

Accessories



**Art. 11228**

Temperature  
detector for  
OXYTEST pulse  
oximeter



Pediatric



**Art. 11224**

Pediatric sensor  
for OXYTEST  
pulse oximeter



Pediatric



**Art. 11222**

Neonatal  
sensor for  
OXYTEST pulse  
oximeter

**Art. 8336**

## SATFINGER II

Finger pulse oximeter

- ✓ Easy-to-use
- ✓ Portable
- ✓ CE certified
- ✓ For home and professional use



Size (LxHxD)	29x35x69 mm
Weight	60 g

# RESUSCITATION

Our resuscitation line includes all the devices aimed at the restoration and support of the respiratory system. The range covers laryngoscopes with optical fibre, masks, cannulas, resuscitation bags and other complementary items related to oxygen administration

Resuscitation

150

Product Catalogue 2023




Laryngoscopes



Resuscitation kits



Resuscitation masks



Resuscitation bags



Oxygen reservoirs and circuits



Aerosol Therapy masks



Oxygen tank



Pressure regulators



# ION



Discover MeBer products on-line.  
Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website.  
You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/reaanimation.aspx>

Resuscitation

151

Product Catalogue 2023



# STANDARD LIGHT LARYNGOSCOPES

## MILLER

Miller standard light laryngoscope in stainless steel, available in 3 or 4 blades set (handle included, Art. 1650)

### Art. 1634

Standard light laryngoscope blade n. 4

### Art. 1633

Standard light laryngoscope blade n. 3

### Art. 1632

Standard light laryngoscope blade n. 2

### Art. 1631

Standard light laryngoscope blade n. 1

### Art. 1630

Standard light laryngoscope blade n. 0



Pediatric

## McINTOSH

McIntosh standard light laryngoscope in stainless steel, available in 3 or 4 blades set (handle included Art. 1650)

### Art. 1615

Standard light laryngoscope blade n. 5

### Art. 1614

Standard light laryngoscope blade n. 4

### Art. 1613

Standard light laryngoscope blade n. 3

### Art. 1612

Standard light laryngoscope blade n. 2

### Art. 1611

Standard light laryngoscope blade n. 1

### Art. 1610

Standard light laryngoscope blade n. 0



Pediatric

### Art. 1638

**Laryngoscope set with 3 standard light blades**  
(Art. 1633 - 1632 - 1631) **with handle and case**

### Art. 1640

**Laryngoscope set with 4 standard light blades**  
(Art. 1633 - 1632 - 1631 - 1630) **with handle and case**

### Art. 1618

**Laryngoscope set with 3 standard light blades**  
(Art. 1614 - 1613 - 1612) **with handle and case**

### Art. 1619

**Laryngoscope set with 3 standard light blades**  
(Art. 1612 - 1611 - 1610) **with handle and case**

### Art. 1620

**Laryngoscope set with 4 standard light blades**  
(Art. 1614 - 1613 - 1612 - 1611) **with handle and case**



**Art. 1620-001**  
**Laryngoscope carrying case**

Laryngoscope case.  
Fit for 3-4 blades kit with handle.



**Art. 1645**

Pediatric handle for standard light laryngoscope



Pediatric



**Art. 1650**

Adult handle for standard light laryngoscope



**Art. 1655**

Spare bulb

RESUSCITATION

# RESUSCITATION KIT



Art. 2000/F | Art. 2000/D

## RIASET

Resuscitating set with carrying case

Resuscitating set in ABS cover with fixing support complete with 0,5 l rechargeable steel oxygen bottle, fixed calibration reduction gear, reanimator bag, mask (size 4), 3 Guedel canulas (sizes 1-2-3), tongue puller, helicoidal mouth opener, tube.

Size: mm 400x260x130



Art. 2100/F | Art. 2100/D | Art. 2100/GB

## RIABOX

Resuscitating set with carrying case

Resuscitating set in ABS cover with fixing support complete with 0,5 l rechargeable steel oxygen bottle, reducer with manometer and dosing unit, reanimator bag, 3 masks (sizes 3-4-5), 3 Guedel canulas (sizes 1-2-3), tongue puller, helicoidal mouth opener, tube

Size: mm 445x320x150

Art. 1090

## GAG


Helicoidal gag

Helicoidal gag in polyethylene with non-slip knurled knob for a secure grip





# CANNULAS AND MASKS



Art. 1409 - 1410 - 1411 - 1412 - 1413 - 1414

## GUEDEL CANULAS

### Guedel pattern airway

Guedel pattern airway, made in nontoxic material, sterile and single-packed. **There are 6 available sizes, in 6 different colours.**

Item code	1409	1410	1411	1412	1413	1414	1415
Size	00	0	1	2	3	4	Set

Resuscitation



Art. 1080

## RESCUE MASK

### Resuscitation mask

Indispensable for all rescue operators, this mask avoids any contact between mouths during emergency reanimations. Made of transparent material, it is composed by a soft cushion and a non-return valve. The mask fits in a soft, shock-proof plastic case. The case colour could change.



Art. 1787

## KEYRING

### Keyring with mouth-to-mouth mask

Mouth-to-mouth safety mask that allows ventilation in complete safety thanks to a transparent sheet and a non-return valve which prevents the patient's exhalation from coming into contact with the rescuer. It is supplied in a handy container with keyring. The case colour could change.

RESUSCITATION

# RESUSCITATION MASKS



Art. 2400 - 2401 - 2402 - 2403 - 2404 - 2405

## SILICONE MASK

Resuscitation and anesthesia face mask

Resuscitation and anesthesia face mask. The anatomical border in silicone guarantees perfect fitting to the face while the shell, made in transparent material, allows the patient to be monitored during the application. Standard connector. Can be autoclave sterilized, available in 6 sizes.

Item code	2400	2401	2402	2403	2404	2405
Size	0	1	2	3	4	5



Resuscitation

156



Art. 2411 - 2413 - 2414 - 2415

## DISPOSABLE MASK

Disposable resuscitation mask

Disposable resuscitation face mask. Size 3-4-5 also feature valve and spider.

Item code	2411	2413	2414	2415
Size	1	3	4	5

## RESUSCITATION

# RESUSCITATION BAGS



Art. 1274 - 1276 - 1278

## RESUSCITATION BAGS

### Silicone resuscitation bags

Transparent, autoclavable, equipped with safety valve with adjustment at 60 mbar for adults and 40 mbar for children and new-borns. Equipped with connections for oxygen tank, reservoir and masks.

Art. 1230 - 1235

## DISPOSABLE BAGS

### Disposable resuscitation bags with mask

Disposable resuscitation bags for adults with disposable mask. Equipped with non-return valve.

Item code	1274	1276	1278	1230	1235
Volume (ml)	1700	500	350	1700	500
Autoclavable	✓	✓	✓	✗	✗
Safety valve	✓	✓	✓	✓	✓

### Accessories



**Art. 12028**  
Blue case for resuscitation bag and mask

Size: 22 x 30 x 12 cm (LxHxD). Weight: 0.8 kg.



**Art. 1056**  
Extension tube for oxygen f/f 2 metres long



**Art. 1057**  
Goggles for oxygen therapy with rear sliding fastener

# OXYGEN RESERVOIRS and **MANUAL CIRCUITS**



Art. 1258

## ADULT DISPOSABLE CIRCUIT

**Manual resuscitation circuit for adults**

Manual BVM resuscitation circuit for adults, equipped with:

- ✓ 1000 ml bag
- ✓ Pressure limit valve (40 cmH2O)
- ✓ Size 4 mask
- ✓ Reservoir
- ✓ Oxygen tube extension



Pediatric



Art. 1257

## DISPOSAL PEDIATRIC CIRCUIT

**Manual resuscitation circuit for children**

Manual BVM resuscitation circuit for children, equipped with:

- ✓ 550 ml bag
- ✓ Pressure limit valve (40 cmH2O)
- ✓ Size 3 mask
- ✓ Reservoir
- ✓ Oxygen tube extension



Pediatric



Art. 1256

## DISPOSAL NEONATAL CIRCUIT

**Manual resuscitation circuit for new-borns**

Manual BVM resuscitation circuit for new-borns, equipped with:

- ✓ 280 ml bag
- ✓ Pressure limit valve (40 cmH2O)
- ✓ Size 1 mask
- ✓ Reservoir
- ✓ Oxygen tube extension



Art. 1233

**Resuscitation bag - Oxygen reservoir with valve for adults**



Pediatric

Art. 1234

**Resuscitation bag - Oxygen reservoir with valve for children**



# AEREOSOL THERAPY MASKS



**Art. 1053**

## **Adult mask for inhalation therapy**

Soft mask for inhalation therapy equipped with ampulla and anti-flattening connection tube



**Art. 1052**

## **Oxygen therapy mask for adults - variable concentration - Venturi system**

Soft mask system equipped with coloured connections kit for oxygen modulation and anti-flattening connection tube



Pediatric

**Art. 1040**

## **Oxygen therapy mask for children**

High concentration oxygen therapy mask for children. Produced with non-toxic material, complete with adjustable elastic, clip for nose and connecting tube of cm 210



Pediatric

**Art. 1045**

## **Oxygen therapy mask for children with bag**

High concentration oxygen therapy mask for children. Produced with non-toxic material, complete with bag, adjustable elastic, clip for nose and connecting tube of cm 210



**Art. 1050**

## **Oxygen therapy mask for adults**

High concentration oxygen therapy mask for adults. Produced with non-toxic material, complete with adjustable elastic, clip for nose and connecting tube of cm 210



**Art. 1055**

## **Oxygen therapy mask for adults with bag**

High concentration oxygen therapy mask for adults. Produced with non-toxic material, complete with bag, adjustable elastic, clip for nose and connecting tube of cm 210



RESUSCITATION

# OXYGEN TANKS

Art.  
1465

Art.  
1468

Art.  
1470

Art.  
1475

Art.  
1478


Art. 1465 - 1468 - 1470 - 1475 - 1478

## OXYGEN TANKS

Tanks for oxygen therapy

Tanks for oxygen therapy in compliance with F.U. norms, produced with stainless steel, refillable and provided with test certificate.

Item code	Capacity	Gaseous capacity	Height	Diameter	Weight
1465	0,5 L	75 L	370 mm	60 mm	2,90 kg
1468	2 L	300 L	400 mm	120 mm	4,10 kg
1470	5 L	1000 L	550 mm	140 mm	9,20 kg
1475	7 L	1400 L	700 mm	140 mm	10,80 kg
1478	10 L	2000 L	955 mm	137 mm	13,25 kg

Art.  
1440

Art.  
1441

Art.  
1442


Art. 1440 - 1441 - 1442

## OXYGEN TANKS WITH PRESSURE REGULATOR

Oxygen tanks with pressure regulator and cover

Tanks for oxygen therapy in compliance with F.U. norms, produced with stainless steel, refillable and provided with test certificate.

Item code	Capacity	Height	Diameter
1440	2 L	435 mm	110 mm
1441	3 L	565 mm	110 mm
1442	5 L	620 mm	140 mm

# PRESSURE REGULATORS



Art. 1452

## STANDARD PRESSURE REGULATOR

Standard pressure regulator with dispenser

Diaphragm pressure reducer with manometer and dispenser to regulate the supply of oxygen



Art. 1455

## FLOW METER REGULATOR

Flow meter pressure regulator

Diaphragm pressure reducer with pressure gauge equipped with rotameter flowmeter 0 - 15 liters / min with hose connection. Safety valve with 7 bar calibration.



Art. 1460

## FLOW METER REGULATOR AND HUMIDIFIER

Flow meter pressure regulator with humidifier

Diaphragm pressure regulator equipped with pressure gauge rotameter flowmeter 0 - 15 liters / min and humidifier in makrolon. Safety valve with 7 bar calibration.

# FIRST AID EQUIPMENTS

8

Thanks to his own equipment, a rescuer can make a difference on the field. Besides the main rescue devices, many more products play a critical role in first aid operations. We have selected tens of products in order to enable any structure, team and rescuer to work everywhere in the best conditions. As far as rescue knapsacks are concerned, we have found numerous solutions suitable to the various existing specializations.

First aid equipments

162

Product Catalogue 2023







Discover MeBer products on-line.  
Shoot the QR code with your tablet or smartphone to land on the dedicated page on the website.  
You might be required to install a free application to navigate.

Alternatively, go to:  
<https://www.meber.it/en/products/first-aid-equipments.aspx>

First aid equipments

FIRST AID EQUIPMENTS

# RESCUE KNAPSACKS

NEW  
PRODUCT

Art. 1318

## ROMEO

### Rescue knapsack

Modular knapsack, compact and simple, suitable for volunteers and rescue professionals. Entirely made in red closed-cell PVC. Tearproof, washable and resistant, with transparent front pocket for customization and internal bands for objects, with double stitching, 10 mm strong zip, high visibility reflective bands on all external sides, PVC base and plastic supports. Equipped with an internal rigid protective shell structure made of thermoplastic material and with a large central compartment with double zip front opening and 2 external front pockets. The central unit compartment is equipped with Velcro® attachments for the anchoring of 4 standard colored cases (Art. 1362-002), and 2 professional bags (Art. 1362-001). Two rubber-coated handles complete the non-slip backpack for carrying by hand, together with its two shoulder straps reinforced with reflective bands.

**Size: 33 x 54 x 36 cm (W x H x D). Weight: 3,8 kg**

NEW  
PRODUCT

Art. 1319

## GIULIETTA

### Emergency bag

Red Emergency bag made in closed-cell PVC, with PVC bottom, plastic supports and upper transparent pocket for customizations. Compact but roomy, thanks to its internal modularity it allows the organization of space in a personalized way in order to optimally house what is necessary for emergencies. 4 pockets, 2 on the sides, one front and one rear, predisposed with supports and rubber bands. It is also equipped with an additional removable bag with transparent front. Thanks to its construction, it allows internal cleaning and disinfection. Portable by hand, shoulder strap or as a backpack.

**Size: 45 x 33 x 34 cm (W x H x D). Weight: 2,2 kg**





FIRST AID EQUIPMENTS

# RESCUE KNAPSACKS

Art. 1325/BR

Art. 1325/RG

Art. 1325

## ANTARTIDE

Complete compact rescue knapsack

The Antarctic model differs from the Artide model in the large front pocket with crystal window for the AED (any model of defibrillator) or other medical material. Made of blue / red or royal blue nylon with double stitching, large high visibility reflective bands, PVC bottom, reinforced shoulder straps for shoulder transport. The large external pocket favors the rescuer by giving a clear view of the content. Internally it is equipped with Velcro® attachments for the anchoring of 3 professional bags (blue, yellow and green), supplied with the knapsack.

**Size: 36 x 44 x 29 cm (W x H x D). Weight: 3,9 kg**

Art. 1326

## ARTIDE

Complete compact rescue knapsack

Rescue knapsack, made in royal blue tear-resistant nylon with double stitching, large high visibility reflective bands, PVC bottom, reinforced and detachable shoulder straps for shoulder transport. Internally it is equipped with Velcro® attachments for the anchoring of 3 professional bags (blue, yellow and green), supplied with the knapsack.

**Size: 36 x 44 x 23 cm (W x H x D). Weight: 3,6 kg**

FIRST AID EQUIPMENTS

# RESCUE KNAPSACKS

NEW  
PRODUCT

Art. 1322

## ATENE

### Red professional rescue knapsack

Athens is the backpack that meets the rescuer's needs: compact, functional and light. The central compartment is equipped with Velcro® attachments for the anchoring of six coloured bags (supplied with the knapsack). Equipped with a front pocket and two small double zip side pockets. The backpack is completed by three high visibility reflective bands (1 front and 2 on the sides) and front and side elastic bands to accommodate sanitary material of various kind.

**Size: 40 x 50 x 24 cm (W x H x D). Weight: 1,6 kg**



Art. 1390

## YANKEE

### Red Cordura Backpack

Compact and functional backpack. Made in tear-proof fabric, it is equipped with a front and a central pocket for medical supplies. Large high visibility reflective strips on the front side and on the shoulder straps. Particularly fit for footed rescuers or nurses for the stashing of medical equipment, thanks to its small size

**Size: 33 x 49 x 8 cm (W x H x D). Weight: 1 kg**





FIRST AID EQUIPMENTS

# RESCUE KNAPSACKS

Art. 1362

## AMBUMED

### Rescue knapsack

New modular knapsack, practical, compact, simple and suitable for rescue volunteers and professionals. Externally it is made of blue tear-resistant double stitched nylon, sturdy 10 mm zips, reflective high visibility strips on all sides and a closed-cell PVC bottom, with plastic supports. On the inner side there is a protective hard shell structure made of thermoplastic material, consisting of a central wide space, characterized by a double-zip opening and 4 external side pockets. The central space is supplied with Velcro® straps for the anchoring of 4 coloured cases, each one supplied with a large transparent window reinforced with double stitching, and 2 professional bags to store and organize different tools and materials. This knapsack is completed by a frontal pocket with inner straps to secure different objects, 3 anti-slip rubber handles (1 on the upper side and 2 on the lateral sides, 2 shoulder bands with reflecting stripes on the back.

**Size: 53 x 53 x 25 cm (W x H x D). Weight: 4,1 kg**



First aid equipments

Art. 1343

## RESCUEMED

### Rescue knapsack

Made of red tear-resistant material, the RESCUEMED pack is ideal for first aid operators thanks to its versatility and robust qualities. Ample front pockets to ensure that all necessary materials are at hand. Internally, it is equipped with smaller bags, with and without handles, permitting instant recognition of the instruments necessary for the intervention. Two yellow reflective strips complete the features of this first aid knapsack.

**Size: 37 x 52 x 29 cm (W x H x D). Weight: 3,5 kg**



Product Catalogue 2023

167

FIRST AID EQUIPMENTS

# RESCUE KNAPSACKS



Art. 1354

**115**

## Red knapsack with side pockets and compartments

Knapsack made of red anti-tear nylon, double seams, reflective strips, reinforced shoulder straps and sturdy 10 mm zips. Supplied with front opening wide space inside of which are 4 rack pockets to store various devices. Possibility to fix 2 litres oxygen tanks. Two wide lateral sections contain four bags in different colours allowing a fast identification of the kits. PVC Water-proof bottom.

**Size: 44 x 56 x 26 cm (W x H x D). Weight: 4 kg**



Art. 1354-001

## Carrying bag for knapsack "115"

Carrying bag for Art. 1354, produced with fire-proof PVC. This bag ensures stowage and transport of the knapsack to preserve its integrity.

**Weight: 1,2 kg**



# ACCESSORIES FOR RESCUE KNAPSACKS

Art. 1373

## Green multi-functional roll bag

Practical and functional accessory, designed for the placement of medical, intubation and resuscitation devices. Internally the roll is equipped with rubber bands that block medical instruments and a transparent zippered pocket. Innovative accessory, adaptable to any type of rescue backpack.



Art. 1362-001

## Square bags

Yellow, green, red and blue square bags for knapsacks

Size: 23 x 23 x 10 cm (W x H x D).

Art.  
1362-001/V

Art.  
1362-001/G

Art.  
1362-001/R

Art.  
1362-001/B



Art. 1362-002

## Rectangular bags

Yellow, orange, red and blue rectangular bags for knapsacks.

Size: 30 x 15 x 10 cm (W x H x D).

Art.  
1362-002/R

Art.  
1362-002/Y

Art.  
1362-002/B

Art.  
1362-002/A



Art. 1362-003

## Rectangular bags

Professional rectangular bags for rescue knapsacks. Available in yellow (Art. 1362-003/G), orange (Art. 1362-003/A), Blue (Art. 1362-003/B) and green (Art. 1362-003/V), also in its pediatric version (Art. 1362-003/PED).

Size: 35 x 19 x 14 cm (W x H x D).

Weight: 0,5 kg

Art.  
1362-003/G

Art.  
1362-003/A

Art.  
1362-003/B

Art.  
1362-003/V

Art.  
1362-003/PED



Pediatric



FIRST AID EQUIPMENTS

# EMERGENCY RESCUE BAGS



Art. 1365

## NEVADA

### Emergency rescue bag

Bag made of royal blue coloured nylon. It is equipped with high visibility reflective stripes, sturdy double zips, a comfortable shoulder strap with metal hooks and handles to carry by hand. The bag is made up of a central space divided into 3 wide compartments with rigid adjustable dividers. Externally it is characterized by a frontal compartment with a pocket and fixed dividers, 3 side pockets with zip and two open compartments for the placement of bottles of various kinds. It is completed by a customizable name plate on the lateral side. The bottom is in water-proof PVC with plastic support feet. The linings can be easily removed to allow internal cleaning and disinfection.

**Size: 50 x 25 x 28 cm (W x H x D). Weight: 1,5 kg**

First aid equipments



Art. 1345

## EMERGENCY RESCUE BAG

### Empty emergency bag



170

Emergency rescue bag made of Nylon 840. It is compact but roomy and, thanks to its internal modularity, it allows personalized organization of the space so that everything required for an emergency can be housed in the best way. It is made up of 4 identical pockets, 2 side ones, a front one and a back one, which are prepared and provided with supports and elastics for holding medical material. It also houses another small removable bag with a transparent front that allows the content to be seen. Thanks to its construction, the linings can be completely removed to allow internal cleaning and disinfection. Can be carried by hand, over the shoulder or as a backpack.

**Size: 45 x 33 x 34 cm (W x H x D). Colour: Red. Weight: 1,8 kg**



FIRST AID EQUIPMENTS

# EMERGENCY RESCUE BAGS



Pediatric

Art. 1342

## BABY BAG

Empty emergency bag



Pediatric emergency rescue bag made in Nylon with waterproof PVC coating. It consists of three inner compartments, a removable bag with transparent front that enables to see the content. It is equipped with 4 external pockets, two at the sides and two at the front. It can be carried by hand, over the shoulder or as a backpack.

**Size: 45 x 30 x 30 cm (W x H x D). Weight: 1,4 kg**

Art. 1370/RB

## DEFIBRILLATOR CASE

Defibrillator case with front pocket

Universal defibrillator case, practical and functional. It is equipped with a double-zip opening, a water-proof PVC bottom and support feet positioned both on the bottom and on the back side, in order to ensure that the medical device is stable and protected. On the inner side there is a section with Velcro® straps to hold the device and a transparent pocket to store spare electrodes. Externally it has a big frontal service pocket, allowing to carry a whole resuscitation kit. It is completed with an anti-slip rubber handle and a shoulder strap with metal hooks. It can be carried both by hand and over the shoulder and it is available in red or blue colour.

**Size: 28 x 32 x 20 cm (W x H x D). Weight: 1 kg**



FIRST AID EQUIPMENTS

# POUCHES and VIAL HOLDERS

Art. 11342

## POUCH

**Multipocket pouch (Blue or red)**

Practical multi pockets pouch designed for professional rescuers and for volunteers. It is divided into different external and internal compartments with 1 practical frontal pocket and 1 little side pocket.

**Size: 18 x 28 x 13 cm (W x H x D).**



Art. 1350

## VIAL HOLDER

**Bag with partitions for vials transport**

Indispensable vial holder for protection and identification of any kind of vials. Made of nylon, it consists of 2 compartments, one split in 3 parts and one in 2 parts. Equipped with elastic bands to fix up to 90 vials and internal transparent windows for customization.

**Size: 18 x 28 x 13 cm (W x H x D).**



# PHARMA-MED



Art. 1375

## PHARMA-MED

### Professional bag with compartments for drugs transport

Made of tearproof, waterproof and fireproof nylon, it is equipped with sturdy zippers, a comfortable handle covered in rubber and 4 support feet in plastic material to isolate the device from the support surface. It consists of three different sections. The sector Liquids and Powders is centrally positioned and is divided into 3 sections, equipped with elastic bands for vials, interspersed with two polyurethane foam membranes (slow-recovery) to protect the vials from breaking. The COLD sector, marked by the red Cordura, is equipped with elastic bands for vials on one side and reusable hot / cold gel bags on the other (supplied with 6 Mini DispoGels). The Capsule compartment consists of a transparent pocket with a zip and a transparent pocket with Velcro® straps for boxes or blisters. Pharma-Med is equipped with a transparent external pocket for the placement of a check list or a medication expiration date.

**Size: 20 x 16 x 28 cm (W x H x D).**

**Number of vials contained: min 108 – max 120. Weight: 1,1 kg**



First aid equipments

173




Accessories



Art. 1375-001

**Re-usable HOT/COLD  
Gel bag 5,5x12 cm for  
drug transport bags**



Product Catalogue 2023



FIRST AID EQUIPMENTS

# INSTRUMENT HOLDERS

Art. 1355/B

## INSTRUMENT HOLDER

MeBer instrument holder

Belt instrument holder made of blue tear-resistant material. It can hold essential first-aid instruments such as sphygmomanometers, phonendoscopes, clothes-cutting scissors, latex gloves, etc.

**Size: 11 x 18 x 10 cm (W x H x D)**



Art. 1356/B

## INSTRUMENT HOLDER

Blue MeBer instrument holder with reflective strip

Handy instrument holder to be fixed to the waist belt. Equipped with three pockets and a reflection band, it is an ideal solution to easily carry any essential first aid tool.

**Size: 18 x 13,5 x 2,5 cm (W x H x D)**





FIRST AID EQUIPMENTS

# TOOLS KIT and ASPIRATORS



Art. 1360

## TOOLS KIT

### Emergency vehicles tools kit

Tools kit for road accidents, ideal for all rescue vehicles.  
All tools needed are included inside of a convenient nylon case.

**Size: 48 x 28 x 12 cm (W x H x D)**

Content:

- ✓ 1 axe
- ✓ 1 crowbar
- ✓ 2 wind-proof torches
- ✓ 2 pairs of professional gloves
- ✓ 1 bolt cutter
- ✓ 1 tow cable
- ✓ 2 smoke candles
- ✓ 1 fire blanket

First aid equipments

Art. 11068

## Manual suction unit

This device was designed to ensure optimal pharyngeal and tracheal aspiration in emergency procedures.  
The secretion collection jar volume is 300 ml.

**Size: 12,5 x 6 x 16 cm. Weight 200 g**



175

Product Catalogue 2023

# CORMOS

PORTABLE ASPIRATOR FOR SECRETIONS

Art. 11010/FA

With autoclavable  
canister



Art. 11014



Art. 11010/LN

With disposable liner

Art. 11010

## CORMOS

Portable aspirator for secretions

Compact and powerful, indicated for emergency services. CORMOS works directly with the 12V of the vehicle or with an inner accumulator with more than 45 min. of autonomy. CORMOS is silent and powerful thanks to its engine of 24 litres/min which enables a constant aspiration, with control manometer, positioned on the front of the aspirator, up to 850 mbar. Easy to carry thanks to its compact dimension, the anatomic handle and its limited weight. New case design, made of ABS/UL high resistance anti-shock material. It is supplied with connection cable to the 12V of the vehicle, carrying case, wall locking device 10G proof as per EN 1789 normative for the use inside emergency vehicles, 1000 ml autoclavable canister with safety valve and antibacterial filter to fix directly onto the cap.

Item complies to ISO 10079-1; IEC 600601-1 3rd ed., approval e1 according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11.

Classification according MDD 93/42/EEC: IIb.



Art. 11128/FA



Art. 11128/LN

Art. 11128

## On board medical suction unit

### On board suction unit with autoclavable canister

- ✓ On board suction unit for rescue vehicles
- ✓ Maximum aspiration flow 30lt/min.
- ✓ Maximum aspiration Power 800 mbar.
- ✓ 1000 ml Autoclavable canister with disposable cartridge (model LN)
- ✓ 1000 ml autoclavable canister with safety valve and protection filter directly fixed onto the cap (model FA)
- ✓ Engine block size: approx. 18x19x9 cm. Built-in regulator block size: 9x13x8.5 cm approx.
- ✓ In compliance with all existing norms.

<p>Supplied with</p>	<p><b>Art. 11028</b> Autoclavable collection canister ml 1000 for disposable liner</p>	<p><b>Art. 11023</b> 10 Proof wall mount with integrated CORMOS and CORMOS PLUS charging system</p>	<p><b>Art. 11016</b> Green case</p>
	<p><b>Art. 11014</b> Rigid aspiration probe</p>	<p><b>11034</b> 220V 2-pole power supply. Output voltage 13.7 V</p>	

<p>Accessories</p>	<p><b>Art. 11162</b> 1 litre disposable vase for secretion collection</p>	<p><b>Art. 11024 (6 pcs.) Art. 11025 (12 pcs.) Art. 11026 (36 pcs.)</b> Disposable 1000 ml bag</p>	<p><b>Art. 11012</b> Autoclavable vase for collection</p>	<p><b>Art. 11013</b> Antibacterial filter</p>
				<p><b>Art. 11022-001</b> 2011 Power supply / 3-pole 12VCC for suction unit</p>

# CORMOS

## PLUS

PORTABLE ASPIRATOR  
FOR SECRETIONS



Art. 11022/FA

Portable aspirator  
for secretions with  
autoclavable canister

Art. 11022

## CORMOS PLUS

Portable aspirator for secretions

The CORMOS PLUS suction units are easy-to-use and compact battery powered medical devices. They differ from the CORMOS series for their higher power supplied by a suction flow of 30 liter/minute engine. The suction rate can be controlled thanks to a 800 mbar manometer. Supplied with: carrying case, 10G proof wall mount as per EN1789 for rescue vehicles operations, 1000 ml autoclavable canister with safety valve and anti-bacterial filter, tube with Jankauer probe and universal battery charger.

Item complies to ISO 10079-1; IEC 600601-1 3rd ed., approval e1 according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11.

Classification MDD 93/42/EEC: IIb.





**Art. 11022/LN**

**Portable aspirator with  
vase and adapter for  
disposable bags**

Supplied with



**Art. 11028**

Autoclavable  
collection  
canister ml 1000  
for disposable  
liner.



**Art. 11023**

10 Proof wall  
mount with  
integrated  
CORMOS and  
CORMOS PLUS  
charging system



**Art. 11032**

Green case



**Art. 11014**

Rigid aspiration  
probe



**11034**

220V 2-pole  
power supply.  
Output voltage  
13.7 V

Accessories



**NEW  
PRODUCT**



**Art. 11162**

1 litre disposable  
vase for secretion  
collection



**Art. 11024 (6 pcs.)  
Art. 11025 (12 pcs.)  
Art. 11026 (36 pcs.)**

Disposable 1000 ml bag



**Art. 11012**

Autoclavable vase  
for collection



**Art. 11013**

Antibacterial  
filter



**Art. 11022-001**

2011 Power supply  
/ 3-pole 12VCC for  
suction unit

FIRST AID EQUIPMENTS

# FIRST AID EQUIPMENTS

Art. 706 - 708

## TRANSFER SHEET

Transfer sheet with 6/8 handles

Transfer sheet made of fire-proof, tear-proof, water-proof nylon.  
Available with 6 or 8 strong handles.

Capacity 160 kg. Size: cm 192 x 75

Art. 706



Art. 708



Art. 712

## TRANSFER SHEET

Obese transfer sheet with 12 handles

Made of fire-proof, tear-proof, water-proof nylon. Equipped with  
12 reinforced handles.

Capacity 350 kg. Size: cm 274 x 150



Bariatric



Art. 1950

## SCISSORS

Emergency scissors

Emergency scissors with thermoplastic handles and steel  
blade for clothes cutting. Available in 3 different colours: black  
(1950/N), yellow (1950/G) and red (1950/R).

Length cm 18

Art. 1950/N

Art. 1950/R

Art. 1950/G



Art. 1980



Art. 1982



Art. 1980 - 1982

## SCISSORS

### Emergency scissors

Multi-purpose emergency scissors, ideal for rescue operators. Produced in stainless steel and plastic-coated handles. Ideal for glass breaking, clothes cutting and safety belts cutting in emergency rescue. Supplied with wall mount.

Art. 1910 - 1905

## ICE BAG

### Instant ice bag

Disposable instant ice bag with a 30 minutes constant-temperature autonomy. Available in PVC bag (art. 1910) or TNT bag (Art. 1905), which produces less outer condensation.

**25 pieces per box. Size: cm 14 x 18.**



Art. 1910

Art. 1905

First aid equipments

Art. 1915

## SPRAY ICE

### Instant spray ice

Instant spray ice in spray can (400 ml), idea for bruises and traumas.




Art. 1960

## TOURNIQUET

### Adjustable tourniquet

Adjustable tourniquet with quick opening for a faster application. Fit for rescue knapsacks and bags



Product Catalogue 2023

FIRST AID EQUIPMENTS

# FIRST AID EQUIPMENTS

Art. 12026

## GLOVE BOX

3-tier glove box for rescue vehicles

Safe ABS box for the storage of 3 disposable gloves. Equipped with safe-lock mechanisms against accidental opening. The smooth surface and rounded corners make this a safe and easy-to-clean piece of equipment.

Can be equipped on any rescue vehicle.

**Height 41 cm, Width 26 cm, Depth 9,5 cm, Weight 1 kg**



Art. 11360/B



Art. 11360/G



Art. 11360/R



Art. 11360

## HELMET

Protective emergency helmet in white, red or yellow ABS

Anti-suffocation throat-strap with automatic release between 15 dN and 25 dN. Inner washable padding, spider shape. Conform to CE EN 397

Art. 11372

Clear visor



Accessories





Art. 1776/R - 1776/AZ

## RESCUE BLANKET

### Pile rescue blanket

Red and blue pile blanket, sewn on all perimeter, supplied with 2 reflecting bands

**Size: cm 140x210.**



Art. 1776/R

Art. 1776/AZ



Art. 1770

## RESCUE BLANKET

### Isothermal blanket

Necessary article for the first aid. Two-sides, one gold and one silver, this blanket protects from severe weather conditions (heat or freezing). Packed in a re-closable bag.

**Size: cm 160x210.**

Art. 1385 - 1388

## PILLOWCASE

### Disposable TNT pillowcase

Disposable TNT pillowcase. Available in green (25 gr/mq – cm 50x70) and white (28 gr/mq – cm 50x80). Packed in 50 pieces

Art.  
1385

Art.  
1388



Art.  
1380

Art.  
1382



Art. 1380 - 1382

## BED STRETCHER SHEET

### Disposable bed stretcher sheet in TNT

Disposable bed stretcher sheet in TNT. Available in green (25 gr/mq – cm 90x230) and white (28 gr/mq – cm 100x220). Packed in 50 pieces.

# FIRST AID EQUIPMENTS

Art. 11264

## Mini first aid box mod. 388/3 annex 2

For companies or production units with less than three workers, not belonging to group A. Fit for small companies and service vehicles. Equipped with wall mount.

**External size: 25x19x9 cm. External size: kg 1,4**

- ✓ 1 copy of Min. Decree 388, 15.07.03
- ✓ 2 pairs of sterile gloves
- ✓ 1 bottle of disinfectant 500 ml IODOPOVIDONE, 10% Iodine PMC
- ✓ 1 bag of sterile saline solution 250 ml EC
- ✓ 1 compressed sterile gauze cm 18x40
- ✓ 3 compressed sterile gauze cm 10x10
- ✓ 1 sterile forceps
- ✓ 1 pack of cotton wool
- ✓ 1 PLASTOSAN, 10 assorted patches
- ✓ 1 Adhesive plaster spool m 5x2,5 cm
- ✓ 1 3.5x10 cm gauze bandage
- ✓ 1 pair of scissors 14.5 cm DIN 58279
- ✓ 1 Tourniquet
- ✓ 1 ICE PACK disposable instant ice
- ✓ 1 bag for medical waste
- ✓ 1 MULTILINGUAL first aid instructions



Art. 11260

## First aid cabinet mod. min. decree 388/1 annex 1

For companies belonging to group A and B, with 3 or more workers

**External size: 46x30x14 cm. External size: kg 7,5**

- ✓ 1 copy of Min. Decree 388, 15.07.03
- ✓ 5 pairs of sterile gloves
- ✓ 1 mask with splash guard visor
- ✓ 2 bottles of disinfectant 500 ml IODOPOVIDONE, 10% Iodine PMC
- ✓ 3 bottles of sterile saline solution 500 ml EC
- ✓ 2 compressed sterile gauze cm 18x40
- ✓ 10 compressed sterile gauze cm 10x10
- ✓ 2 Sterile sheets cm. 40x60 DIN 13152-BR
- ✓ 2 sterile forceps
- ✓ 1 pack of cotton wool
- ✓ 2 PLASTOSAN, 10 assorted patches
- ✓ 1 Box tubular elastic bandage
- ✓ 2 Adhesive plaster spool m 5x2,5 cm
- ✓ 1 pair of scissors 14.5 cm DIN 58279
- ✓ 3 Tourniquets
- ✓ 2 ICE PACK disposable instant ice
- ✓ 2 bags for medical waste
- ✓ 1 PERSONAL sphygmomanometer with stethoscope
- ✓ 1 CE clinical thermometer with case
- ✓ 1 MULTILINGUAL first aid instructions

Art. 11270  
Spare pack for  
art. 11260



Art. 11262

## Abs first aid box mod. min. decree 388/2 annex 1

For companies belonging to group A and B, with 3 or more workers. Equipped with wall mount.

**External size: 44x34x15 cm. External size: kg 8,3**

- ✓ 1 Copy of Min. Decree 388, 15.07.03
- ✓ 5 Pairs of sterile gloves
- ✓ 1 Mask with splash guard visor
- ✓ 3 Bottles of sterile saline solution 500 ml. CE
- ✓ 2 Bottles of disinfectant 500 ml. IODOPOVIDONE 10% iodine PMC
- ✓ 10 compressed sterile gauze cm. 10x10
- ✓ 2 compressed sterile gauze cm. 18x40
- ✓ 2 Sterile sheets cm. 40x60 DIN 13152-BR
- ✓ 1 Packet cotton wool
- ✓ 2 Sterile forceps
- ✓ 1 Box tubular elastic bandage
- ✓ 2 Reels adhesive plaster 5m x 2.5cm.
- ✓ 2 PLASTOSAN 20 assorted plasters
- ✓ 1 Pair bandage cutting scissors cm. 14.5 DIN 58279
- ✓ 3 Tourniquets
- ✓ 2 ICE PACK disposable instant ice
- ✓ 2 Bags for sanitary waste
- ✓ 1 Clinical thermometer CE with case
- ✓ 1 PERSONAL Sphygmomanometer with phonendoscope
- ✓ 4 Hemmed bandages 3.5m x 10cm.
- ✓ 1 PIC 3 box containing:
- ✓ 3 liquid soap, 3 disinfectant wipes, 2 ammonia wipes
- ✓ 1 PINOCCHIO+VENTO, complete kit 1 Isothermal blanket, gold/silver cm. 160x210
- ✓ 1 Triangular sheet non-woven fabric cm. 96x96x136
- ✓ 2 Preparations 3.5 g sterile GEL FOR BURNS
- ✓ 2 Elastic bandages 4m x 6cm. DIN 61634
- ✓ 1 Sterile compressed bandages mm. 80x100 DIN13151 M
- ✓ 1 MULTILINGUAL first aid instruction booklet



# DISPOSABILE GLOVES



## LATEX GLOVES

### Disposable latex gloves

High quality natural latex gloves. Both powdered or powder-free, non-sterile, ambidextrous.

**Available in 4 sizes, packed in boxes of 100 pieces**

#### POWDERED

Art. 1930 SIZE S	Art. 1932 SIZE M	Art. 1934 SIZE L	Art. 1935 SIZE XL
---------------------	---------------------	---------------------	----------------------

#### POWDER-FREE

Art. 1936 SIZE S	Art. 1937 SIZE M	Art. 1938 SIZE L	Art. 1939 SIZE XL
---------------------	---------------------	---------------------	----------------------



#### NITRILE

Art. 1940 SIZE S	Art. 1942 SIZE M	Art. 1944 SIZE L	Art. 1945 SIZE XL
---------------------	---------------------	---------------------	----------------------

## NITRILE GLOVES

### Powder-free disposable nitrile gloves

Disposable, non-sterile, nitrile medical glove suitable for people allergic to natural latex rubber and sensitive to type I allergic reactions (immediate hypersensitivity). Light blue color, odorless, uniform thickness and free from bubbles, excellent tactile sensitivity and dexterity, increased comfort even during prolonged operations. Anatomical shape, ambidextrous, with anti-roll-up reinforced elastic cuff, absence of latex proteins and chemical processing residues, free from lubricating powder.

**Available in 4 sizes, packed in boxes of 100 pieces**



#### DOUBLE THICK LATEX

Art. 1925 SIZE S	Art. 1926 SIZE M	Art. 1927 SIZE L	Art. 1928 SIZE XL
---------------------	---------------------	---------------------	----------------------

## DOUBLE THICK LATEX GLOVES

### Powder-free, double thick, disposable blue latex gloves

Disposable high thickness natural rubber latex gloves, increased length, powder-free, high biocompatibility, micro-rough, controlled and certified quality.

Subjected to halogenation process. Increased resistance compared to that of standard disposable gloves.

**Available in 4 sizes, packed in boxes of 50 pcs**

# FASTENING SYSTEMS

9

Fastening systems



Discover MeBer products on-line  
Shoot the QR code with your tablet  
or smart-phone to land on the  
dedicated page on the website.  
You might be required to install a  
free application to navigate.

Alternatively, go to: <https://www.meber.it/en/products/fixing-systems.aspx>



# FASTENING SYSTEMS FOR STAIR CHAIRS

Art. 975

## ERGOLIFT Fastening system

### Extra Ergolift Stair chair fastening system

Designed to secure the Extra Ergolift Stair Chair to an emergency vehicle. Built in stainless steel, this device is corrosion-resistant and easy-to-clean. Its innovative mono-block configuration ensures a quick and easy locking and unlocking of the chair. This fastening device is compatible with all Extra chairs as well as Zena: 670/BR - 676 - 678 - 679. **Device certified according to the European harmonized safety standard UNI EN 1789.** Size: 58 x 111 x 23 cm (WxHxD)



Art. 960

## EXTRA Fastening system

### Extra stair chair fastening system

Designed to secure the Extra stair chair to an emergency vehicle. Built in stainless steel and thermoplastics, corrosion-resistant and easy-to-clean. This device ensures a quick, easy and efficient fastening. The operator can easily move the stair chair before and after locking, thanks to the innovative mono-block configuration. **Device certified according to the European harmonized safety standard UNI EN 1789.** Size: 23 x 42 x 20 cm (WxHxD)



Art. 977

## KUMA Fastening system

### Kuma stair chair fastening system

Designed to secure the Kuma stair chair to an emergency vehicle. Corrosion-resistant and easy-to-clean. Size: 55 x 120 x 30 cm (WxHxD)



Art. 618/SUP

Art. 618/INF



### Art. 618 Fastening system for carry chair

Designed to secure all Maya and Quick carry chairs to an emergency vehicle. 2-piece, built in stainless steel. **Device certified according to the European harmonized safety standard UNI EN 1789.**




### Art. 988 Extra Ergolift FX (Art. 675) and Proxima Pro (Art. 683) fastening system

Floor fastening system with integrated unlocking device for Art. 675 Extra Ergolift FX and Art. 683 Proxima Pro. **Device certified according to the European harmonized safety standard UNI EN 1789.**



• Also available separately

# FASTENING SYSTEMS FOR **AMBULANCE COTS**



Art. 956 PROOF

Art. 955 PROOF

## FASTENING SYSTEM

## FASTENING SYSTEM

Designed to ensure maximum stability to the stretcher on the emergency vehicle, these fastening systems are entirely built in stainless steel to offer increased durability and corrosion-resistance. Thanks to their structure, these easily-installed fastening systems ensure perfect containment, in any direction, of the cot or trolley as well as of the patient. The Art. 955 PROOF is fit for all 10G ambulance cots, while the Art. 956 PROOF is fit for the "TITANUS" trolley. **Device certified according to the European harmonized safety standard UNI EN 1789.**

Art.	Front height	Rear height	Width	Length	Weight
955 PROOF	13 cm	2,6 cm	48 cm	177 cm	8,5 kg
956 PROOF	13 cm	2,6 cm	42 cm	175 cm	8,5 kg

### Accessories

Art. 955-001



Front strike plate for 955 - 956 Proof certified fastening systems

Art. 955-002



Rear strike plate for 955 - 956 Proof certified fastening systems

Art. 980

**"MEBER-X" fastening kit for multi-level stretcher, including front lock, lateral lock and rear safety lock.**



**Front lock**  
Length: 54 cm  
Width: 68 cm  
Height: 14 cm



**Lateral lock**  
Length: 75 cm  
Width: 8 cm  
Height: 24 cm



**Rear safety lock**  
Length: 10 cm  
Width: 2 cm  
Height: 5 cm

# FASTENING SYSTEM FOR SCOOP STRETCHERS AND FIXED STRETCHERS

## Art. 959

### Wall mount for Variable geometry ERGON scoop stretcher

Wall mount for the ERGON scoop stretcher, made in thermo-plastics. This device ensures a safe and quick installation on any emergency vehicle, vertically. **Device certified according to the European harmonized safety standard UNI EN 1789.**



## Art. 614

### Horizontal fastening system for scoop stretchers

Aluminium support with safe quick-locking device for wall mount or horizontal storage of a scoop stretcher. **Device certified according to the European harmonized safety standard UNI EN 1789.**



## Art. 616

### Certified Vertical fastening system for scoop stretchers



## Art. 616/A

### Certified Vertical fastening system for scoop stretchers



Aluminium support with safe quick-locking device for wall mount or vertical storage of a scoop stretcher. The stretcher can be folded and then locked. **Device certified according to the European harmonized safety standard UNI EN 1789..**



## Art. 963

### Front fastening system for fixed stretchers Art. 1510 PROOF and 1515 PROOF

Max Width 18,5 cm  
Length 9,5 cm  
Height 5 cm



# FASTENING SYSTEMS

Art. 962

## EVOX PROOF Fastening system EN 1789

Fastening system for EVOX art. 7112 Proof, built entirely in stainless steel to offer corrosion-resistance and durability. Thanks to its structure, this device can withstand intense stress up to 10G in all directions, as per UNIEN 1789.



Art. 964

### Front lock for art. 638 - 639 - 640



Art. 965

### Front lock for ambulance cot

Blocking device, to be installed on the floor of an ambulance, enabling the locking of an ambulance cot towards its front part. Built in stainless steel.



Art. 966

### Front lock for ambulance trolley

Blocking device, to be installed on the floor of an ambulance, enabling the locking of an ambulance trolley Art. 7020 towards its front part. Built in stainless steel.



Art.  
950/1  
H. 188 mm



Art.  
950/2  
H. 160 mm



Art.  
950/3  
H. 148 mm



Art.  
950/4  
H. 129 mm

### Rear locking systems for ambulance cots

Built in steel, these systems are installed on the ambulance floors to safely lock an ambulance cot. Available in 4 different heights.



Art. 1720

### Block mm 30 with lever

Aluminium block, height 30 mm, with security lever.



Art. 1722

### Block mm 30 with lever and little plate

Aluminium block, height 30 mm, equipped with rotating security lever and a little plate for the wall mount.



Art. 1725

### Block mm 73 with lever and little plate

Aluminium block, height 73 mm, equipped with rotating security lever and a little plate for wall mount.



Art. 970

### Complete din blocking device

Device for the locking on stainless steel supports of the rigid stretchers (Art. 1510-1515), DIN sizes. Set with 4 pieces.



# FASTENING SYSTEMS

## COMPARATION CHART FOR FASTENING SYSTEMS

Art.	Suitable for	Art.	Suitable for
 <b>614</b>	Maxima [630, 631] - 6200	 <b>975</b>	Extra Ergolift [671] - Extra [670, 670/BR] - Extra Lite [676, 678] - Extra Lite Zena [679]
 <b>616</b>	Maxima [630, 631]	 <b>960</b>	Extra [670, 670/BR] - Extra Lite [676, 678] - Extra Lite Zena [679]
 <b>616/A</b>	Maxima [630, 631]	 <b>962</b>	Evov Proof [7112]
 <b>618</b>	Maya [650, 653] - Quick [658, 659]	 <b>963</b>	1510 Proof - 1515 Proof
 <b>988</b>	Ergolift fx [675] - Proxima Pro [683]	 <b>964</b>	638 - 639 - 640
 <b>950/1</b>	Winner [910, 914]	 <b>965</b>	Frog [7210] - Frog Plus [7240] - Winner [910, 914] - Gecko [7122/V, 7122/AZ]
 <b>950/2</b>	Winner [910, 914]	 <b>966</b>	Titanus [F7020]
 <b>950/3</b>	Frog [7200, 7210] - Frog plus [7230, 7240] - Winner [910, 914] - Gecko [7122/V, 7122/AZ]	 <b>970</b>	1510 Proof - 1515 Proof
 <b>950/4</b>	Titanus [F7020] - 638 - 639 - 640	 <b>977</b>	Kuma [668]
 <b>955 PROOF</b>	Mercury [7070, 7072, 7076, 7078] - Mercury Lite [7080, 7082] - Mercury 5 [7075/4RG, 7077/4RG] - Mercury 5 Lite [7095/4RG, 7097/4RG] - Frog [7210 Proof] - Frog Plus [7240 Proof] - Frog Lite [7260 proof] - Winner [910 Proof, 914 Proof] - 7270	 <b>1720</b>	Galaxy [9010, 9012, 9016] - Duck [9022, 9024]
 <b>956 PROOF</b>	Titanus [7060, 7062]	 <b>1722</b>	Galaxy [9010, 9012, 9016] - Duck [9022, 9024]
 <b>959</b>	Ergon [16150/G, 16150/V] - Ergon Junior [16160/G, 16160/V]	 <b>1725</b>	Multi purpose fixing system
		 <b>980</b>	Meber-X [7300]

# RESTRAINT SYSTEMS

10

Restraint systems

192

Product Catalogue 2023



Discover MeBer products on-line  
Shoot the QR code with your tablet  
or smart-phone to land on the  
dedicated page on the website.  
You might be required to install a  
free application to navigate.

Alternatively, go to: <https://www.meber.it/en/products/restraint-systems.aspx>

# SPIDER STRAPS FOR SPINAL BOARDS AND ERGON



Art. 9156

## X-FIX

**Straps restraint system with carabiners**

Standard restraint straps with carabiners: innovative system designed to ease rescue operations, ensuring high comfort and spinal stability.



Art. 690

## SPIDER

**Spinal board straps**

This strap system ensures a precise patient restraint. Thanks to the transversally sliding VELCRO® straps, this device can fit any spinal board and any patient.



Art. 692

## SPIDER BABY

**Pediatric spinal board straps**

This strap system ensures a precise patient restraint on a "DUCK" pediatric board. Thanks to the transversally sliding VELCRO® straps, this device can fit any spinal board and any patient.

# STAIR CHAIRS AND AMBULANCE COTS **BELTS AND STRAPS**



Art. 610

**Black Thoracic-abdominal belt with 4-buckles system for padded chair.**



Art. 612/N-2/MEB

**Black/orange or full black adjustable chest belt with 4-buckles system and MEBER logo for ambulance cots and stair chairs**



Art. 612/AN-2/MEB

# BELTS



**Art. 605/C** 1 piece orange belt mm 50, plastic buckle and MEBER logo for mattresses proof, length 115 cm



**Art. 605/Nmeb** 1 piece black belt with plastic buckle and MEBER logo, 50 mm



**Art. 605/MEB** 1 piece orange belt with plastic buckle and MEBER logo, 50 mm



**Art. 606-A/MEB** 1 piece orange belt with metal buckle and MEBER logo (Length 190 cm)



**Art. 606-N/MEB** 1 piece black belt with metal buckle and MEBER logo (length 190 cm)



**Art. 607/A** 2-pieces orange belt with plastic buckle for stretchers, 50 mm



**Art. 607/A-MEB** 2-pieces orange belt with plastic buckle and MEBER logo, 50 mm



**Art. 607/B** Blue belt, 2-pieces with plastic buckle for stretchers, 50 mm



**Art. 608-A/MEB** 2-pieces orange belt with metal buckle and MEBER logo (Length 155 cm)



**Art. 608/N** 2-pieces black belt with metal buckle (Length 155 cm)



**Art. 608-N/MEB** 2-pieces black belt with metal buckle and MEBER logo (Length 155 cm)



**Art. 572** 2-pieces blue belt with quick-release metal buckle and carabiners, MEBER logo



**Art. 568** 2-pieces orange belt with plastic buckle and MEBER logo, 40 mm



# BELT SETS WITH CARRYING POUCH



**Art. 566** Black pouch with 3 belts Art. 607/Nmeb, 2-pieces black belt with orange MEBER logo for EXTRA stair chairs



**Art. 574** Blue pouch with 3 belts Art. 572, 2-pieces blue belt with quick-release metal buckle and carabiners, MEBER logo



**Art. 584** Black pouch with 3 belts art. 568 for ERGON stretchers



**Art. 604/MEB** Orange pouch with 2 belts 605/MEB



**Art. F602/B** Blue pouch with 3 blue belts with plastic buckle



**Art. 602/MEB** Orange pouch with 3 belts 605/MEB



**Art. 590/N** Black nylon pouch with 3 black belts Art.605/Nmeb for stair chairs, 1 piece belts with plastic buckle



**Art. 603/MEB** Nylon pouch with 3 belts Art. 607/Ameb



**Art. 603** Nylon pouch with 3 belts Art. 607/A



**Art. 600** Pouch with 3 belts Art. 606



**Art. 600/N** Pouch with 3 belts Art. 606-N/MEB



**Art. 601** Pouch with 3 belts Art. 608



**Art. 601/N** Pouch with 3 belts Art. 608-N/MEB


**NEW**  
PRODUCT

**NEW**  
PRODUCT

# AMBULANCE COTS **SUPPORTS**

11

Ambulance cots supports



Discover MeBer products on-line.  
Shoot the QR code with your tablet  
or smart-phone to land on the  
dedicated page on the website.  
You might be required to install a  
free application to navigate.

Alternatively, go to: <https://www.meber.it/en/products/stretchers-supports.aspx>



## AMBULANCE COTS SUPPORTS

# MTX

SINGLE SPACE SUPPORT WITH LATERAL MOVEMENT



Art. 19152-130

## MTX-AR-130

Reduced height ambulance cot support with lateral movement (no space underneath), with 130 mm ramp



Art. 19130-955

## MTX-955

Single-space cot support with lateral movement, 180 mm ramp and Art. 955 fastening system – **Certified EN1789**

Art. 19130-956

## MTX-956

Single-space cot support with lateral movement, 180 mm ramp and Art. 956 fastening system – **Certified EN1789**

Art. 19130

## MTX

Single-space cot support with lateral movement, 180 mm ramp – **Certified EN1789**

# MTX-2

DOUBLE-SPACE AMBULANCE COT SUPPORT WITH LATERAL MOVEMENT



Double-space ambulance cot support with lateral movement and 350 mm ramp. The underneath space may house both a scoop stretcher and a spinal board. Fastening systems 955 and 956 can be integrated.



Art. 19140

## MTX-2

Double-space ambulance cot support with lateral movement and 350 mm ramp – **Certified EN1789**

Art. 19140-955

## MTX-2-955

Double space ambulance cot support with lateral movement and 350 mm ramp, fastening system Art. 955 – **Certified EN1789**

Art. 19140-956

## MTX-2-956

Double space ambulance cot support with lateral movement and 350 mm ramp, fastening system Art. 956 – **Certified EN1789**

Art.	Size with folded ramp (Width x Height x Length)	Size with ramp (Width x Height x Length)	Weight	Load capacity	Lateral movement
19130	60,5 x 19,5 x 214,5 cm	60,5 x 37 x 202 cm	58 kg	400 kg	30 cm
19130-955	60,5 x 19,5 x 214,5 cm	60,5 x 37 x 202 cm	66,5 kg	400 kg	30 cm
19130-956	60,5 x 19,5 x 214,5 cm	60,5 x 37 x 202 cm	66,5 kg	400 kg	30 cm
19152-130	60,5 x 13,5 x 208 cm	60,5 x 31 x 202 cm	42 kg	400 kg	30 cm
19140	63 x 27,5 x 229 cm	63 x 62 x 202 cm	76 kg	400 kg	30 cm
19140-955	63 x 27,5 x 229 cm	63 x 62 x 202 cm	84,5 kg	400 kg	30 cm
19140-956	63 x 27,5 x 229 cm	63 x 62 x 202 cm	84,5 kg	400 kg	30 cm

# MX | MX-2

SINGLE-SPACE AND DOUBLE-SPACE FIXED SUPPORT  
FOR AMBULANCE COTS

Art. 19136

Art. 19146



Art. 19136

## MX

Single space fixed support with 180 mm ramp - **Certified EN1789**

Art. 19146

## MX-2

Double-space fixed support with 350 mm ramp - **Certified EN1789**

Art. 19136-955

## MX-955

Single space fixed support with 180 mm ramp and fastening system Art. 955 - **Certified EN1789**

Art. 19146-955

## MX-2-955

Double-space fixed support with 350 mm ramp and fastening system Art. 955 - **Certified EN1789**

Art. 19136-956

## MX-956

Single space fixed support with 180 mm ramp and fastening system Art. 956 - **Certified EN1789**

Art. 19146-956

## MX-2-956

Double-space fixed support with 350 mm ramp and fastening system Art. 956 - **Certified EN1789**

Art.	Size with folded ramp (Width x Height x Length)	Size with ramp (Width x Height x Length)	Weight	Load capacity
19136	60,5 x 13,5 x 208,5 cm	60,5 x 37 x 202 cm	41,5 kg	400 kg
19136-955	60,5 x 13,5 x 208,5 cm	60,5 x 37 x 202 cm	50 kg	400 kg
19136-956	60,5 x 13,5 x 208,5 cm	60,5 x 37 x 202 cm	50 kg	400 kg
19146	63 x 21,5 x 223 cm	60,5 x 31 x 202 cm	59,7 kg	400 kg
19146-955	63 x 21,5 x 223 cm	63 x 62 x 202 cm	68,2 kg	400 kg
19146-956	63 x 21,5 x 223 cm	63 x 62 x 202 cm	68,2 kg	400 kg
19156-130	62 x 24 x 205 cm	63 x 62 x 202 cm	67 kg	400 kg
19156-130-AL	62 x 24 x 205 cm	63 x 62 x 202 cm	47 kg	400 kg
19162	77,5 x 27 x 201 cm		75 kg	400 kg



# MTX-2R

REDUCED HEIGHT DOUBLE-SPACE COT SUPPORT  
WITH LATERAL MOVEMENT



Art. 19156-130

Art. 19156-130-AL



Art. 19156-130

## MTX-2R

Reduced height double space-cot support with  
130 mm ramp

Reduced height double space cot support with 130 mm ramp. The space underneath may house both a scoop stretcher and a spinal board, despite the reduced height (~40 mm) compared to the MTX-2. **Certified EN1789**. Max lateral movement 30 cm.

Art. 19156-130-AL

## MTX-2R AL

Reduced height double-space cot support in  
aluminium with 130 mm ramp

Sharing the same features of the MTX-2R, this support presents an aluminium structure, reducing its weight of approximately 20 kg. **Certified EN1789**. Max lateral movement 30 cm.

Ambulance cots supports

# MTX-2-X



DOUBLE SPACE COT SUPPORT WITH LATERAL MOVEMENT  
FOR MEBER-X STRETCHER

Art. 19162

## MTX-2-X

Double-space cot support  
with lateral movement with  
special MEBER-X system with  
holes, no ramp



# HERCULES

BARIATRIC AMBULANCE COT SUPPORT



Ambulance cots supports



Bariatric

Art. 19170

## MTX-2RH-AL

Double-space ambulance space cot support in aluminium, with reduced height, lateral movement, HERCULES system (stainless steel) and 180 mm ramp with roller



Bariatric

Art. 19172

## MTX-2RPH-AL

Double space ambulance cot support in aluminium, with reduced height, comb lateral movement, HERCULES system (stainless steel) and 180 mm ramp with roller

200

The HERCULES is the new ambulance cot support with motor-powered loading and unloading.

The integrated position detectors and small pushbutton ensure a safe and easy use of this support. The loading takes place by means of a quick-release buckle system with belt, actioned by an integrated coaxial motor. No outer volume. Equipped with a double-safety system, to ensure total safety during loading and unloading). Can be installed on MTX-2R and MTX-2R-AL. 130 mm ramp. The underneath space may house both a scoop stretcher and a spinal board, despite the reduced height (~40 cm) compared to the MTX-2.

**Certified EN1789. 30 cm lateral movement.**

Art.	Length	Width	Height	Weight	Load Capacity	Lateral movement
19172	216,5 cm	57,8 cm	29 cm	65 kg	400 kg	30 cm (15 + 15)
19170	216,5 cm	57,8 cm	29 cm	65 kg	400 kg	30 cm (15 + 15)



Art. 19170 + 7270 + 12084 +  
12170 + 12172 + 7013/N

### High efficiency motor

For an efficient operator  
support during the  
loading of bariatric  
patients.

Double-space  
Ergon/spinal board/  
scoop stretcher



Pushbutton

Ambulance cots supports

201

Product Catalogue 2023

# ACCESSORIES

**MTX** AND **MX** AMBULANCE COT SUPPORTS



**Art. 19110**

**Stainless steel plate mm 1917 x 572 with connections for MTX, MTX-2 and MX.**



**Art. 19112**

**100 mm ramp with support connections.**



**Art. 19114**

**130 mm ramp with support connections.**



**Art. 19116**

**180 mm with support connections.**



**Art. 19118**

**350 mm ramp with connections and gas spring.**



**Art. 787L**

**Gas spring for support's side.**

**NEW  
PRODUCT**

## AMBULANCE COT SUPPORTS WITH LATERAL MOVEMENT



**Art. 19115**

**180 mm ramp with lateral sides and roller for lateral movement**



**Art. 19123**

**130 mm ramp with lateral sides and roller for lateral movement**

**NEW  
PRODUCT**

**NEW  
PRODUCT**

Art.	Length	Width	Height	Weight
19110	194 cm	59 cm	2,5 cm	16,3 kg
19112	10 cm	62 cm	1,6 cm	1,4 kg
19114	13 cm	61 cm	1,7 cm	1,7 kg
19115	18 cm	62,5 cm	1,7 cm	2,8 kg
19116	18 cm	63 cm	1,7 cm	2,3 kg
19118	35 cm	62 cm	5,4 cm	4,3 kg
19123	13 cm	60,5 cm	1,7 cm	2,3 kg




# VETERINARY



**Art. 6100**

## **Foldable stretcher (long) for veterinary use**

Foldable stretcher with foot supports and handles, built in aluminium, with tearproof, water-proof PVC canvas. Ideal for transfers. Easily stored. Two available sizes



**Art. 6150**

## **Foldable stretcher (short) for veterinary use**

Foldable stretcher with foot supports and handles, built in aluminium, with tearproof, water-proof PVC canvas. Ideal for transfers. Easily stored. Two available sizes



**Art. 6200/V**

## **Scoop stretcher for veterinary use**

Entirely built in aluminium, this stretcher can be detached, ensuring a safe loading of the patient, reducing movement. Equipped with two quick-release buckle straps

Art.	Length	Canvas length	Width	Height	Folded height	Depth	Weight
6100	167 cm	140 cm	50 cm	13 cm	83 cm	13 cm	5,5 kg
6150	127 cm	100 cm	50 cm	13 cm	63 cm	13 cm	5,5 kg
6200/V	130 cm	-	42 cm	7 cm	-	-	7,8 kg
6300/B	100 cm	-	65 cm	-	-	-	3,5 kg



**Art. 6300/B**

## **Vacuum mattress for veterinary use**


This article ensures a precise immobilization of the patient, thus easing transfers. Built in PVC, this article does not absorb moisture and is easy-to-clean. Equipped with 4 handles for transport



**Art. 896/DF**

## **Suction pump**

ABS suction pump for vacuum mattress. Can be used to inflate



**Art. 6315**

## **6-handles transfer sheet**

Animal transfer sheet, made in water-proof, fireproof, tear-proof nylon.  
**Size 75 x 100 cm.**





# GENERAL SALE CONDITIONS

## ORDERS

---

For a faster and more accurate service please send orders by fax, letter or e-mail and make sure that all your data, including V.A.T. Registration Number, are included. Please also specify code and description of the articles.

## PRICE

---

Prices are subject to variation at any time without previous notice. Prices are ex works Fontanini – Italy and do not include V.A.T.

## DELIVERY

---

Goods travel at customer's risk.

## PACKING

---

Standard packing is supplied free of charge.

## WARRANTY

---

The articles are guaranteed for 12 months. The warranty covers damages and/or breakages by defects in the materials used, defects in the working or imperfect assembly, furthermore it expires if the products are not used in the correct way, if they are altered and/or repaired by personnel not authorized by Me.Ber.

## PAYMENT

---

Payments must always be in advance or on delivery. Other terms of payment may be arranged with our sales department or with our agents in your area. In the event of outstanding payments, the current bank interest plus a fixed amount of 26,00 € will be added to the due amount.

## CLAIMS

---

Claims are accepted within 8 days from receipt of the goods.

## RETURNS

---

Returns will be accepted only upon authorization of our sales department and free of transport charges..

## CONTROVERSIES

---

In case of any controversy, Parma is the competent court..

Texts, articles, technical data and pictures are subject to variation at anytime without notice.







Stretchers

Self loading stretchers and trolleys

Sedan chairs

Accessories

Immobilization

Diagnostics

Resuscitation

First aid equipments

Fastening systems

Restraint systems

Supports

Veterinary line

Pediatric products

10g certified products

# MEBER<sup>+</sup>

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy

Ph. +39 0521 648770

Fax +39 0521 648780

e-mail: [info@meber.it](mailto:info@meber.it)

