



TECHNICAL DATASHEET

Item code	12084	12086	12088
Name	Set	Bariatric wing for self loading stretchers with fixing kit	
		Right	Left
Dimensions [mm]			
Length	-	1360	1360
Width	-	270	270
Table width	-	195	195
Height	-	105	105
Used materials	Stainless steel, Al		
Weight [kg]	14	7	7
Declared use	Aid to transport, in horizontal position, of bariatric patients temporarily or permanently unable, whose pathologies are compatible; only on ME.BER stretchers MBS-01, Frog and Mercury		
Conformity	MDR 2017/745 - Class I Medical Device		



2305

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it