



TECHNICAL DATASHEET

Item code	13018	
Name	CRONOS - Multi-level stretcher with backrest and trendelenburg	
Dimensions [mm]	Length	1980
	Width	590
	Height	480 ÷ 860
Used materials	AISI 304, Al alloy, HDPE, ABS	
Wheels	n°2 Ø200 mm (rotating w/brake) n°2 Ø200 mm (fixed)	
Weight [kg]	37	
Max load [kg]	170	
Declared use	Manual transportation of recumbent patients, with compatible pathologies, in hospitals, nursing homes and rest homes	
Conformity	MDR 2017/745 – Class I Medical Device	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it