



TECHNICAL DATASHEET

Item code	16100	
Name	TOBOGA Basket stretcher with belts and footrest	
Dimensions [mm]	Length	2160
	Width	610
	Height	170
Used materials	Aluminium, ABS, EVA foam	
Weight [kg]	15	
Max load [kg]	300	
Temperature range [°C]	-25 ÷ 60	
Declared use	Transport patients to conventional cots or stretchers from locations that are difficult to access	
Conformity	MDR 2017/745 – Class I Medical Device	



2212

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it