



## TECHNICAL DATASHEET

Item code	16102		
Name	TOBOGA Detachable basket stretcher with belts and footrest		
Dimensions [mm]	Length	2160	
	Length divided	1080 (1140)	
	Width	610	
	Height	170	
Used materials	Aluminium, ABS, EVA foam		
Weight [kg]	15		
Max load [kg]	300		
Temperature range [°C]	-25 ÷ 60		
Declared use	Transport patients to conventional cots or stretchers from locations that are difficult to access		
Conformity	MDR 2017/745 – Class I Medical Device		



2212

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)