




## TECHNICAL DATASHEET

<b>Item code</b>	<b>6200</b>	
<b>Name</b>	<b>MAXIMA Pediatric anti trauma scoop stretcher</b>	
<b>Dimensions [mm]</b>	Length	1300
	Width	420
	Height	70
<b>Used materials</b>	Al alloy	
<b>Weight [kg]</b>	7,8	
<b>Max load [kg]</b>	170	
<b>Declared use</b>	Transport, in horizontal position, of pediatric patients temporarily or permanently unable, whose pathologies are compatible	
<b>Conformity</b>	MDR2017/745 – Class I Medical Device	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)