

TECHNICAL DATASHEET

Item code	Name	
631	Grey anodized anti trauma scoop stretcher with head/trunk in 1 pc	
631/G	Yellow painted anti trauma scoop stretcher with head/trunk in 1 pc	
Dimensions [mm]	Length	1670 ÷ 2200
	Folded length	1200
	Width	420
	Height	70
	Folded height	90
Used materials	Aluminium alloy	
Weight [kg]	9,5	
Max load [kg]	170	
Declared use	Transport, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible	
Conformity	MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 Certified	
Standard equipment	Bag with n°3 two pieces belts (603/MEB)	



2104

All rights are reserved. Any variations on images and subjects can be done without notice.
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it