

TECHNICAL DATASHEET

Item code	638
Name	Folding rescue stretcher with adjustable backrest
Dimensions [mm]	Length 1870
	Folded length 980
	Width 490
	Height 140
	Thickness folded 170
Used materials	Al alloy, Nylon+PVC
Weight [kg]	9,2
Max load [kg]	170
Equipment	n°2 wheels Ø100 mm n°2 feet
Declared use	Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.
Conformity	MDR 2017/745 – Class I Medical Device



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it