


## TECHNICAL DATASHEET

<b>Item code</b>	<b>640</b>
<b>Name</b>	<b>Rescue stretcher without backrest</b>
<b>Dimensions [mm]</b>	Length 1860
	Folded length 920
	Width 480
	Height 120
	Thickness folded 170
<b>Used materials</b>	Al alloy, Nylon+PVC
<b>Weight [kg]</b>	6,5
<b>Max load [kg]</b>	170
<b>Equipment</b>	n°4 feet n°4 handles
<b>Declared use</b>	Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device



2207

All rights are reserved. Any variations on images and subjects can be done without notice.  
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: [info@meber.it](mailto:info@meber.it)