

## TECHNICAL DATASHEET

Item code	643	
Name	Duo fold rescue stretcher	
Dimensions [mm]	Length	1950
	Folded length	970
	Thickness folded	170
	Width	490
	Folded width	120
	Height	120
Used materials	Al alloy, Steel, Nylon+PVC	
Weight [kg]	5	
Max load [kg]	170	
Equipment	n°4 feet n°4 handles	
Declared use	Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.	
Conformity	MDR 2017/745 – Class I Medical Device	



2207

All rights are reserved. Any variations on images and subjects can be done without notice.  
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: [info@meber.it](mailto:info@meber.it)