

TECHNICAL DATASHEET

Item code	644
Name	Stiff fixed stretcher for stowage
Dimensions [mm]	Length 1860
	Width 480
	Height 170
Used materials	Al alloy, steel, Nylon+PVC
Weight [kg]	4,9
Max load [kg]	170
Equipment	n°4 feet
Declared use	Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.
Conformity	MDR 2017/745 – Class I Medical Device



2207

All rights are reserved. Any variations on images and subjects can be done without notice.
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it