


TECHNICAL DATASHEET

| | | |
|-----------------|---|------|
| Item code | 646 | |
| Name | Folding stiff stretcher | |
| Dimensions [mm] | Length | 2010 |
| | Folded length | 1010 |
| | Width | 500 |
| | Height | 140 |
| | Thickness folded | 170 |
| Used materials | Al alloy, steel, Nylon+PVC | |
| Weight [kg] | 6,2 | |
| Max load [kg] | 170 | |
| Equipment | n°4 feet | |
| Declared use | Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible. | |
| Conformity | MDR 2017/745 – Class I Medical Device | |



2207

All rights are reserved. Any variations on images and subjects can be done without notice.
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it