

## TECHNICAL DATASHEET

<b>Item code</b>	<b>660</b>
<b>Name</b>	<b>4 pieces folding stretcher with telescopic handles</b>
<b>Dimensions [mm]</b>	Length 1920 ÷ 2160
	Folded length 940
	Width 570
	Folded width 110
	Height 180
	Thickness folded 210
<b>Used materials</b>	Al alloy, steel, Nylon+PVC
<b>Weight [kg]</b>	9,3
<b>Max load [kg]</b>	170
<b>Equipment</b>	n°4 feet n°4 handles
<b>Declared use</b>	Transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device



2207

All rights are reserved. Any variations on images and subjects can be done without notice.  
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: [info@meber.it](mailto:info@meber.it)