

TECHNICAL DATASHEET

Item code	7062/4RG PROOF	
Name	TITANUS - Stainless steel variable height trolley certified EN 1865 high version with 4 rotating wheels	
Dimensions [mm]	Length	1970
	Folded length	2000
	Width	440
	Rear axle width	560
	Max height	890
	Intermediate height	610
	Folded height	250
	Loading height	680 ÷ 720
Used materials	AISI304, PVC, HDPE	
Wheels	n°4 Ø200 mm (2 w/brake) n°2 Ø125 mm n°2 Ø100 mm n°2 Ø60 mm	
Weight [kg]	37	
Max load [kg]	220	
Declared use	Transport of patients in a relaxed position and loading onto emergency vehicles of patients with compatible pathologies, together with fixed stretchers manufactured by ME.BER.	
Conformity	MDR 2017/745 – Class I Medical Device EN1865-1 – EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice.
Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it