



TECHNICAL DATA SHEET

Item code	7080/4RG PROOF	
Name	MERCURY LITE – Aluminium alloy self loading stretcher with variable height with 4 swivel wheels	
Dimensions [mm]	Length	1970
	Length with folded trolley	1650
	Width	570
	Max height	950
	Intermediate height	720
	Folded height	310
	Loading height	640 ÷ 660
Weight [kg]	40	
Max load [kg]	250	
Used materials	Al alloys, HDPE	
Wheels	n°4 Ø200 · n°2 Ø150 · n°2 Ø60	
Characteristics	<ul style="list-style-type: none"> • Backrest inclination 0 ÷ 75° (h. p.) • Trendelenburg: inclination 0 ÷ 20° 	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled.	
Conformity	MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified	



2405

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it