



## TECHNICAL DATASHEET

<b>Item code</b>	<b>7092 PROOF</b>	
<b>Name</b>	<b>MERCURY EL - Black painted steel cot certified EN 1865</b>	
<b>Dimensions [mm]</b>	Length	1970
	Width	570
	Max height	950
	Folded height	310
	Loading height	640 ÷ 660
<b>Used materials</b>	S235JR, ABS, HDPE	
<b>Wheels</b>	n°4 Ø200 mm (2 fixed+2 swivel) n°2 Ø150 mm n°2 Ø60 mm	
<b>Weight [kg]</b>	40	
<b>Max load [kg]</b>	250	
<b>Declared use</b>	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device EN 1865-1 - EN 1789 certified	



2301

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)