



TECHNICAL DATASHEET

Item code	7094 PROOF	
Name	MERCURY EL - Black painted steel cot high version certified EN 1865	
Dimensions [mm]	Length	1970
	Width	570
	Max height	990
	Folded height	310
Dimensions [mm]	Loading height	690 ÷ 710
Used materials	S235JR, ABS, HDPE	
Wheels	n°4 Ø200 mm (2 fixed+2 swivel) n°2 Ø150 mm n°2 Ø60 mm	
Weight [kg]	41	
Max load [kg]	250	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
Conformity	MDR 2017/745 – Class I Medical Device EN 1865-1 - EN 1789 certified	



2301

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it