



TECHNICAL DATASHEET

Item code	7095/4RG PROOF	
Name	"MERCURY LITE 5" Aluminium self loading stretcher, 5 positions with 4 rotating wheels EN 1865 certified	
Dimensions [mm]	Length	1970
	Length with folded trolley	1670
	Width	570
	Height	950 – 870 – 790 – 685 – 390
	Folded height	310
	Loading height	640 ÷ 660
Used materials	Aluminium alloys	
Wheels	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm	
Weight [kg]	41	
Max load [kg]	250	
Declared use	Transporting, in the lying down position, of temporarily or permanently disabled patients with pathologies that allow this. The stretcher is suitable for transporting patients on vehicles appropriately prepared with fastening systems	
Conformity	MDR 2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified	



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it