



## TECHNICAL DATASHEET

Item code	7122/AZ	7122/V
Name	GECKO - Self loading stretcher entry-level	
	Blue	Green
Dimensions [mm]	Length	1970
	Width	560
	Height	905
	Folded height	300
	Loading height	620
Used materials	S235JR, Al alloys, HDPE	
Wheels	n°4 Ø200 mm (2 fixed+2 rotating) n°2 Ø125 mm	
Weight [kg]	37	
Max load [kg]	170	
Declared use	Transporting, in the lying down position, temporarily or permanently disabled patients with pathologies that allow this	
Conformity	MDR 2017/745 – Class I Medical Device	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)