



TECHNICAL DATASHEET

Item code	7210 PROOF		
Name	FROG – Yellow self-loading stretcher with Trendelenburg certified		
Dimensions [mm]	Length	1950	
	Length with folded trolley	1660	
	Width	560	
	Max height	920	
	Folded height	325	
	Loading height	600 ÷ 660	
Used materials	S235JR, ABS, HDPE		
Wheels	n°4 Ø200 mm (2 fixed+2 rotating) n°4 Ø125 mm n°2 Ø100 mm		
Weight [kg]	44		
Max load [kg]	250		
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled		
Conformity	MDR 2017/745 – Class I Medical Device EN 1865-1 – EN 1789 certified		



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it