



TECHNICAL DATASHEET

Item code	7210
Name	FROG – Yellow self-loading stretcher with trendelenburg
Dimensions [mm]	Length 1950
	Length with folded trolley 1690
	Width 560
	Max height 920
	Folded height 325
	Loading height 600 ÷ 660
Used materials	S235JR, ABS, HDPE
Wheels	n°4 Ø200 mm (2 fixed+2 rotating) n°4 Ø125 mm n°2 Ø100 mm
Weight [kg]	43
Max load [kg]	250
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled
Conformity	MDR 2017/745 – Class I Medical Device



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it